


NeoBellyBand

Goal: a boundary on the tummy will help her use her tummy muscles and give them support to prevent further expansion. This can potentially decrease her belly size.

Wear: Please open during cares to check skin, and then refasten it.

FIT: Please place with lower border at plane of hip points (ASIS) and upper border at xiphoid process.



Place your hand under her low back to be sure it is long and rounded before closing. **The belt is fastened at her girth.** Draw the sides in by pulling straight up and then lay over. **Do not squeeze her in or overtighten.** The goal is support NOT compression.

It can be over the diaper. If soiled, there is a second in the room, hand wash and hang to dry, stretch to reshape. We will try the belt for a few days and re-evaluate next week, often it stays in use until coming off CPAP.

**Please leave a message with therapy if the belt is being removed for any reason with the time and reason for removal.*