

MULTIPULSE HOPLUS

WHEN POWERFUL VERSATILITY MEETS
UNMATCHED PRECISION



JENA SURGICAL
LASER AT YOUR SIDE

POWER IN YOUR HANDS

MULTIPULSE HOPLUS

The **MultiPulse HoPLUS** is a high power holmium laser device characterized by great flexibility and efficiency. Ideal for both, surgical treatment of benign prostatic hyperplasia (BPH) and endoscopic treatments for lithotripsy it is highly recommended for long-lasting procedures requiring high power thanks to its laser source reaching 150 W.

The **MultiPulse HoPLUS** guarantees superior surgical performance thanks to high levels of energy and pulse repetition rate, guaranteeing reduced operating time and at the same time, effective and delicate cutting procedures with an excellent hemostatic effect.

ADVANTAGES FOR THE SURGEON

- High power laser
- High pulse repetition rate
- Fast vaporization and cutting with excellent hemostasis
- Ability to immediately remove the enucleated tissue with the integrated morcellator available as an optional accessory
- High versatility in endourology
- Dual pedals to select the required mode without changing the laser fiber and without looking away from the operating field, e.g. for cutting and coagulation
- Easy to install and operate

ADVANTAGES FOR THE PATIENT

- Minimally invasive surgery
- Low probability of bleeding: recommended also for patients treated with anticoagulants
- Low rate of side effects
- Shorter catheter application time in the post-operative phase
- Shorter hospital stay and recovery time
- Low rate of recurrences
- Immediate relief of symptoms

MULTIPULSE HOPLUS

ASPIRATION PUMP
for integrated morcellator
(optional)

PORT
for the blast shield holder

IMPROVED COOLING
for low noise
Holmium laser


RFID
recognition
system

FACILITATED MAINTENANCE
due to sliding covers

2 KW
less power
consumption



MULTIPULSE HOPLUS WORLDWIDE



OUR SERVICE REGIONS

Asclepion creates the latest generation workstations used in urology, gynecology, oncology and otolaryngology for surgical teams all around the world.

For our portfolio of high quality "Made in Germany" laser devices, Asclepion ensures high availability worldwide as well as high performance service care through our global distributor network.

VERSATILITY

150 W / 100 HZ
no upgrade
230V ~50/60Hz; 32A
integrated Morcellator
optional



MULTIPULSE HOPLUS

PREMIUM

POWERKIT⁵⁰

additional 50W of power for the high power Holmium laser with 100W and 100Hz easy and fast to upgrade



100 W / 100 HZ
upgradable to 150W/100Hz
230V ~50/60Hz; 32A
integrated Morcellator optional

POWERKIT²⁰

additional 20W of power for the high power Holmium laser with 80W and 75Hz easy and fast to upgrade



100 W / 75 HZ
no upgrade
230V ~50/60Hz, 32A
no integrated Morcellator

80 W / 75 HZ
upgradable to 100W/75Hz
230V ~50/60Hz, 32A
no integrated Morcellator

UPGRADES

BASIC

THE MULTICUT

MULTICUT

The MultiPulse HoPLUS laser system contains an **integrated morcellation module***, which allows the morcellation of enucleated adenoma within the bladder to be performed very quickly. Thanks to the peristaltic pump, the aspiration is highly efficient and additionally ensures a maximum of safety and protection against unintended bladder wall injuries.

* for the 100 Hz system

ADVANTAGES

- Ergonomic design, independent right or left handed use
- Light and easy to clean handpiece
- No consumables except for blades
- Excellent blade design, two different sizes
- New and facilitated coupling – *Plug and Cut*
- Optimized cutting teeth geometry guarantees enough tissue for a histopathological examination
- Suction of the peristaltic pump is highly efficient and additionally ensures a maximum of safety and protection against unintended bladder wall injuries
- Possibility to always perform the histopathological examination of enucleated tissue samples
- Unique combination of laser and morcellation system in one device (optional)

**CUTTING TEETH
GEOMETRY**
with rotating Blades



INTEGRATED MORCELLATION MODULE

Accessories | Tube set, morcellation handpiece with blades

Dimensions | 15 Fr x 350 or 400 mm (diameter x free length)

Aspiration pump flow rate | Adjustable


Oscillation frequency | 0 or 2 – 5 Hz

Rotation speed | 500 – 3000 rpm

LIGHT & EASY
ergonomic design

NEW
Coupling system

CLINICAL CASES



MORCELLATION, 250 CC ADENOMA, APPROX. 6 MINUTES

[Courtesy of Dr. Edmundo Guarnizo Olivera - Clinica Monterrico - Perú]

APPLICATIONS AND CLINICAL CASES

The **MultiPulse HoPLUS** is undoubtedly the most versatile laser system for urology. It offers a wide set of applications ranging from treatment of strictures (urethrotomy) or bladder neck incision (BNI), but also being universally recognized as the gold standard in endoscopic lithotripsy for stone fragmentation and pulverization, regardless of their composition and anatomical location.



UROLOGY (HOLEP, BNI, URETHROTOMY, LITOTRIPSY)

The MultiPulse HoPLUS is best suited for the surgical treatment of benign prostatic hyperplasia and endoscopic lithotripsy. It is recommended for long-lasting use thanks to its power of up to 150 watts.



GYNECOLOGY (TREATMENT OF VVF)

The Holmium (Ho:YAG) laser wavelength is intended for ablation, vaporization, incision, excision, welding and coagulation of soft tissue. Treating VVF (vesicovaginal fistula) can be done with the MultiPulse HoPLUS.





ENT (FRAGMENTATION OF SALIVARY DUCT STONES)

Sialolithiasis can cause blockage of salivary ducts resulting in painful inflammation. Most salivary stones occur in the submandibular gland, followed by the parotid gland and infrequently in the sublingual or minor salivary glands. While smaller stones may pass on their own, larger stones generally require medical or surgical intervention.

Minimally invasive and nonsurgical techniques of treating salivary stones have been evolving rapidly. One of the possible treatments is laser assisted minimal invasive therapy. The Holmium Laser Lithotripsy offers a useful tool to destruct those salivary stones.



GENERAL SURGERY (GASTROENTEROLOGY)

Bile duct stones (Cholelithiasis) are gallstones that are present in the bile ducts, sometimes within the cystic duct of the gallbladder, but more frequently in the common bile duct. Ho:YAG laser lithotripsy can be done via percutaneous or peroral access to the biliary tree under an appropriate direct visual control. Once the stones are fragmented, the debris is flushed out the access catheter or may pass through the sphincter of Oddi into the duodenum.


UROLOGY

Prostate enucleation (BPH treatment) | Urinary lithotripsy | Urethrotomy

BENIGN PROSTATIC HYPERPLASIA

The MultiPulse HoPLUS enables the enucleation of the prostate as a minimally invasive alternative to the traditional techniques of prostatectomy such as transurethral resection (TURP) and open



surgery. HoLEP, together with the presence of the integrated mechanical morcellator (optional), makes it possible to reduce operating time and treat prostates of any size.



The MultiPulse HoPLUS does not produce carbonization or necrosis, leaving the surgical field free and clearly visible. The holmium laser emits a light of 2,100 nm wavelength, which is well absorbed in water, resulting in photothermal and photomechanical effects which lead to effective vaporization and good hemostasis within the tissue.

DIFFERENT STEPS OF HOLEP FOR THE TREATMENT OF BENIGN PROSTATIC HYPERPLASIA (BPH), PERFORMED USING THE MULTIPULSE HOPLUS LASER SYSTEM.

[Courtesy of S. Piesche, MD - Clinic of Urology and Urological Oncology, Sana Klinikum Hof – Germany]



DIFFERENT STEPS OF HOLEP (BPH TRILOBULAR, BIG VOLUME) PERFORMED USING THE MULTIPULSE HOPLUS LASER SYSTEM.

[Courtesy of Dr. Edmundo Guarnizo Olivera - Clinica Monterrico - Perú]



DIFFERENT STEPS OF HOLEP (BPH, 360 CC) PERFORMED USING THE MULTIPULSE HOPLUS LASER SYSTEM.

[Courtesy of Dr. Edmundo Guarnizo Olivera - Clinica Monterrico - Perú]



ADVANTAGES

- Use of saline not glycine – long operations on large prostates (>200 gm) possible
- No TURP syndrome
- Identical surgical cavity as TURP
- Urodynamic outcomes identical to TURP
- Enough tissue for histology
- No electromagnetic interaction with patient
- Minimal bleeding
- Minimal risk of transfusion
- Shorter catheter time, shorter hospital stays (outpatient/24 hour stay)
- Fewer nursing requirements, lower costs


UROLOGY

Prostate enucleation (BPH treatment) | Urinary lithotripsy | Urethrotomy

LITHOTRIPSY

The MultiPulse HoPLUS laser has the capability of fragmenting calculi but can also ablate tissue. The Ho:YAG laser (2,100 nm) is very rapidly absorbed in water and thus has a thermal effect. It produces a very small steam bubble which quickly expands and contracts. The water content of uroliths and other calculi turned out to be high enough for adequate energy absorption; therefore calculi can be destroyed regardless of their color, hardness or composition.

RETROGRADE INTRARENAL SURGERY for staghorn calculus occurring in the renal pelvis, with branches extending into the medium and inferior calyces. RIRS was performed with flexible ureterorenoscopy and using MultiPulse HoPLUS laser with 200 µm fiber.



[Courtesy of Prof. M. S. Minervini, M.D. - Head of Urological Division - Hospital of Sondrio, Italy]

BLADDER STONE TREATMENT (pulverization and fragmentation) with the MultiPulse HoPLUS.



[Courtesy of S. Piesche, MD - Clinic of Urology and Urological Oncology, Sana Klinikum Hof – Germany]

DIFFERENT STEPS OF REMOVING BLADDER STONES performed by using the MultiPulse HoPLUS laser system.



[Courtesy of Dr. med. Martin Kanne – Königin Elisabeth Herzberge Krankenhaus – Germany]

LITHOTRIPSY BENEFITS WITH THE MULTIPULSE HOPLUS

- Small fibers used with either flexible or rigid scopes
- Minimal stone movement
- Virtually no bleeding
- Effective on all stone compositions - unlike EHL, ultrasound or lithoclast
- Efficient in removing embedded stones





UROLOGY

URETHROTOMY - TREATMENT OF URETHRAL STRICTURES

Urethrotomy is an operation for men who have a decreased urinary stream and problems passing urine because of a bladder neck stenosis or a urethral stricture. The aim of the urethrotomy is to cut through these bands of tissue to resolve


the narrowing and allow for an improved urinary stream and relief of symptoms. The narrowing caused by a urethral stricture can cause bladder obstruction that leads to incomplete bladder emptying.

ADVANTAGES

- Minimally invasive and safe procedure suitable for short urethral strictures and for more significant strictures or for the treatment of urethral atresia
- Ablation of damaged tissue with a clean cut and without penetrating too deeply
- Virtually no bleeding procedure
- Reduces the risk of lateral thermal damage, complications, recurrences or the formation of post-surgery fibrosis
- Resulting in complete restoration of normal urinary flow

TECHNICAL SPECIFICATIONS

| MultiPulse HoPLUS | |
|---|---|
| Laser Source | Ho:YAG |
| Wavelength | 2,100 nm |
| Emission Mode | Pulsed Wave (pw) |
| Power | up to 150 W (depending on the configuration) |
| Pulse Energy | 0.25 - 6 J |
| Repetition Rate | Up to 100 Hz |
| Pulse Duration | 150 - 1,700 μ s |
| Cooling | Internal water cooling |
| Control Panel | 10.1" LCD display 16:9 |
| Integrated Morcellation Module (Optional) | Integrated auxiliary morcellator Ergonomic handpiece Exchangeable blades |
| Device Accessories | Fiber Handpieces and Cannulas Bare fibers (reusable and single use) available in following diameters: 200, 272, 365, 400, 550, 600, 800, 1000 μ m |
| Dimensions and Weight | 116 (H) x 105 (D) x 46 (W) cm - 285 kg |



JENA SURGICAL
LASER AT YOUR SIDE

JenaSurgical is the brand of the surgical
business unit of Asclepion Laser Technologies.
Specifications are subject to change without notice.
Brochure not for the USA



JENA SURGICAL
LASER AT YOUR SIDE

MADE IN GERMANY
ALWAYS THE LATEST PRODUCT INFORMATION



FOLLOW US

Asclepion Laser Technologies GmbH
Brüsseler Str. 10
07747 Jena | Germany

www.jenasurgical.com