


Abstract

Mesenchymal stromal cells (MSCs) has been identified in the 60s and since then, they stole the spotlight, especially in the field of cellular therapy. MSCs are able to communicate with their microenvironment through both juxtacrine and paracrine pathways. Numerous studies show that beneficial effects may mainly be due to the secretome and mostly the extracellular vesicles (EVs) cells derived(1). This opens the door to a new approach which can be named: acellular therapy. MSCs-derived EVs are classically obtained either in a chemically defined media or in a starved condition. Hence, we focus on the determination of better EVs production conditions compatible with therapeutic applications. Our approach is to work both on the media composition and on the culture platform. We first produce EVs-depleted fetal bovine serum (FBS) through tangential flow filtration(2). We combine this EVs-free FBS with the use of wheat peptone as coating while cultivating the MSCs. Our in-house method led us to recover an increased EVs amount compared to more classical methods. Besides, EVs produced with our method exhibit better efficiency on two functional assays: a migration assay and a wound healing assay. Our technique may be useful for upscaling EVs production process which is a mandatory step for therapeutic EVs production.

Key word :MSC / peptone/ TFF/ EV production/ coating

Problematic



Objectives


- Optimize the production of EVs by MSC using peptone which has the capacity to :
- To attach the cells to the surface and stimulate cell growth and biomolecules synthesis.
 - To reduce the utilization of FBS (4)
 - To be a less expensive alternative to FBS
- Comparison between the EVs produced in the presence/absence of peptone.

Experiments

A/EVs Production:


Medium used for the production of EVs:

- Medium 1: MEM No-Coating
- Medium 2: MEM coating with wheat peptone (C-CELL W106)
- Medium 3: MEM + FBS 5% EVs free + No-Coating
- Medium 4: MEM + FBS 5% EVs free + Coating with C-CELL W106



Results

A/Microscopic observations of MSCs:




B/MSCs Characterization:

FACS results

	CD31 FITC	CD34 PE	CD45 FITC	CD73 APC	CD90 FITC	CD105 PE	CD144 PE	CD309 PE
Control	-	-	-	+	+	+	-	-
αMEM NC	-	-	-	+	+	+	-	-
αMEM C	-	-	-	+	+	+	-	-
FBS - NC	-	-	-	+	+	+	-	-
FBS - C	-	-	-	+	+	+	-	-


C/EVs Quantification:




C/ EVs Concentration/Purification:

	Before TFF (EVs/mL)	After TFF (EVs/mL)
αMEM No Coating (MC1 4)	9,00E+08	1,20E+10
αMEM With Coating (MC1 4)	1,10E+09	1,20E+10
αMEM+FBS 5%+No Coating (MC1 4)	4,00E+09	4,00E+10
αMEM+FBS 5%+Coating (MC1 4)	4,20E+09	5,20E+10
αMEM+FBS 5%+No Coating (MC5-6)	5,00E+09	2,70E+10
αMEM+FBS 5%+Coating (MC5-8)	4,40E+09	4,00E+10


D/EVs Characterization:



D/EVs Functionalities



2/Migration assay (ECFC)




B/MSCs Characterization:

FACS: Using the markers:

	CD31 FITC	CD34 PE	CD45 FITC	CD73 APC	CD90 FITC	CD105 PE	CD144 PE	CD309 PE
Control	-	-	-	+	+	+	-	-

C/EVs Concentration/Purification:

Using Tangential Flow Filtration (TFF)



D/EVs Characterization:


- Transmembrane proteins associated to EVs: CD9 – CD63 – CD81
- Cytosolic protein recovered in EVs: ALIX
- Endoplasmic reticulum protein : Calnexin

E/EVs Functionalities:

1/Wound healing (HaCaT: keratinocytes)

These assays involve creating a scratch or a gap in a cell monolayer and monitoring the movement and proliferation of cells to close the wound over time.

2/Migration assay (ECFC)



Discussion

- The presence of FBS in the medium (αMEM) helps the cells to produce more EVs.
- The coating with wheat peptone stimulates the cells to produce more EVs and to stand a longer period on the flask.
- The characteristics of MSC and EV have not changed with wheat peptone (C-CELL W106) coating, moreover it improves the functionalities that has been showed on our two models (i.e., migration & wound healing)

References

- Kou, M., Huang, L., Yang, J. et al. Mesenchymal stem cell-derived extracellular vesicles for immunomodulation and regeneration: a next generation therapeutic tool? *Cell Death Dis* 13, 580 (2022). <https://doi.org/10.1038/s41419-022-05034-x>
- Lorenzini, B., et al. Producing vesicle-free cell culture additive for human cells extracellular vesicles manufacturing. *Journal of controlled release : official journal of the Controlled Release Society*, 355, (2023). <https://doi.org/10.1016/j.jconrel.2023.01.073>
- Van der Valk, J., Bieback, et al. "Fetal bovine serum (FBS): Past – present – future", *ALTEX - Alternatives to animal experimentation*, 35(1), (2018). doi: 10.14573/altex.1705101.
- Lezin, C., Mauduit, P., Uzan, G., & Abdelgawad, M. E. An Evaluation of Different Types of Peptone as Partial Substitutes for Animal-derived Serum in Vero Cell Culture. *Alternatives to laboratory animals : ATLA*, 50(5), (2022). <https://doi.org/10.1177/02611929221122780>