

TRATEK INDUSTRIES INC.

604-820-7754 or 888-850-9691

www.tratek.ca

COURT ORDERED ALCOHOL MONITORING



GUIDELINES TO FOLLOW REGARDING YOUR SCRAM REMOTE BREATH MONITOR



NEVER consume or put in your mouth anything that may contain alcohol, including but not limited to mouthwash, food made with alcohol, or any beverage that you yourself did not prepare. Do not consume anything but water for the 15 minutes leading up to a test, if possible.

NEVER allow the battery of the remote breath device to become depleted, as this is a violation, and will be reported.

NEVER cover your face, or hold the straw during a test.

NEVER submerge the Remote Breath device in any water. This is a violation, and will be reported.

ALWAYS be aware of your environment. It is your responsibility to never expose the Remote Breath device or yourself to alcohol.

ALWAYS take your scheduled or random tests on time. If you miss a test, turn the device on and take a client-initiated test right away.

ALWAYS follow the prompts on the device. If the device asks for a retest, perform the retest immediately when prompted to do so. The device will give you directions and timing for a retest. If a timed retest is necessary, the timer will count down on the screen and audibly prompt for a retest.

ALWAYS be aware of where your Remote Breath device is located. It is your responsibility to have the device available during any prompted test.

- All data of the device is recorded. Deviation from any of the above **MUST** be reported. We do not and cannot make determinations on what is reported or are we responsible for decisions regarding your program participation.
- In the event you cannot complete or take tests due to a health-related issue, it is your responsibility to discuss the issue with your lawyer, and to obtain any necessary medical documentation to support your claim.

Our website is located at www.tratek.ca , feel free to look and learn more about us.

PLEASE CONTACT US FOR ADDITIONAL INFORMATION, OR TO FIND A REP NEAR YOU.

