

FastGrafter™

Autograft Harvesting System



Fully-sterile system for quick and efficient
harvest of autogenous bone

TREACE[®]
Medical Concepts, Inc.

The Future of Hallux Valgus Surgery™

FastGrafter™

Autograft Harvesting System

The FastGrafter™ Autograft Harvesting System is a sterile-packed, single-use device for quick and efficient harvest of cancellous autogenous bone from the calcaneus, distal tibia, and other harvest sites through a minimal incision approach.

- **Fully-sterile system** for quick and efficient harvest of autograft bone
- **Dual-purpose cutting tip** penetrates cortex and morselizes cancellous bone
- **Single-piece harvester** designed to reduce instrumentation and system complexity

FastGrafter™ Harvester



Dual-Purpose Cutting Tip
penetrates cortex and morselizes
autogenous bone

FastGrafter™ Pusher

Calcaneal Autograft Harvest

Key surgical steps*



1.

Incision & Dissection

Make a small incision over the lateral aspect of the calcaneus, posterior and inferior to the peroneal tendon and sural nerve, approximately 3cm below and 2cm posterior to the distal fibula. Use blunt dissection to expose the bone.



2.

Bone Graft Harvest

Insert the **FastGrafter™ Harvester** into a 1/4" Jacobs chuck on a powered driver. Place the cutting tip of the **Harvester** onto the exposed bone surface. Beginning at low speed in the forward direction, advance the **Harvester** to the desired depth.



3. (Optional)

Additional Bone Graft Harvest

To make space for the additional bone graft, insert the **FastGrafter™ Pusher** through the distal tip of the **Harvester** and push the harvested bone towards the proximal end of the cannula. Reinsert the **Harvester** tip into the original harvest site and make a pass approximately 20-30 degrees from the original harvest path.




4.

Removal of Morselized Bone Graft

Disassemble the **Harvester** from the 1/4" Jacobs chuck. Over a sterile container, insert the **Pusher** through the distal tip of the **Harvester** to expel the morselized graft through the proximal opening of the **Harvester**.

Additional Harvest Sites

The FastGrafter™ Autograft Harvesting System can be used percutaneously through multiple harvest sites including the calcaneus, proximal tibia, and distal tibia.



Proximal Tibia

Distal Tibia

Ordering Information

SK27 FastGrafter™ Autograft Harvesting System (7mm)

The Future of Hallux Valgus Surgery™

TREACE[®]
Medical Concepts, Inc.

203 Fort Wade Rd. | Suite 150 | Ponte Vedra, FL 32081 | 904.373.5940 | www.treace.com

See treace.com/patents ©2019 Treace Medical Concepts, Inc. All rights reserved. M547A