



PaceAid

Improving the Quality of the Life You Save



**PaceAid PA03 - INFANT / PEDIATRIC
1-MAN / 2-MAN
VISUAL CPR GUIDANCE DEVICE**

Device Name: PaceAid

Model Number: PA03

Classification: Class 1, 510(k)-exempt

Sterility Status: Nonsterile

Single-use/Reusable: Single-Use

Material Disclosures: Latex-free

Intent for Use: **PaceAid PA03** is a visual guidance device designed to aid anyone in performing CPR in emergency situations by offering real-time visual cues for chest compressions and rescue breaths, aligned with established CPR guidelines, including those published by the American Heart Association (AHA).

Indications for Use: For use on infants and pediatric patients during CPR, in any setting where CPR may be required.

Contraindications: Do not use **PaceAid PA03** in situations where visual guidance is not appropriate or the device cannot be observed by the rescuer.

Do not use PaceAid PA03 in MRI environments. The device is NOT MRI SAFE and must be removed before any MRI procedures.

Warnings

- Contains a button battery. If swallowed, seek immediate medical attention.
- PaceAid is an aid and does not replace CPR training.
- Survival depends on multiple factors beyond device use.

FDA: PaceAid is a Class 1 medical device registered and listed with the FDA. FDA registration means the manufacturer has provided required information to the FDA, but the FDA has not reviewed or approved this product for safety to current good manufacturing practices (CGMPs).

INFORMATION FOR USE CONTINUED:

Directions for Use

Activating PaceAid PA03:

- Locate the overhanging liner at the bottom of the device; activation is automatic.



Applying the Device:

- Place on victim's upper chest, shoulder, arm, or any area within rescuer's direct line of sight. Ensure secure placement before starting CPR.

Adjusting the Mode:

- Use the toggle button to switch between modes:
 - 1-Man CPR (30 compressions: 2 rescue breaths)
 - 2-Man CPR (15 compressions: 2 rescue breaths)

LED Signal Indicators:

- Red Light (Heart Icon): Compression needed
- Green Light (Head Icon): Deliver breath
- Yellow Light: Current operational mode.

Default Mode:

- Device defaults to 1-Man CPR mode when activated. Use toggle to switch modes as needed.

Device Description

- Visual CPR guidance device for infants/pediatrics
- Contains LED indicators for real-time prompts.
- Does not provide clinical decision-making.
- Materials: Latex-free, contains button battery

Storage & Handling: Store in a cool, dry place. Avoid exposure to extreme temperatures.

Disposal Instructions: Dispose of as per local regulations for single-use medical devices and batteries.

Lot Number / Tracking Information: Refer to packaging for lot number and tracking details.

Adverse Event Reporting: Report adverse events to FDA MedWatch.

Labeling Symbols: Definitions for symbols are provided on the packaging.

Additional Information:

- **PaceAid** can be used alongside AED; shock can be delivered without removing.
- May remain in place during transport or follow-up care if victim remains unstable.

 Manufactured by:
Butler Technologies, Inc.
231 West Wayne, Butler, PA 16001
www.butlertechnologies.com

Exclusively for:

P.O. Box 267
Matewan, NJ 07747

 Medical Device



Single Use Only



MR Unsafe



Button Battery

WARNING! This product contains a button battery. If swallowed, a lithium button battery can cause severe or fatal injuries within 2 hours. Keep batteries out of reach of children. If you think battery may have been swallowed or placed inside any part of the body, seek immediate medical attention.



Website: www.paceaid.com

email: info@paceaid.com