

---

# **QUICK START HERE**

## **Wiz by Pocket Neurobics**

### **(i) Product Suite**

- Uwiz-4: 1/2/3/4 eeg channels, pIR HEG
- Q-wiz: 1/2/3/4 eeg chs, eeg Cap, pIR/nIR HEG
- U-wiz: 1/2 eeg channels
- X-wiz: 1/2 dc-coupled chs: nIR near-infrared blood oxygenation biofeedback, plus pIR passive-infrared forebrain temperature (metabolism) biofeedback

### **(ii) Software**

- Windows, MAC, Android:
  - BioSignal.Network web app in browser
- Windows
  - EEGer
  - BioExplorer
  - BioEra
  - Better Fly

### **(iii) General Installation and Startup**

#### *USB Driver Installation*

- On Windows, see [www.biosignal.network/WizDriverForWindows.html](http://www.biosignal.network/WizDriverForWindows.html)
- On MAC, Android, Wiz just works.

#### *Operation*

- Uwiz-4, U-wiz & X-wiz do not have any buttons, plugging them into a USB will cause them to start transmitting. Any mode changes are commanded by the software
- Q-wiz has mode buttons, but in general again, mode changes will be commanded by the software

#### *Modes*

Current Mode is defined by the number of blinks of the Mode LED:

- No blinks - mode to update firmware
- 1 blink - HEG mode: 2 ch, 32sps
- 2 blinks - EEG mode: 2 ch, 128sps
- 3 blinks - ExG wideband (eg EMG): 2 ch, 512sps
- 4 blinks - EEG4 mode: 4ch, 256sps

## *Device & Signal Integrity*

- At startup, a test signal is transmitted on each channel for five seconds
- EEG Electrode Impedance can be measured
  - off-line (ie momentary) on 4 channels
  - in real-time on 2 channels (Ch1 & Ch2)
- EEG variance is tested and evaluated as OK/not OK.
- EEG signals evaluated as in/out of range
- Signal evaluations are reported by
  - reversal on Mode LED blink on Wiz
  - sent to PC application

## *In case of fire*

- *U4/U/X-wiz have pin-hole button on side to sequence through Modes: Program/1/2/3/4 blinks*
- *Pressing pin-hole button (or Mode button on Q-wiz) during Power-up will do factory reset. Latest firmware version can then be uploaded using facility at [BioSignal.Network](http://BioSignal.Network)*

*Questions?*

[BioSignal.Network](http://BioSignal.Network)