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COROAS OU ABUTMENT PERSONALIZADO

DISSILICATO

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COROAS OU ABUTMENT PERSONALIZADO

DISSILICATO

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1. Abi-Rached FO, Fonseca RG, Haneda IG, Almeida-Junior AA, Adabo GL. The effect of different surface treatments on the shear bond strength of luting cements to titanium. *J Prosthet Dent* 2012;108:370-376.
2. Baratto SSP, et al. Silanated Surface Treatment: Effects on the Bond Strength to Lithium Disilicate Glass-Ceramic. *Brazilian Dental Journal* (2015) 26(5): 474-477.
3. Blatz MB, Bergler M, Holst S, Block MS. Zirconia abutments for single tooth implants - rationale and clinical guidelines. *J Oral Maxillofac Surg* 2009;67:74-81.
4. Lim BS, Heo SM, Lee YK, Kim CW. Shear bond strength between titanium alloys and composite resin: sandblasting versus fluoride-gel treatment. *J Biomed Mater Res* 2003;64B:38-43.
5. Long M, Rack HJ. Titanium alloys in total joint replacement – a materials science perspective. *Biomater* 1998;19:1621-1639.
6. Nejatidanesh F, Savabi O, Ebrahimi M, Savabi G. Retentiveness of implant-supported metal copings using different luting agents. *Dent Res J* 2012;9:13-18.
7. Guilherme N, et al. Effect of surface treatments on titanium alloy bonding to lithium disilicate glass-ceramics. *J Prosthet Dent* 2016;116:797-802.
8. Özcan M, Matinlinna J. Surface conditioning protocol for the adhesion of resin-based cements to base and noble alloys: how to condition and why? *J Adhes Dent* 2015;17:372-373.
9. Özcan M, Volpato CA. Adhesion Protocol for Bonding Abutments or Fixed Dental Prostheses on Titanium Bases in Implant-borne Reconstructions: How and Why? *The journal of adhesive dentistry*, 2016, Vol.18(3), pp.268-9.
10. Özcan M, Pekkan G, Khan A. Does rinsing following particle deposition methods have a negative effect on adhesion to titanium? *J Adhes Dent* 2013;15:307-310.
11. Özcan M, Valandro LF. Effect of silano coupling agents and alloy primers on adhesion to titanium. *Minerva Stomatol* 2011;60:427-434.
12. Özcan M, Volpato CA. Surface conditioning protocol for the adhesion of resin-based materials to glassy matrix ceramics: how to condition and why? *J Adhes Dent* 2015;17:292-293.
13. Özcan M. Air abrasion of zirconia resin-bonded fixed dental prosthesis prior to adhesive cementation: Why and how? *J Adhes Dent* 2013;15:394.
14. Rasetto FH, Driscoll CF, Prestipino V, Masri R, von Fraunhofer JA. Light transmission through all-ceramic dental materials: a pilot study. *J Prosthet Dent* 2004;91:441-446.
15. Rocha SS, Adabo GL, Spinola SG, Fonseca RG, Ferreira AR. Effect of metal conditioners on the adhesive bonding of resin cements to cast titanium. *Quintessence Int* 2007;38:506-510.
16. Taira Y, Yoshida K, Matsumura H, Atsuta M. Phosphate and triphosphate primers for bonding prosthodontic luting materials to titanium. *J Prosthet Dent* 1998;79:384-388.
17. Yanagida H, Matsumura H, Astuta M. Bonding of prosthetic composite material to Ti-6Al-7Nb alloy with eight metal conditioners and surface modification technique. *Am J Dent* 2001;14:291-294.
18. Adhesion Protocol for Bonding Abutments or Fixed Dental Prostheses on Titanium Bases in Implant-borne Reconstructions: How and Why? Mutlu Özcan / Cláudia Ângela Maziero Volpato *J Adhes Dent* 2016;18:268-269.

MDP (metacrilóiloxidecil dihidrogenofosfato)