

MedPark

S1

NEW BONE  
NEW LIFE

**DENTAL**  
Clinical case



# CONTENTS



<b>CLINICAL CASE (+BIOPSY)   1</b>	<b>Alveolar Ridge Preservation without Membrane</b>	p. 3
<b>CLINICAL CASE (+BIOPSY)   2</b>	<b>Alveolar Ridge Preservation</b>	p. 5
<b>CLINICAL CASE (+BIOPSY)   3</b>	<b>Alveolar Ridge Preservation without Membrane</b>	p. 8
<b>CLINICAL CASE (+BIOPSY)   4</b>	<b>Alveolar Ridge Preservation</b>	p. 11
<b>CLINICAL CASE (+BIOPSY)   5</b>	<b>Sinus Floor Elevation (Lateral Approach)</b>	p. 13
<b>CLINICAL CASE (+BIOPSY)   6</b>	<b>Vertical &amp; Horizontal Bone Augmentation</b>	p. 15
<b>CLINICAL CASE (+BIOPSY)   7</b>	<b>Regeneration of the Odontogenic Cyst Area</b>	p. 17
CLINICAL CASE   8	Horizontal Alveolar Bone Augmentation	p. 20
CLINICAL CASE   9	Minor Bone Augmentation without Membrane	p. 22
CLINICAL CASE   10	Moldable Augmentation in Anterior Area	p. 23
CLINICAL CASE   11	Alveolar Ridge Preservation	p. 25
CLINICAL CASE   12	Periodontal Defect Management	p. 27
CLINICAL CASE   13	Socket Management in Posterior Area	p. 28
CLINICAL CASE   14	Dehiscence Defect in Anterior Area	p. 29
CLINICAL CASE   15	Open Socket Management in Posterior Area	p. 30
CLINICAL CASE   16	Narrow Ridge Augmentation in Anterior Area	p. 32
CLINICAL CASE   17	Easy Augmentation of Dehiscence Defect	p. 33
CLINICAL CASE   18	Simplified Vertical Augmentation for Advanced Bone Resorption	p. 34
CLINICAL CASE   19	Dough Staged Sticky Bone Graft	p. 35
CLINICAL CASE   20	Immediate Implant Placement	p. 36
CLINICAL CASE   21	Easy 3D Ridge Augmentation	p. 38
CLINICAL CASE   22	Labial Fenestration Defect	p. 39
CLINICAL CASE   23	Minimally Invasive Bone Grafting in Anterior Area	p. 40
CLINICAL CASE   24	Horizontal Bone Augmentation	p. 41
CLINICAL CASE   25	Labial Dehiscence in Upper Anterior Area	p. 42
CLINICAL CASE   26	Labial Fenestration in Lower Anterior Area	p. 43
CLINICAL CASE   27	Peri-implantitis Treatment with Laser Therapy	p. 44
CLINICAL CASE   28	Extraction Socket Management	p. 45
CLINICAL CASE   29	Moldable Augmentaion ☒in Posterior Area	p. 46
CLINICAL CASE   30	Immediate Implant Placement	p. 49
CLINICAL CASE   31	Immediate Implant Placement	p. 52
CLINICAL CASE   32	Sinus Augmentation	p. 55
CLINICAL CASE   33	Socket Management in Posterior ☒Area & Immediate Implant Placement	p. 57
CLINICAL CASE   34	Minor Bone Augmentation☒(Maxillary Sinus Lift Revision Surgery)	p. 59
CLINICAL CASE   35	Horizontal Alveolar Bone ☒Augmentation	p. 61
CLINICAL CASE   36	Horizontal Alveolar Bone ☒Augmentation	p. 62
CLINICAL CASE   37	Horizontal Ridge Augmentation	p. 64
CLINICAL CASE   38	Horizontal Guided Bone Regeneration	p. 65
CLINICAL CASE   39	Narrow Ridge Augmentation ☒in Posterior Area	p. 66

# Alveolar Ridge Preservation without Membrane


**Initial assesment** A 78- year-old female patient has suffered from gingival swelling & pain due to a vertical root fracture of Lt. 2nd premolar.  
**Products** **S1 Bone graft material** (Powder type, 0.2~1.0mm) | **Use of membrane** ☐ Yes ☒ No

## Objectives

- › Extensive bone destruction was observed, including the apical and buccal aspects, as seen in the CT scan. Therefore, an alveolar ridge preservation procedure with S1 was required before placing an implant.

## Conclusions

- › Despite using S1 alone without membrane coverage for alveolar ridge preservation (ARP), extra-fine soft tissue healing and full coverage were observed within one month. The ridge contour was preserved for 5 months of healing time. Histologic evaluation revealed extra-fine new bone formation in both the quantity and quality of the regenerated bone.



Preoperative X-ray



Tooth extraction



Application of S1 bone graft material and suture (Open socket)



Post OP (2 weeks)



Post OP (2 months)




Post OP (5 months)



Post OP (5 months), Implant placement




Post OP (5 months), Suture




Post OP (5 months), Implant placement



Post OP (7 months), Custom abutment



Post OP (7 months), Final prosthesis




Alveolar bone changes on CBCT

**Biopsy time** 5 months after using S1 bone graft material

**Biopsy method** Collected at a depth of 6 mm using a trephine bur from the implant placement area of #13 (universal system)

**Findings** Despite the single use of S1 without membrane coverage for ARP, histologic evaluation revealed extra-fine new bone formation in both the quantity and quality of the regenerated bone.

› Hematoxylin & Eosin staining / Masson's Trichrome staining



Total Region of Interest Area		
	Area (mm <sup>2</sup> )	Percent(%)
Bone graft area	0.65	8.96
<b>New bone area</b>	<b>3.43</b>	<b>47.24</b>
Others	43.80	43.80
Total bone area	4.08	56.20



# Alveolar Ridge Preservation

**Initial assesment** Advanced peri-implantitis was found on the implant in the left maxilla, with a buccal bone defect, in a 41-year-old male patient who is suspected to be a heavy smoker.



**Products** **S1 Bone graft material** (Powder type, 0.2~1.0mm) | **Use of membrane** ☒ Yes ☐ No

## Objectives


- › After fixture removal, severe buccal dehiscence and poor soft tissue conditions allowed for only alveolar ridge preservation in the left premolar area. S1 bone graft was applied to the dehiscence defect, and a collagen membrane was used for coverage.

## Conclusions

- › Not only was a volumetric change found in the CT scan, but the clinical situation for implant placement was also acceptable. An additional bone graft was performed simultaneously with implant placement five months after ARP. Fine bone regeneration was observed in the histologic evaluation.



# Treatment Steps




Preoperative x-ray



Peri-implantitis condition




Defect after removal of the old fixture



Significant defect after removal of the old fixture



Application of S1 bone graft material



Application of S1 bone graft material




Bone covered with granulation tissue and a collagen membrane



Post-OP (5 months)



Post-OP (5 months)




Post OP (5 months), Implant placement



Post OP (5 months), Implant placement



Post OP (9 months), Prosthesis



Alveolar bone changes on CBCT




**Biopsy time** 5 months after using S1 bone graft material

**Biopsy method** Collected at a depth of 6 mm using a trephine bur from the implant placement area of #12 (universal system)

**Findings** Although it was a large defect, fine bone regeneration was observed in the histologic evaluation

› Hematoxylin & Eosin staining



Total Region of Interest Area		
	Area (mm²)	Percent(%)
Bone graft area	1.47	26.77
<b>New bone area</b>	<b>1.46</b>	<b>26.64</b>
Others	2.55	46.59
Total bone area	2.93	53.41

# Alveolar Ridge Preservation without Membrane

**Initial assesment** A 47-year-old male patient has a floating tooth on the right mandibular canine with advanced apical periodontitis. Alveolar ridge preservation was required due to a two-wall extensive bony defect and poor soft tissue condition. A healing period of more than six months was expected before implant placement.


**Products** **S1 Bone graft material** (Powder type, 0.2~1.0mm) | **Use of membrane** ☐ Yes ☒ No

## Objectives

› After tooth extraction, the thick inflamed soft tissue was carefully peeled off from the underlying bone. S1 bone was then applied to the large defect and covered with granulation tissue, which served as a protective membrane for the graft material.


## Conclusions

› For five months, the CT scans showed that the soft tissue and bony contours were well-preserved. However, despite histologic results indicating acceptable bone regeneration, the regenerated bone exhibited a soft quality. For implant placement in cases of 1-wall or 2-wall defects, a healing period of more than six months would be required.






# Treatment Steps




Preoperative X-ray



Huge defect after tooth extraction



Application of S1 bone graft material



Bone covered with granulation tissue without a membrane



Post-OP (2 weeks)



Post-OP (6 weeks)



Post-OP (5 months)




Post-OP (5 months)




Post OP (5 months), Implant placement



Post OP (5 months), Implant placement



Post OP(5months)



Alveolar bone changes on CBCT



Post OP (8 months), Custom abutment



Post OP (8 months), Final prosthesis




Post OP (9 months)

**Biopsy time** 5 months after using S1 bone graft material

**Biopsy method** Collected at a depth of 6 mm using a trephine bur from the implant placement area of #27 (universal system)

**Findings** In cases of 1-wall or 2-wall defects, a healing period of more than six months is required for implant placement. While the regenerated bone shows a soft quality, the histologic results are acceptable for implantation.

› Hematoxylin & Eosin staining



Total Region of Interest Area		
	Area (mm <sup>2</sup> )	Percent(%)
Bone graft area	1.21	16.00
<b>New bone area</b>	<b>1.52</b>	<b>20.16</b>
Others	4.835	63.83
Total bone area	2.73	36.17



# Alveolar Ridge Preservation

**Initial assesment** A patient came in for treatment of peri-implantitis affecting the mandibular first and second molars. The existing implant was removed, and inflammation was managed through curettage and guided bone regeneration surgery. Once the healing process is assessed, a new implant will be placed.


**Products** **S1 bone graft material** (Powder type, 0.2~1.0mm) | **Use of membrane** ☐ Yes ☒ No

## Objectives


- › When there is insufficient support space for implant placement due to the creation of a wide socket range after removing peripheral inflammatory tissue caused by peri-implantitis.
- › On the day of GBR surgery, after soft tissue elevation, traces of the bone graft material used for the previous implant placement were observed.

## Conclusions

- › Six months after using the S1 bone graft material, the boundary between the S1 graft and the natural alveolar bone naturally connected, generating new bone.
- › Biopsy results of the new bone formed at the implant placement site confirmed that osteocytes had developed around the graft material and connected to form mature bone.



Preoperative X-ray (#36, #37)



Implant removal and inflammation curettage



Application of S1 bone graft material



Application of S1 bone graft material



Intraoral photo taken 6 months after the application of S1 bone graft material



Incision of soft tissue for implant placement



Implant placement in newly formed bone



Suturing after attaching the healing abutment




X-ray taken 6 months after the application of S1 bone graft material and implant placement

**Biopsy time** 6 months after using S1 bone graft material


**Biopsy method** Collected at a depth of 6 mm using a trephine bur from the implant placement area of #36

- Findings**
- › After 6 months, the bone formation rate was confirmed to be 60.29%, with most of it forming a network with mature bone.
  - › The graft material is confirmed to be safe, with no signs of inflammation or immune rejection.
  - › Osteogenesis occurred around the graft material, with uniform distribution of osteocytes observed, confirming the mature stage of ossification.

› Hematoxylin & Eosin staining



› Masson's Trichrome staining



Total Region of Interest Area		
	Area (mm <sup>2</sup> )	Percent(%)
Bone graft area	1.44	20.99
New bone area	4.13	60.29
Others	1.28	18.72
Total bone area	5.57	81.28

※ Biopsy analysis : OBen tissue analysis



# Sinus Floor Elevation (Lateral Approach)

**Initial assesment** A patient who had been using dentures for a long time visited for implant placement. Vertical and horizontal bone augmentation was performed at the base of the maxillary sinus using a lateral approach, as there was only 1 mm of residual bone in the maxillary left posterior region near the sinus.

**Products** **S1 bone graft material** (Powder type, 0.2~1.0mm) | **Use of membrane** ☐ Yes ☒ No

## Objectives

- › Insufficient condition with only 1mm of vertical residual bone in the upper left molar area
- › Access to the maxillary sinus by creating a window through a lateral approach
- › Implanting S1 bone graft material after sinus elevation

## Conclusions

- › This case demonstrates the effectiveness of using S1 bone graft material to create a stable space by lifting the sinus floor membrane.
- › Hard bone formation was observed, with noticeable resistance when using Trephine burs to collect a biopsy from the S1 bone graft area.
- › Biopsy results indicated that new bone had formed well around the graft material, with mature bone tissue evident.



Preoperative X-ray



Intraoral photo taken before the grafting procedure in the left maxilla.




Application of S1 bone graft material from the buccal side of the remaining bone in left maxilla




X-ray taken immediately after applying the S1 bone graft material




Implant placement in the #26 area four months after using S1 bone graft material



Biopsy request for posterior region of #26 at 11 months after using S1



Final prosthesis



X-ray taken 26 months later




After 33 months of using S1 bone graft material

**Biopsy time** 11 months of using S1 bone graft material

**Biopsy method** Collected a sample from the posterior side of the implant placement area at a depth of 6 mm using a trephine bur

- Findings**
- › The graft material is considered to be safe, with no signs of inflammation or immune rejection
  - › High-quality bone has formed uniformly from the host bone to the maxillary sinus mucosa
  - › Most of the new bone is integrated with the graft material, creating a very dense network

› Hematoxylin & Eosin staining



› Comparison of S1 and Product 'B' for bone formation rate

	S1	Product 'B'
Graft area	Sinus	Sinus
Observation period	11 months	12 months
Bone formation rate	<b>28.22 %</b>	26.60 %

※ Biopsy analysis : Knotus, Tissue analysis by Wonkwang University

**Reference**

- (1) Son WK, Shin SY, Yang SM, Kye SB. Maxillary sinus floor augmentation with anorganic bovine bone: Histologic evaluation in humans. J Korean Acad Periodontol. 2009;39(1):95-102.
- (2) Lee YM, Shin SY, Kim JY et al. Bone reaction to bovine hydroxyapatite for maxillary sinus floor augmentation: Histologic results in humans. Int J Periodontics Restorative Dent. 2006;26:471-481.



# Vertical & Horizontal Bone Augmentation

**Initial assesment** The situation required extracting the entire maxillary tooth and removing inflamed tissue due to periodontitis. Because extensive horizontal and vertical bone reconstruction of the maxillary alveolar bone was needed, implant placement was planned for 4 months after the bone augmentation.

**Products** **S1 bone graft material** (Powder type, 0.2~1.0mm) | **Use of membrane** ☒ Yes (Non-absorbable) ☐ No

## Objectives

- › Alveolar bone loss due to periodontitis had progressed extensively, resulting in large bone defects in the maxillary right canine and first premolar areas.
- › Horizontal and vertical guided bone augmentation was performed using S1 bone graft material and a non-absorbable membrane.

## Conclusions


- › Using a widely performed guided bone regeneration technique, the biopsy revealed excellent ossification and bone quality, with no visible boundary between the bone graft and the new bone.
- › After elevating the soft tissue, it was found that 4 months after using the bone graft material, there was sufficient bone quality and width for implant placement.



Preoperative X-ray



Postoperative X-ray taken after using S1 bone graft material



CBCT taken after using S1 bone graft material



Intraoral photo of the lateral right maxilla taken 3 months later




Intraoral photo taken 3 months later



Soft tissue incision for implant placement after 4 months



Abutment connection for the final prosthetics



X-ray taken after the final prosthetics




Final prosthetics

**Biopsy time** 4 months of using S1 bone graft material

**Biopsy method** Collected at a depth of 6 mm using a trephine bur from the implant placement area of #36

- Findings**
- › Ossification was outstanding, with no visible boundary between the bone graft material and the new bone, and the bone quality was high.
  - › Mature ossification has been confirmed, with osteoblasts observed in the bone graft material, indicating a stage beyond the initial bone formation phase where osteoblasts are actively involved.
  - › The graft material is considered to be safe, showing no signs of inflammation or immune rejection.
  - › Bone formation is very successful, with a dense connection between the graft material and the new bone, reflecting a strong affinity between them.

› Hematoxylin & Eosin staining



› Comparison of S1 and Product ‘B’ for bone formation rate

	S1	Product ‘B’
Graft area	Alveolar ridge	Alveolar ridge
Observation period	4 months	6 months
Bone formation rate	25.23 %	18.3 %

※ Biopsy analysis : OBen tissue analysis

**Reference**

(1) Lee YM, Shin SY, Kim JY et al. Bone reaction to bovine hydroxyapatite for maxillary sinus floor augmentation: Histologic results in humans. Int J Periodontics Restorative Dent. 2006;26:471-481.

(2) Nicola U. Zitzmann, Dr Med Dent\* Peter Schärer, Prof Dr Med Dent, MS\*\* Carlo P. Marinello, Prof Dr Med Dent, MS\*\*\* Peter Schüpbach, Dr Sc Nat, PhD\*\*\*\* Tord Berglundh, DMD, PhD\*\*\*\*\*

(3) Zitzmann NU, Schärer P, Marinello CP, Schüpbach P, Berglundh T. Alveolar ridge augmentation with Bio-Oss: a histologic study in humans. Int J Periodontics Restorative Dent. 2001 Jun;21(3):288-95. PMID: 11490406.



# Regeneration of the Odontogenic Cyst Area

**Initial assesment** To reconstruct the bone defect caused by an inflammatory cyst at the root of the mandibular right first molar, guided bone reconstruction surgery using S1 bone graft material will be performed first. Implant placement will follow after confirming bone tissue formation.



**Products** **S1 bone graft material** (Powder type, 0.2~1.0mm) | **Use of membrane** ☐ Yes ☒ No

## Objectives

- › Implant placement is planned for the right mandibular first and second molars
- › A large cyst and inflamed tissue were observed beneath the apex of the first molar
- › Implant placement will proceed after extraction and guided bone reconstruction in the affected area

## Conclusions


- › Three months after using S1 bone graft material, a dense new bone formation was observed in the radiographs
- › Biopsy results from the implant placement site confirmed successful new bone induction
- › Healthy new bone, including blood vessels, has formed even in areas with large defects



# Treatment Steps



Preoperative X-ray (before the extraction of #46)



X-ray taken after the extraction of the mandibular right first molar



Intraoral photo taken after the extraction of the mandibular right first molar



Residual bone condition after elevating the soft tissue flap




Application of S1 bone graft material (#46)




Application of S1 bone graft material (#47)



X-ray taken 3 months after the application of S1 bone graft material



Postoperative CBCT cross-section




Implant placement 4 months after the use of S1 bone graft material



Postoperative intraoral photo after the final prosthetics



Final prosthetics (#46, #47)



X-ray after the final prosthetics




**Biopsy time** 4 months after using S1 bone graft material


**Biopsy method** Collected at a depth of 6mm using a trephine bur from the implant placement area of #36

- Findings**
- › 22.3% of the new bone was preserved in the large defect socket area, which is a significant finding
  - › A substantial amount of new bone formed between the safe graft material and the bone graft material, with no notable inflammatory reaction
  - › Osseointegration anticipated to be excellent
  - › Forming a pattern indicating very healthy new bone with blood vessels

› **Hematoxylin & Eosin staining**



› **Masson's Trichrome staining**



Total Region of Interest Area		
	Area (mm <sup>2</sup> )	Percent(%)
Bone graft area	2.10	11.46
<b>New bone area</b>	<b>4.09</b>	<b>22.34</b>
Others	12.13	66.18
Total bone area	6.19	33.81

※ Biopsy analysis : OBen tissue analysis

# Horizontal Alveolar Bone Augmentation

Dr. Dong-Wook Chang

**Initial assesment** A 63-year-old woman visited the dental clinic because of the discomfort with her dentures and expressed a desire for implants.


**Products** **S1 Bone graft material** (Powder type, 0.2~1.0mm) | **Use of membrane** ☒ Yes ☐ No

## Objectives

- › Extraction of the remaining canine teeth and implantation of implants in both the anterior and posterior areas
- › Bone grafting was performed due to resorption on the labial side of the anterior teeth

## Conclusions

- › Grafting S1 into the thin labial of the anterior teeth, followed by examining the augmented alveolar bone and healed soft tissue after 3 months using
- › CBCT





# Treatment Steps




Preoperative panoramic radiograph




Preoperative intraoral photo



Preoperative CBCT



Extraction of the remaining upper right canine tooth and partial incision of the anterior tooth




Transplantation of S1 bone graft material



Use of COLLA Membrane



Suture



Panoramic radiograph taken 3 months after surgery




Intraoral photograph showing soft tissue healing 3 months after surgery



Final panoramic radiograph



Final prosthesis



Postoperative CBCT

# Minor Bone Augmentation without Membrane

**Initial assesment** The edentulous right maxillary premolar area of a 67-year-old female patient did not have sufficient horizontal volume for implant placement. An onlay bone graft was required, but the patient opted against extensive GBR due to diabetes.

**Products** **S1 Bone graft material** (Powder type, 0.2~1.0mm) | **Use of membrane** ☐ Yes ☒ No

## Objectives

› After implant placement, a buccal fenestration defect was discovered on the first premolar implant. S1 bone was applied to the buccal defect without the use of a membrane or additional fixation for space maintenance. Minor flap advancement was performed, and primary closure was achieved using only 5-0 nylon sutures.

## Conclusions

› Abundant bone regeneration was observed after three months, with the regenerated bone even covering the cover screws. The volume of regenerated buccal bone remained stable over 8 months.



Preoperative X-ray



Implant placement




Buccal fenestration defect



Suture



Post OP (3 weeks)



Post OP (3 months), 2nd OP



Post OP (3 months), 2nd OP




Post OP (3 months), 2nd OP




Post OP (4 months)



Post OP( 4 months), Final prosthesis



Post OP (4 months), Final prosthesis



Alveolar bone changes on CBCT



# Moldable Augmentation in Anterior Area

**Initial assesment** A 66-year-old male patient presented with complaints of a strong odor around his maxillary anterior prosthesis (a 4-unit bridge) and expressed a desire for implant treatment. Severe alveolar bone contraction had progressed due to root inflammation, with a residual buccal bone width of 3 mm observed on CBCT.



**Products** **S1 bone graft material** (Powder type, 0.2~1.0mm) | **Use of membrane** ☐ Yes ☒ No

## Objectives

- › Extensive alveolar bone augmentation is needed due to the buccal width of the alveolar bone being only 3mm, necessitating the extraction of the remaining maxillary teeth (#11, #21, #22) and the removal of inflamed granulation tissue.
- › Upon visual inspection, a defect in the form of fenestration and dehiscence of the alveolar bone was observed after soft tissue elevation.

## Conclusions

- › Good healing of the soft tissue around the bone graft material was observed, despite the wide surgical area.
- › CBCT revealed that the new bone around the implant is naturally integrated with the existing alveolar bone.
- › Although no membrane was used, the position of the bone graft material remained stable, and the bone augmentation surgery was performed successfully and with relative ease.






# Treatment Steps



Preoperative panoramic radiograph




Preoperative CBCT



Preoperative intraoral photo



Implant placement at positions #11, #21, and #22



Dehiscence defect and significant bone loss



Application of S1 bone graft material



Application of S1 bone graft material




Suture



6 weeks after grafting S1



X-ray taken after final prosthetics



CBCT taken after final prosthetics



Final prosthetics

# Alveolar Ridge Preservation

**Initial assesment** Visited for implant placement in the tooth loss area, with a plan for vertical bone augmentation using bone graft material simultaneously with the implant placement in the defect of the mandibular left first molar

**Products** **S1 bone graft material** (Powder type, 0.2~1.0mm) | **Use of membrane** ☐ Yes ☒ No


## Objectives

- › The defect area at #36 displays both vertical and horizontal alveolar resorption.
- › The plan is to perform vertical bone augmentation simultaneously with implant placement to address the discrepancy between the adjacent teeth and the excessive bone loss.


## Conclusions

- › The vertically augmented alveolar bone around the #36 implant is identified by comparing preoperative and postoperative X-rays taken three months after placing the final prosthesis.
- › Because of the successful bone and soft tissue augmentation, a final prosthesis was fabricated that aligns well with the adjacent teeth, which is expected to benefit future periodontal management.

Before



After



Before



After



Before




After





# Treatment Steps



Preoperative X-ray



Preoperative intraoral photo



Implant placement following soft tissue flap elevation and inflammation curettage



Application of S1 bone graft material to the implant bone loss area



Vertical and horizontal bone augmentation using S1



Postoperative X-ray




10 weeks after the surgery



Observation of new bone formation during the second surgery, 10 weeks after the initial procedure




Abutment connection for tooth restoration



Final prosthetics 3 months after S1 grafting



Final prosthetics 3 months after S1 grafting



Alveolar bone changes on CBCT  
Before surgery    Right after surgery    4months



# Periodontal Defect Management

Dr. Chingu kim

**Initial assesment**    Advanced periodontal disease was observed in the left mandibular molar, along with a buccal bone defect, in a 57-year-old patient


**Products**                **S1 bone graft material** (Powder type, 0.2~1.0mm)    |    **Use of membrane** ☐ Yes ☒ No

**Objectives**

› After removing the inflammatory granulation tissue, destruction of the alveolar bone connected to the buccal side of the two molars was observed, along with a large, puddle-shaped defect around the mesial root of the second molar. It is planned to use only S1 bone graft material in the defect area, where membrane placement is challenging.


**Conclusions**

› Although no membrane was used, the volume of the bone graft material was well-maintained, and new bone formation was evident on the CT scan.  
› New bone formation was observed as radiopacity on the CBCT, and after two months, it had integrated with the surrounding bone tissue.




1

Preoperative X-ray




2

Flap operation




3

Inflammation curettage




4

Defect following the removal of inflamed tissue




5

Ready-to-use S1 bone graft material




6

Application of S1 bone graft material




7

Application of S1 bone graft material




8

Suture




9

Post OP




10

Post OP (3 months)



11

Post OP (3 months)



12

Alveolar bone changes on CBCT

# Socket Management in Posterior Area

**Initial assesment** A patient with severe periodontitis in the maxillary left first molar had minimal or no remaining bone at the site for implant placement. Bone augmentation was performed in the tooth extraction area.


**Products** **S1 bone graft material** (Powder type, 0.2~1.0mm) | **Use of membrane** ■ Yes (PRF) □ No

## Objectives

- › Alveolar bone condition after extraction of the maxillary first molar
- › Assessment of the extraction site and defect area following soft tissue elevation
- › Hydration of S1 with PRF followed by application of PRF membrane

## Conclusions


- › Despite the large defect area, the shape was well maintained without a membrane, and new bone formed effectively.
- › Excellent results were achieved using only residual bone and S1 graft material, without the need for maxillary sinus elevation.
- › Bone formation was successful, with the height of the marginal bone around the implant harmonizing well with the adjacent teeth.



Preoperative X-ray




Preoperative intraoral photo



Hydration of S1 using PRF



Application of S1 bone graft material after implant placement




X-ray 1 month after using S1 bone graft material




Bone formation observed during the second surgery, 7 months after placement



Abutment placement before tooth restoration



X-ray taken after final prosthetics



X-ray taken 10 months after using S1 bone graft material



After several months of using S1 bone graft material



# Dehiscence Defect in Anterior Area

**Initial assesment** Horizontal bone augmentation surgery and implant placement were planned after extracting the maxillary left and right lateral incisors, with a remaining net lingual width of 2.8mm

**Products** **S1 bone graft material** (Powder type, 0.2~1.0mm) | **Use of membrane** ☐ Yes ☒ No

## Objectives

- › Extraction of the maxillary left and right lateral incisors due to apical inflammation
- › Horizontal bone augmentation was performed, with the remaining bone having a lingual width of 2.8 mm
- › S1 graft material was applied after placing implants with a diameter of 3.5mm

## Conclusions


- › Despite the insufficient bone width and extensive range of the defect area, the desired shape was achieved using S1 bone graft material, and the graft material remained stably positioned.
- › The S1 bone graft material maintained its position without the use of a membrane, promoting new bone regeneration.
- › The periodontal tissue in the S1 graft area is aesthetically pleasing and healed stably.



Preoperative X-ray



Preoperative intraoral photo (temporary crown on #12~22)




X-ray taken after the extraction of the maxillary left and right lateral incisors




Bone loss at the site of the maxillary right lateral incisor



Molding S1 bone graft material




X-ray taken after implant placement and the application of S1 bone graft material



Bone formation pattern 4 months after using S1 bone graft material



Final prosthetics



X-ray taken 11 months after using S1 bone graft material



# Open Socket Management in Posterior Area

**Initial assesment** Periodontitis in the maxillary right first and second premolar and first molar areas has progressed over a long period with inflammation, affecting the surrounding alveolar bone and significant bone loss observed on CBCT. The patient is seeking implant treatment.


**Products** **S1 bone graft material** (Powder type, 0.2~1.0mm) | **Use of membrane** ☐ Yes ☒ No

## Objectives


- › Confirmation of alveolar bone loss and remaining bone morphology due to periodontitis in the maxillary right first and second premolar and first molar areas
- › Immediate implant placement planned following the extraction of #14 and #16, with a 3-unit prosthesis scheduled for production
- › S1 bone graft material will be used to address the bone defect at the extraction sites

## Conclusions

- › Despite being sutured in an open wound state immediately after implant placement with bone graft, secondary healing occurred without significant loss of the S1 graft material.
- › To account for soft tissue contraction and the potential separation of graft particles when implanted immediately after extraction, a transmucosal bone graft was used to extend the graft material up to the top of the soft tissue.
- › Well-formed new bone around the implant was observed on the X-ray 3 months after surgery.



# Treatment Steps



Inflammation and alveolar bone morphology on preoperative CBCT



Preoperative intraoral photo



Identification of bone loss area following the extraction of #14, #15, and #16



Hydration of S1 bone graft material



Implant placement following S1 grafting in the extensive bone loss area



Implant placement following S1 grafting in the extensive bone loss area



Suture



Intraoral photo taken 2 months after surgery



Final prosthetics



X-ray taken 3 months after using S1 bone graft material



CBCT taken 3 months after using S1 bone graft material

# Narrow Ridge Augmentation in Anterior Area

Dr. Choong Noh

**Initial assesment** The patient visited with peri-implantitis diagnosed at another clinic. After removing the existing implant, S1 bone graft material was applied to address the bone defect.

**Products** **S1 bone graft material** (Powder type, 0.2~1.0mm) | **Use of membrane** ☐ Yes ☒ No

## Objectives

- › Implant removal due to peri-implantitis in the mandibular anterior area (#31)
- › Horizontal bone augmentation surgery was performed on the labial side of the residual bone

## Conclusions

- › Final prosthetic placement was completed 3 months after S1 implantation.
- › Aesthetic periodontal tissue was observed after the second implant placement.



Preoperative X-ray



Preoperative intraoral photo



Application of S1 bone graft material for horizontal bone augmentation following implant placement (single-stage surgery)



Suture



Abutment placement for prosthetics 3 months after S1 grafting



Abutment placement before prosthetics



Condition of the gingival tissue



Final prosthetics



X-ray taken after final prosthetics



# Easy Augmentation of Dehiscence Defect

**Initial assesment** A 67-year-old woman visited the hospital for the reconstruction of a defect in the maxillary left lateral incisor

**Products** **S1 bone graft material** (Powder type, 0.2~1.0mm) | **Use of membrane** ☒ Yes (CGF) ☐ No

## Objectives

- › Labial resorption of the surrounding alveolar bone was noted after the extraction of the maxillary left lateral incisor.
- › To address aesthetic concerns, horizontal bone augmentation was planned. Bone graft material was required to increase the height of the alveolar bone in the marginal area. With the patient's consent, blood was collected to process CGF, which would be used as a membrane.

## Conclusions

- › Despite the insufficient bone width in the defect area and the challenging anatomical structure, the moldable characteristics of S1 allowed us to easily achieve the desired shape.
- › The position of the bone graft material was stably maintained without using a membrane, which helped preserve the shape.



X-ray of the defect area of alveolar bone before surgery



Incision and flap elevation



Dehiscence defect condition after implant placement



Application of S1 bone graft material



Molding to fit the shape of the surrounding alveolar bone



Formation of a membrane using CGF



Postoperative X-ray

# Simplified Vertical Augmentation for Advanced Bone Resorption

**Initial assesment** The patient, who had the second premolar and the first and second molars extracted from the mandibular right posterior region due to severe periodontitis, is planning to undergo implant placement along with vertical bone augmentation of the defect area.


**Products** **S1 bone graft material** (Powder type, 0.2~1.0mm) | **Use of membrane** ☒ Yes (COLLA) ☐ No

## Objectives

- › After removing the teeth and inflammatory granulation tissue due to severe periodontitis, the residual bone has an irregular shape, necessitating vertical bone augmentation
- › New bone regeneration will be induced using S1 bone graft material following implant placement in the #45, #46, and #47 regions

## Conclusions

- › Despite the extensive bone loss in the defect area, S1 bone graft material could be applied stably to the desired region.
- › The S1 bone graft material retained its position without the use of a non-absorbable shape-fixing membrane, facilitating new bone regeneration.
- › The periodontal tissue in the S1 graft area healed both aesthetically and stably.



X-ray of the defect area of alveolar bone before surgery




Alveolar bone defect after implant placement




Application of S1 bone graft material to the bone defect area



Use of a COLLA membrane for vertical bone augmentation



Postoperative X-ray



X-ray taken 4 months after bone graft material transplantation



# Dough Staged Sticky Bone Graft

**Initial assesment** A 54-year-old woman visited for the reconstruction of defects in the mandibular right second premolar and first molar areas

**Products** **S1 bone graft material** (Powder type, 0.2~1.0mm) | **Use of membrane** ■ Yes (COLLA) □ No

## Objectives


- › Since the height of the alveolar bone in the mandibular left defect area is uneven, the shape of the marginal area needs to be reconstructed.
- › Intramarrow penetration is conducted to encourage the formation of new blood vessels for the bone graft material, as considerable time has passed since the extraction. Meanwhile, the remaining bone surface is healed with dense bone.

## Conclusions

- › Radiographs taken 3 months after surgery showed that the shape of the marginal bone integrated harmoniously with the adjacent teeth, and the new bone density was satisfactory.
- › The moldable property of S1 was confirmed as advantageous in maintaining the one-wall defect at the distal site of the #46 implant placement.



Preoperative X-ray



X-ray of the defect area of alveolar bone before surgery




Incision and flap elevation



Alveolar bone defect after implant placement



Application of S1 bone graft material to the bone defect area



Formation of a membrane using COLLA



Final coverage using CGF



Suture



X-ray taken 3 months after surgery

# Immediate Implant Placement

**Initial assesment** Visited because the maxillary left central incisor was loose

**Products** **S1 bone graft material** (Powder type, 0.2~1.0mm) | **Use of membrane** ☒ Yes (COLLA) ☐ No

## Objectives

- › Immediate implant placement was performed to preserve the alveolar bone and soft tissue.
- › Due to labial bone loss, dehiscence defect surgery was carried out, with the incision kept minimal to account for potential edema.


## Conclusions

- › The use of S1 and COLLA resulted in better preservation of the periodontal tissue compared to before extraction, and aesthetically satisfactory outcomes were achieved.
- › The procedure involved minimizing the incision on the periodontal tissue, and the membrane was securely fixed, leading to a favorable prognosis with strong sutures.

Before



After




Before



After



Before



After






# Treatment Steps



Preoperative X-ray



Preoperative intraoral photo



Incision and flap elevation



Application of S1 bone graft material after Implant placement



Application of S1 bone graft material



Use of COLLA membrane



Suture



Postoperative X-ray



Intraoral photo taken 2 months after surgery



Final prosthetic connection 3 months after surgery



Intraoral photo after prosthetics



X-ray 3 months after surgery

# Easy 3D Rigde Augmentation

**Initial assesment** A 75-year-old woman visited to restore the defect area in the mandibular right first premolar and molar


**Products** **S1 bone graft material** (Powder type, 0.2~1.0mm) | **Use of membrane** ☒ Yes (CGF) ☐ No

## Objectives

- › A socket-shaped, three-wall defect was identified in the #44 area after soft tissue elevation during surgery
- › Vertical bone augmentation was carried out in the #44 to #46 regions using S1 bone graft material

## Conclusions

- › Despite the narrow bone width and wide defect area, S1 bone graft material was applied stably to the desired region.
- › After the surgery, the shape of the marginal bone integrated harmoniously with the adjacent teeth, and the new bone density was satisfactory.



Alveolar bone shape in the defect areas #44, #45, #46, Implant placement and #47




Application of S1 bone graft material to the bone defect area



Application of S1 bone graft material to the bone defect area



Use of CGF for vertical bone augmentation surgery



Postoperative X-ray



# Labial Fenestration Defect

**Initial assesment** A 57-year-old female patient visited for implant placement of the mandibular right central incisor

**Products** **S1 bone graft material** (Powder type, 0.2~1.0mm) | **Use of membrane** ■ Yes (COLLA) □ No



Preoperative X-ray



Preoperative intraoral photo (#41)



Defect area after implant placement



Application of S1 bone graft material



Application of COLLA membrane



Suture



Postoperative X-ray



Final prosthetic connection 3 months after surgery



Postoperative intraoral photo after final prosthetics

# Minimally Invasive Bone Grafting in Anterior Area

**Initial assesment** A 67-year-old female patient using partial mandibular anterior dentures visited for implanting a fixed prosthesis. The plan included performing narrow alveolar bone augmentation while minimizing the incision area during surgery.

**Products** **S1 bone graft material** (Powder type, 0.2~1.0mm) | **Use of membrane** ■ Yes (COLLA) □ No



Preoperative X-ray



Preoperative intraoral photo



Application of a COLLA membrane after incision (#42)



Application of a COLLA membrane after incision (#42)



Application of S1 bone graft material between the membrane and the alveolar bone (#32)



Suture



Intraoral photo taken after prosthetics



Intraoral photo taken after prosthetics



X-ray taken after final prosthetics



# Horizontal Bone Augmentation

**Initial assesment** A 53-year-old female patient visited for the reconstruction of the mandibular second premolar and molar areas. A significant time has elapsed since the tooth extraction, leading to buccal bone resorption. The plan is to perform horizontal guided bone regeneration following implant placement.

**Products** **S1 bone graft material** (Powder type, 0.2~1.0mm) | **Use of membrane** ■ Yes (COLLA) □ No



Preoperative X-ray (#45, #46, and #47)



Preoperative intraoral photo



Implant placement



Application of S1 bone graft material and COLLA membrane



Suture



Intraoral photo after prosthetics



Final prosthetics

# Labial Dehiscence in Upper Anterior Area

**Initial assesment** Visited for the restoration of a lost maxillary left canine, with implant placement planned

**Products** **S1 bone graft material** (Powder type, 0.2~1.0mm) | **Use of membrane** ■ Yes (COLLA) □ No



Preoperative intraoral photo



Elevation after soft tissue incision



Application of 0.15 g of S1 bone graft material



Molding of S1 bone graft material to fit the defect area



Application of COLLA membrane



Suture



Intraoral photo after prosthetics



Final prosthetics



# Labial Fenestration in Lower Anterior Area

**Initial assesment** Visited for the restoration of missing mandibular left central incisor and right lateral incisor. Due to the loss of buccal width in the anterior residual bone, guided bone regeneration is planned to be performed simultaneously with implant placement.

**Products** **S1 bone graft material** (Powder type, 0.2~1.0mm) | **Use of membrane** ■ Yes (COLLA) □ No



Preoperative intraoral photo (#31, #42)




Elevation after soft tissue incision



Implant placement



Application of S1 bone graft material



Application of COLLA membrane



Suture

# Peri-Implantitis Treatment

## with Laser Therapy

Dr. Kacper Koryzna

**Initial assesment** A 62-year-old female patient, with inflammation next to the prosthesis at #37, where an implant was placed 10 years ago

**Products** **S1 bone graft material** (Powder type, 0.2~1.0mm) | **Use of membrane** ☐ Yes ☒ No

### Objectives

- › Assess the bone defect around the implant in #37 (the left mandibular second molar) using a panoramic X-ray
- › Treating the defect with laser therapy and filling the bone with S1

### Conclusions

- › Removing granuloma from the implant surface, and filling the bone defect with S1
- › S1 adhered well to the implant surface
- › New bone formation is expected to occur effectively between the implant threads



Preoperative X-ray (#37)



Laser-assisted incision



Cyst and inflamed tissue



Application of S1 bone graft material to the bone defect area



Postoperative X-ray taken after using S1 bone graft material



# Extraction Socket Management

Dr. Alesio Bocari

**Initial assesment** Toothache and swollen gum in the site #27 caused by vertical root fracture (VRF)

**Products** **S1 bone graft material** (Powder type, 0.2~1.0mm) | **Use of membrane** ☐ Yes ☒ No

## Objectives

› After extracting the #27 tooth due to a vertical root fracture (VRF) and removing the inflammation, an immediate implant is placed, and S1 is applied.

## Conclusions

› Alveolar Ridge Preservation (ARP) was performed on #27 using S1 bone graft material without a membrane, and the gingival soft tissue healed without inflammation.  
› After a year, the site remained very stable with no marginal bone resorption observed.



Preoperative (#27)



Measuring the socket size after tooth extraction



Immediate implant placement



Hydration of S1



Socket management by applying S1



Healing abutment placement



#27 Suture



After 1 week



After 1 month & soft tissue healing



After 1 year

# Moldable Augmentation

## in Posterior Area

**Initial assesment** Visited due to chronic inflammation and mobility of the mandibular left first molar and mandibular right first molar



**Products** **S1 bone graft material** (Powder type, 0.2~1.0mm) | **Use of membrane** ☐ Yes ☒ No

### Objectives

- › Extraction planned due to inflammation and vertical mobility of both first molars
- › Three-month observation period after removal of inflammatory tissue following extraction

### Conclusions

- › Both mandibular first molars were close to the nerve and had extensive bone defects, but stability in new bone regeneration and healing were observed.
- › Vertical bone augmentation was performed without a membrane, and CBCT results confirmed that the alveolar bone shape remained stable with fine new bone regeneration.






# Treatment Steps



3 months after the extraction of both mandibular first molars



Intraoral photo of the left mandibular first molar before surgery



Preoperative CBCT of the mandibular right and left first molars



Incision for the mandibular left first molar



Implantation of the mandibular left first molar




Transplantation of S1 bone graft material shaped to fit the defect of the mandibular left first molar



Suture



Postoperative (4 weeks)



Revision surgery (4 weeks later)



Mandibular left first molar, 10 weeks after surgery: healing abutment tightened and sutured



Final prosthesis for the mandibular left first molar



Intraoral photo of the mandibular right first molar before surgery

# Treatment Steps



Gum incision for the right mandibular first molar



Implant placement for the right mandibular first molar



Transplantation of S1 bone graft material, shaped to fit the defect in the right mandibular first molar




Suture



Post OP (2 months)



2nd surgery



CBCT taken after surgery on bilateral mandibular first molars



Right mandibular first molar prosthesis



X-ray taken after the final prosthetics



# Immediate Implant Placement

**Initial assesment** Visited for a fractured maxillary central incisor that had undergone root canal treatment



**Products** **S1 bone graft material** (Powder type, 0.2~1.0mm) | **Use of membrane** ☐ Yes ☒ No

## Objectives

- › Immediate implant placement following the extraction of the fractured maxillary central incisor
- › Bone grafting planned for the extraction socket and thin labial bone

## Conclusions

- › Bone grafting using the tunneling technique on the thin labial bone in the apical region of the maxillary central incisor
- › 3 months after surgery, the bone remained well-maintained and showed no resorption, even with the use of a membrane.




## Treatment Steps



1  
Preoperative panoramic radiograph



2  
Preoperative intraoral photo



3  
#22  
Preoperative CBCT



4  
After the extraction



5  
Grafting S1 bone material



6  
Implant placement in the optimal position



7  
Tunneling technique – soft tissue incision



8  
Tunneling technique – grafting the bone material onto the apical region



9  
Tunneling technique – Suture



10  
Postoperative panoramic radiograph



11  
1 month after surgery



12  
#22  
CBCT taken after surgery



## Treatment Steps



Soft tissue healing 2 months after surgery



X-ray taken 2 months after surgery




Soft tissue healing 3 months after surgery



Panoramic radiograph of the final prosthesis



Final prosthesis of the maxillary left lateral incisor



CBCT of the final prosthesis

# Immediate Implant Placement

**Initial assesment** Visited due to discomfort caused by inflammation and mobility of the bridge prosthesis on the maxillary central and lateral incisors



**Products** **S1 bone graft material** (Powder type, 0.2~1.0mm) | **Use of membrane** ☐ Yes ☒ No

## Objectives

- › Extraction of the maxillary anterior bridge teeth, and removal of the granulation tissue
- › Planning the implant placement and the horizontal bone augmentation simultaneously by observing the absorbed shape on the labial side

## Conclusions

- › Bone grafting and immediate implant placement following the removal of granulation tissue from the maxillary central and lateral incisors
- › 3 months after surgery, the gum and bone graft volume was well-maintained, even without the use of a membrane






# Treatment Steps



Preoperative panoramic radiograph



Preoperative intraoral photo



Preoperative CBCT



Extraction



Molded S1 bone graft material




Grafting S1 bone material



Implant placement in the optimal position




Suturing and bone grafting using the tunneling technique in the apical region of the upper right lateral incisor



Bone grafting using the tunneling technique at the site of the left maxillary central incisor



Postoperative panoramic radiograph



Postoperative CBCT



Temporary prosthesis placement

## Treatment Steps

---



Post OP (6 weeks)




Post OP (3 months)



Final prosthesis



Panoramic radiograph of the final prosthesis



CBCT of the final prosthesis



# Sinus Augmentation

**Initial assesment** Visited for implant placement in the maxillary posterior region

**Products** **S1 bone graft material** (Powder type, 0.2~1.0mm) | **Use of membrane** ☐ Yes ☒ No

## Objectives

- › Inflammation in the left maxillary first and second molars requires removal
- › Planning a maxillary sinus lift and immediate implant placement, with the sinus perforation area sealed using PRF

## Conclusions

- › Inflammation in the left maxillary molar was removed using a syringe
- › 3 months after surgery, the size of the inflammation had decreased, and new bone regeneration and implant stability were well maintained.
- › At one year and one year and six months post-surgery, no inflammatory changes were observed, and mature bone tissue was present



Preoperative panoramic radiograph



Preoperative intraoral photo



Preoperative CBCT



Removal of the inflammation using a syringe



Hydration and mixing of S1 bone graft powder



S1 bone graft material placed at the site of the maxillary left first and second molars



Implant placement in the maxillary left first and second molars



Fastening of the healing abutment




Postoperative panoramic radiograph

# Treatment Steps




Final prosthesis placed after 3 months




Postoperative CBCT




CBCT of the final prosthesis



CBCT of the final prosthesis



Post OP (1 year)



Post OP (1 year and 6 months)



# Socket Management in Posterior Area & Immediate Implant Placement

Dr. Young Jin Cho

**Initial assesment** An 80-year-old man visited for rapid reconstruction due to food impaction, mobility, and inflammation odor in the mandibular right first molar area




**Products** **S1 bone graft material** (Powder type, 0.2~1.0mm) | **Use of membrane** ☐ Yes ☒ No

## Objectives

- › Due to the vertical mobility of the mandibular first molar being at 3 degrees, an immediate implant placement and bone grafting are scheduled following the tooth extraction.

## Conclusions


- › After extracting the mandibular right first molar, inflammation was removed, revealing thin buccal-lingual cortical bone and a large defect.
- › Bone grafting was performed after implant placement to address the large defect at the extraction site.
- › 3 months after implant placement, the final prosthesis was completed, showing well-maintained shape, successful osseointegration, and stable healing of the gingival soft tissue.



# Treatment Steps



Preoperative panoramic radiograph



Preoperative X-ray



Preoperative intraoral photo



Extraction



Immediate implant placement



Placement of moldable S1 bone graft material




Insertion of the bone graft material



Suture




Post OP (4 weeks)



Post OP X-ray (4 weeks)



Post OP (3 months)



Post OP X-ray (3 months)



# Minor Bone Augmentation (Maxillary Sinus Lift Revision Surgery)

Dr. Young Jin Cho

**Initial assesment** A 57-year-old woman needed implant placement and removal of inflammatory tissue because bone regeneration was unsuccessful following right maxillary sinus lift surgery due to inflammation. Planning for implant placement will proceed after evaluating the results of the reoperation for maxillary sinus elevation.

**Products** **S1 bone graft material** (Powder type, 0.2~1.0mm) | **Use of membrane** ☒ Yes ☐ No

## Objectives


- › Removal of the implant and granulation tissue due to lack of bone regeneration on the right side of the maxilla and alveolar bone destruction from inflammation following maxillary sinus lift surgery
- › Planning for new bone regeneration and implant placement using S1 sticky bone graft material for the large defect in the maxillary sinus

## Conclusions

- › During the removal of granulation tissue due to a maxillary sinus infection, perforation occurred, which was closed with a membrane, and S1 bone grafting was performed for the large defect.
- › Despite the challenging case of maxillary sinus perforation and large defects, the use of sticky bone allowed easy shaping of the defect without causing particle separation and infection of the bone graft material.
- › 6 months after maxillary sinus elevation and implant placement in the large defect, a significant reduction in maxillary sinus inflammation was observed, with bone levels and implant fixation stably maintained.




# Treatment Steps



Previous) Maxillary sinus lift CBCT (#15, #16)



Implant removal



Panoramic radiograph after implant removal



Large defects observed after removal of granulation tissue from the maxillary sinus



Insertion and coverage of a membrane at the perforation site in the maxillary sinus



S1 bone graft material preparation




Molding of the bone graft material to fit the large defects



Maxillary sinus lift revision surgery with bone graft




Suture



CBCT taken 4 days after surgery, confirming limited dispersion of the bone graft material



Intraoral photo demonstrating bone regeneration via gum incision



Post OP (6 months)



Implant placement



Identification of gingival soft tissue healing



Final prosthesis



# Horizontal Alveolar Bone Augmentation

Dr. Dong-Wook Chang

**Initial assesment** Visited for mandibular first molar implantation due to prolonged tooth loss

**Products** **S1 bone graft material** (Powder type, 0.2~1.0mm) | **Use of membrane** ☒ Yes ☐ No

## Objectives

› Both adjacent teeth have tilted due to long-term tooth loss, and horizontal augmentation surgery is required due to buccal margin resorption.

## Conclusions

› S1 bone grafting and horizontal alveolar augmentation with a membrane were performed, and after 4 months, the alveolar bone was confirmed to be well-formed.  
› Due to sufficient implant width achieved through horizontal augmentation, the final prosthesis was stably restored after implantation.



Panoramic radiograph at the first visit



Incision at the tooth site #36



Implant placement



Insertion of S1 bone graft material followed by membrane coverage



Suture

Panoramic radiograph after the 1<sup>st</sup> surgery

Post OP (4 months and 2 weeks)



Gingival incision to check for bone regeneration

2<sup>nd</sup> surgery, healing abutment placement

Suture



Final prosthesis of the mandibular left first molar (#36)



Panoramic radiograph after the final tooth restoration

# Horizontal Alveolar Bone Augmentation

Dr. Dong-Wook Chang

**Initial assesment** Visited for mandibular 2nd molar implantation

**Products** **S1 bone graft material** (Powder type, 0.2~1.0mm) | **Use of membrane** ☒ Yes ☐ No

## Objectives

- › Observation after bone grafting for a furcation lesion and a large defect in the extraction site due to chronic periodontitis
- › Planning the implant placement for the mandibular second molar after assessing bone regeneration progress

## Conclusions

- › Horizontal alveolar augmentation was performed using S1 bone graft material at the furcation lesion site of the mandibular right first molar and the implant placement site of the mandibular right second molar. After 5 months, it was confirmed that the alveolar bone had regenerated well.
- › The results were aesthetically satisfactory.



Panoramic radiograph at the first visit



Preoperative intraoral photo



Ascertainment of bone resorption after incision



Mandibular right second molar (#47) implant placement

Mandibular right 1<sup>st</sup> molar (#45) and mandibular right 2<sup>nd</sup> molar (#46) absorbed downward to the bone level

Moldable and sticky bone graft material after hydration



Insertion of S1 to fit the defect



Membrane Insertion



Suture



## Treatment Steps



Post OP (1 month)



Post OP (1 month) panoramic radiograph



Post OP (5 months)



2<sup>nd</sup> OP, healing abutment placement



Final prosthesis



Panoramic radiograph of the final prosthesis

# Horizontal Ridge Augmentation

**Initial assesment** A vertical root fracture (VRF) occurred in tooth #27, resulting in the patient experiencing inflammation and pain.

**Products** **S1 bone graft material** (Powder type, 0.2~1.0mm) | **Use of membrane** ☐ Yes ☒ No

### Objectives

› The buccal alveolar bone of the maxillary right first premolar (#14) was very thin, so implant placement was accompanied by alveolar bone augmentation.

### Conclusions

› Horizontal ridge augmentation was performed using S1 bone graft material on the resorbed buccal side, and after 4 months, it was verified that alveolar bone and volume has been successfully regenerated.

› Healing abutment placed after the second surgery



Horizontal marginal resorption on the buccal side of the maxillary right first premolar (#14)



Confirmation of marginal resorption via incision and elevation



S1 bone graft material preparation



Bone grafting on the buccal margin following implant placement in the maxillary right first premolar



Horizontal ridge augmentation



Suture



Post OP (10 days)



Post OP (4 months)



Placement of the healing abutment after the second surgery



Panoramic radiograph after the second surgery



# Horizontal Guided Bone Regeneration

Dr. Goh

Initial assesment

Products      S1 bone graft material (Powder type, 0.2~1.0mm)    |    Use of membrane ☐ Yes    ☒ No

Objectives

› GBR was performed for the upper resorption of the maxillary left central incisor, and implant placement was carried out concurrently with the surgery.

Conclusions

› Horizontal guided bone regeneration was performed using S1 bone graft material on the resorbed buccal alveolar bone, and it was confirmed that alveolar bone and volume had been successfully regenerated 4 months later.



Intraoral photo of the upper left maxillary lateral incisor



Confirmation of gingival resorption after gum incision



Confirmation of implant penetration into the buccal soft tissue



Buccal GBR performed after implant placement in the left maxillary lateral incisor



Suture



Suture



Post OP (3 days)



Post OP (3 days)



CBCT (3 months)



Intraoral photo (4 months)

# Narrow Ridge Augmentation

Dr. Jaebum Lee

**Initial assesment** After extracting the existing tooth from the narrow buccal bone, alveolar ridge preservation was performed using S1

**Products** **S1 bone graft material** (Powder type, 0.2~1.0mm) | **Use of membrane** ☒ Yes ☐ No



## Objectives

› After tooth extraction, S1 was placed in the socket and open technique was proceeded using a PRF membrane.

## Conclusions

› Even though the open technique was proceeded without a membrane, the soft tissue healed well.

› Bone formation was splendid after 6 months.





# Treatment Steps



Preoperative X-ray



Identification of the bone defect after extraction



Application of S1 bone graft material



Covering bone with a PRF membrane



Suture



Suture



Post OP




Post OP



Post-surgery healing status



Post-surgery healing status



X-ray taken after using S1 bone graft material

# MedPark



**Tel** +82 2 554 8776 **E-mail** [biz@medpark.net](mailto:biz@medpark.net)  
**[www.medpark.net](http://www.medpark.net)**