



Protocol Description for Tatt2Away and Tatt2Revise

Please read the following carefully. Ask your Tatt2Away technician questions about any aspect concerning the Tatt2Away or Tatt2Revise protocols you do not understand.

Tatt2Away® and Tatt2Revise® protocols employ Trans Epidermal Pigment Release® or TEPR® to entirely remove or selectively modify tattoos and permanent makeup. Invented by Rejuvatek Medical Inc (RMI) and distributed by RMI since 2011, TEPR protocols are only performed by highly trained and RMI certified technicians in authorized Tatt2Away Centers.

Tatt2Away removal and Tatt2Revise modification targets all colors, ages, sizes, and styles of tattoo. Only natural products are used. TEPR protocols are designed to minimize risks of scarring and dyspigmentation. Because healing is key to removal, to obtain the best results you **MUST** follow the Post-Treatment Instructions precisely and diligently.

Oscillating needle points, held vertically in a handpiece similar to a rotary tattoo machine, are scrubbed over tattooed skin within small circular disks (called tegulae). The scrubbing and rapidly oscillating needles locally abrade epidermis, disrupt inked tissue, and suffuse Teprsol fluid to a precise depth. To aid healing, all treated tegula are isolated from one another in a nominally hexagonal treatment pattern. Within each tegula, tattoo pigment, disrupted tissue, and blood products agglomerate into separate debris scabs, which are pushed up by the reforming skin below. In two to three weeks, the ink-laden scabs will fall off to reveal fresh skin, clear of tattoo ink. These scabs **MUST** be allowed to fall off naturally.

TEPR removes tattoos piecemeal, like removing unconnected pieces from a picture puzzle. The first three treatments remove most of the tattoo. One or two more treatments are required (in Tatt2Away protocols) to scavenge remnant ink. Before treating intervening skin, prior treated tegulae must be given time to fully heal and mature. This requires a minimum of eight weeks (increasing by one or two weeks with each treatment) and possibly more, depending on how rapidly your skin heals.

If you have already attempted removal of your tattoo, we will carefully evaluate the condition of your skin. Prior laser treatment may eliminate you as a TEPR candidate.

Removal cost is calculated per session and tattoo area (usually, by the number of tegulae treated). Because a minimum cost applies to every session, it may be advantageous to treat several small tattoos at the same time. The maximum areal size treated in one session is limited to two standard 3×4 inch templates (2×138 tegulae). The maximum arc length treated around an extremity is limited to half its circumference. Sixty minute sessions are usually required to treat one template of 138 tegulae.

Geanes Aesthetic Spa
29650 Bradley Rd, Ste A. Menifee CA 92586
(951) 672-8226
www.GeanesAestheticSpa.com