

HUMA Media Technology

GMP-Certified Universal Multi-purpose Media For Adult Stem Cells



ISO 13485 Certified

HUMA Media Custom Production is available



Email: info@rokit.co.kr | Tel: +82 1899 7296 | Web: rokithealthcare.com

The Unprecedented Media for Human Primary MSC culture

Meticulously Phased Smart Media

- * Cell attachment (Phase I)
- * Stabilization/proliferation (Phase II-III)
- * Proliferation & Stemness Maintenance (Phase IV)

Product	HUMA-Media Tech Xeno-Free GMP-certified MSC Media Kit		
Application	Cell culture and expansion	Storage	Shelf-life 6 Months at 4 °C
Cell Type	Human Primary/Stem Cells; Neonatal Dermal Fibroblasts, Bone-Marrow MSC, Adipose derived MSC, Smooth Muscle Cells, Dermal Papilla Cells, Warton's Jelly MSC, Dental Pulp Stem Cells	Notes	No plate coating needed
		Grades	500 ml RUO (custom manufacturing available)

HUMA Media Tech Xeno-Free MSC Media Kit: Media I & II

※ Phase I Media & Phase III Media can be custom manufactured

HUMA-Media Recommended Protocol 1 (Media I: 500ml , Media II: 500ml)

Phase I (Day 0) & Phase II (Day 2)

Mix Media I & Media II at a 1:1 ratio to make a Hybrid Media. Use desired volume (i.e. 10 ml) for seeding. On day 2, without discarding existing Media, just add 75% of the media volume at seeding (i.e. 7.5 ml).


Phase III (Day 4~)

On day 4, discard entire existing Media (i.e. 17.5 ml) and add fresh 10 ml of Hybrid Media. From day 4, replace Hybrid Media for every 2 days.

HUMA-Media Recommended Protocol 2 (Media I: 500ml , Media II: 500ml)

Phase I (Day 0) & Phase II (Day 2)


Media I	300ml
Media II	75ml
Final Vol(ml)	375ml



Mix Media I (300 ml) & Media II (75 ml) at a 4:1 ratio to prepare Phase I Media. Use desired volume (i.e. 10 ml) for seeding. On day 2, without discarding existing Media (i.e. 10 ml), just add Phase I Media with 75% of initial seeding volume (i.e. 7.5 ml).

Phase III (Day 4~6)

Media I	120ml
Media II	120ml
Final Vol(ml)	240ml



On day 4, prepare Phase III Media by mixing Media I (120 ml) & Media II (120 ml) at a 1:1 ratio. On day 4, discard entire Media in the dish and freshly add Phase III Media same as initial seeding volume (i.e. 10 ml).

Phase IV (Day 6~)

Media I	-
Media II	305ml
Final Vol(ml)	305ml



From day 6, Media II is Phase IV Media. Discard entire Media in the dish and add the same volume of Phase IV Media (i.e. 10 ml). From day 6, change Media for every 2 days.

*10ml of Media per dish as an example volume

HUMA-Media Tech Xeno-Free GMP-certified MSC Media Kit Performance




Fig.1 HUMA Media growth rate comparison with commercially available serum free media | HUMA media growth rate and accumulated total cell number were compared with two commercially available serum free media products, Stem*** and Mesen*** up to passage 8. PDT (Population Doubling Time) was around 20-24 hrs, at least 2 folds faster than competing products throughout the entire culturing and passaging periods. Accumulated total cell number at passage 8 was at least 3 folds higher than competing commercial serum free media products.

Proliferation




Fig.1 | Population doubling time is tightly maintained during extended passaging period for all three MSC types.




Fig.2 | Population double level (PDL) showed high performance; All cell types showed doubling within 20 hrs range during late passage stages. BMSC and WJSC showed average 10 hrs PDL.

Cell size





Fig.3 | Endogenous MSC size is known as near 10 µm range and primary MSC size after first passage has known to be 14 µm~20 µm. Cell size is directly related to cell senescence and performance as an adult stem cell. MSC culture in HUMA-media over 7 passages shows excellent cell size maintenance in all three stem cell types; Bone Marrow MSC, Adipose derived MSC, and Wharton's Jelly MSC.

Cell Viability




Fig.4 | In HUMA-media, cell viability has maintained over 98% in relatively old passaging stage for all three MSC cell types.

Cell surface marker (FACS)




Fig.5 | Highly distinctive MSC surface marker has been identified for passage 7 cell types (BMSC, ADSC, and WJSC).

Differentiation




Fig.6 | Expansion with HUMA-Media retains tri-lineage passage 5, 6, 7 differentiation potential.