

How to use the SoberSense Application: iPhone

1. Open SoberSense Application

2. Enter login credentials

- a. Username- firstnamelastname
- b. Password- sober
- c. ALL LOWERCASE WITH NO SPACES
- d. You are now looking at the sessions for the day. The times showing are the end of the window. Any test submitted after will be reported as late.

3. We can now turn on the breathalyzer

- a. *Hold the gray button down for a couple of seconds until the blue light turns on. (This button only turns the device on and off. Hold again to turn off.)*

4. After the blue light is on at the top you are ready to select your session

- a. **Always take your time at this point*. **It should always be the top test.***

5. Tap the middle of the screen

- a. This will start the syncing process.

6. Tap on "Take Test"

7. You will see a countdown on the top right of the screen. Once it is at 0, it will say, "Blow Now"

- a. *You will see a battery voltage notification show up. Disregard that. It will still show even with full battery.*

8. Blow steadily for about 6 seconds until you feel a click inside the breathalyzer.

9. Submit Results!!

10. Tap on the tap right refresh arrows to confirm that the session has disappeared.

