

Characterization of Irradiance in Light Microscopes that use the Villuminator™ Module

Ronian Siew*

*Optical Consultant

Venture Biotech Modules Business (V-BMB) Pte. Ltd.



Acknowledgements



- We are grateful to Nathalie Gaudreault* and Laurent Gelman** for their kind invitation
- Thanks also to Stanley Schwartz‡ for introducing us to QUAREP-LiMi

* Director Microscopy Pipeline at [Allen Institute for Cell Science](#)


** Head of [Facility for Advanced Imaging and Microscopy](#)

‡ Consultant, (Former) Senior Advisor at [Nikon Instruments](#)

Irradiance measurement at the specimen plane of an objective (example)



This is for the 480 nm illumination (surface plot is across 250 x 250 pixels, each pixel is 2.4 μm square)




1 Measure total flux at the specimen plane using power meter

2 Capture beam profile at specimen plane using an image sensor (focus by first reducing iris at the field stop)

3 Scale pixel digital intensities using Matlab, Excel, etc.
(Here, Excel was used)


Scaling of pixel digital intensities (method 1)

Resized raw image into 250 x 250 pixels (easier to handle in Excel)




1 Save image file from ImageJ into .CSV text (power meter measured total flux $\approx 34 \text{ mW} \pm 0.3\%$)

MS Excel® Spreadsheet




2 Open text file with Excel, scale cell values to flux values (total pixel shot noise $\ll \pm 0.3\%$)

MS Excel® Spreadsheet




3 Divide flux values by pixel area to yield flux per unit area

Scaling of pixel digital intensities (method 2)




MS Excel® Spreadsheet



	A	B	C
1	Value	Count	Value x Count
255		253	0
256		254	0
257		255	0
258		SUM of Value x Count →	5850508
259		Total Flux (mW) →	34
260		Scaling Factor →	5.81146E-06
261			

ImageJ 3D Plot



1 Display histogram of pixel digital intensities in ImageJ, click “List” to output histogram data

2 Column C = Value x Count; Sum all, then divide flux by sum to yield scaling factor

3 Digital intensity values in 3D plot must multiply scaling factor, yielding irradiance in mW/pixel

Caveat 1

Field number divided
by magnification

Flux in the FOV

This is a
relative term

We assume that the signal to background ratio is “high”

Mean flux from light scatter
that's been integrated by
power meter and image sensor

Caveat 2


What if you don't have space to mount a camera at the specimen plane?

In this case, one option is to try this experimental setup




The flat mirror can be any high quality 'first surface' mirror. To focus onto the flat mirror, close the iris but leave sufficient opening so as to see the iris's inner rim, then focus until the rim appears sharp. Let the camera record the irradiance distribution. Note that the irradiance profile at the camera may not necessarily be the same as the irradiance at the mirror. In this case, perhaps you can 'baseline' your measurement by first mounting a Villuminator module as the light source to obtain a top hat irradiance profile with > 95% uniformity at the mirror plane. Then, let $V(x, y)$ be the irradiance at the mirror from the Villuminator, and $V'(x, y)$ be its image in the camera. If $T(x, y)$ is the 'transmittance function' of the imaging system, then $V'(x, y) = V(x, y)T(x, y)$. So, if $I'(x, y)$ is the image of your original irradiance profile and $I(x, y)$ is your original irradiance at the mirror, then $I'(x, y) = T(x, y)I(x, y)$. Thus, your original irradiance at the mirror is $I(x, y) = [V'(x, y)/V(x, y)]I'(x, y)$. But since $V(x, y)$ is a top hat, it is spatially constant (within 5% error). Therefore, $I(x, y) \approx V'(x, y)I'(x, y)$.

SNR* is independent of the beamsplitter's %R/%T split ratio



*SNR \equiv Signal to noise ratio


The Villuminator™ modules produce very flat irradiance distributions



So, it is of interest to characterize the irradiance profile in terms of analytic formulas


Analytic expressions relating irradiance distribution with total flux and 'width'

Rotationally Symmetric Gaussian



$$E(\rho) = \frac{2P}{\pi w^2} \exp\left[-2\left(\frac{\rho^2}{w^2}\right)\right]$$

Rotationally Symmetric 'Super-Gaussian'




$$E(\rho) = \frac{nP}{2\pi w^2 \Gamma(2/n)} \exp\left[-\left(\frac{\rho}{w}\right)^n\right]$$

Gamma function

In both formulas, P is the total flux in the distribution

Properties of the 'super-Gaussian top-hat'



➔ Average irradiance in a circle of area $\pi\rho^2$ is $= \frac{2\pi}{\pi\rho^2} \int_0^\rho E(\rho)\rho d\rho = \frac{2}{\rho^2} \int_0^\rho \frac{nP}{2\pi w^2 \Gamma(2/n)} \exp\left[-\left(\frac{\rho}{w}\right)^n\right] \rho d\rho$

➔ In the limit that $n \rightarrow \infty$, the distribution  approaches 


➔ $\lim_{n \rightarrow \infty} E(0) = \frac{P}{\pi w^2} = \frac{P}{\pi R^2}$ There is no distinction between R and w when the distribution is an ideal top hat

Why irradiance matters to us



- A source's 'brightness' is proportional to the irradiance* in the image of that source
- Fluorescence emission is proportional to the irradiance from the illumination
- Given a spatial irradiance distribution $E(x, y)$, the integral $\iint_{\text{Area}} E(x, y) dx dy$ yields flux in the 'Area'
- The flux delivered to an area is at least partially absorbed by matter in that area
- If a distribution is 'top hat', then the irradiance times any area yields flux in that any area ✓

This is what we're after


Globally non-uniform irradiance profile
(but locally uniform)




Globally uniform irradiance profile
(but locally non-uniform)



*R. Siew, [Eur. J. Phys.](#) **29**, 1105 (2008); R. Siew, [Eur. J. Phys.](#) **43**, 035304 (2022)




I JUST NEED A
LITTLE SHADE.
A LITTLE 'SLOPE'
ON THE SIDES
IS OK.




Lublu

Producing top hat illumination using the Villuminator™ module: principle of operation




Example

Experimental Setup




Measured Profile




Fluorescence Microscope Setup

(Nikon 10x Plan Fluorite, 0.3 NA)




Fluorescence Microscope Setup

(Edmund Optics 10x 120i Plan Apo, 0.28 NA)



Tunability of the irradiance distribution



By varying the working distance, the irradiance profile can be made to change from a dome to a top hat and to an inverted-dome.



This means that the FINAL profile at the camera can be made flat even if the 'response' of a fluorescent sample (or the transverse transmittance function of the imaging system) has either a dome-shaped or inverted-dome shaped profile.

Example of tuning the profile

(Villuminator with Edmund Optics 10x 120i Plan Apo, 0.28 NA)

Working Distance = 18 mm

Working Distance = 36 mm



White Papers at <https://v-bmb.com>



WHITE PAPER
Dec 2023, Rev. 2

Introduction to Using Villuminator™ Modules for Producing Uniform “Top Hat” Illumination

Villuminator modules are compact modular subassemblies that produce exceptionally uniform illumination across a surface. This property benefits many applications involving the need to illuminate planar regions of interest, such as spatial light modulators in projection systems and specimen planes being observed under a microscope. In this white paper, we explain how to use Villuminator modules for a variety of illumination applications.



1. INTRODUCTION

In the Jan/Feb 2023 issue of *Biophotonics*, Advanced Products Corporation (APC) described how biotechnology applications are benefiting from using special devices called Villuminator™ modules, which are compact optical systems that produce output beams with



Figure 1. Compact standard Villuminator modules.

exceptionally uniform light distributions [1]. Soon after, the Venture Biotech Modules Business (V-BMB) was formed, offering a variety of precision control solutions, and expanding the Villuminator module product line, providing analytical instrument developers and researchers with unparalleled uniform illumination.

Villuminator modules have been offered in a variety of sizes and designs. For example, the three modules shown at the top right of this page have diameters between 50 and 60 mm, and length of roughly 130 mm (inclusive of a heat sink for the source). Their output

beams are in the red, green, and blue. Other designs are smaller, such as the four modules shown in figure 1. These compact standard designs have dimensions of 40 mm (diameter) x 78 mm (length). The four modules in figure 1 produce output beam colors in green, blue, red, and white. Other colors (wavelengths) are also available, ranging from the near ultraviolet through the visible (including white) and near infrared.

2. PRIMARY COMPONENTS OF A VILLUMINATOR MODULE

Figure 2 shows the components of a Villuminator module. A module consists of a LED, heatsink, power cord, and a patent-pending optical system that projects a non-divergent beam out of the module.




Figure 2. Components of a Villuminator module.


Introduction to Using Villuminator Modules for Producing Uniform “Top Hat” Illumination
©2023 V-BMB | <https://v-bmb.com>

1



WHITE PAPER
Feb 2024, Rev. 2

Selecting Commercial Off-The-Shelf Coupler Lenses for Villuminator™ Modules Used in Fluorescence and Bright-Field Microscopy



Epifluorescence and bright-field microscope systems benefit from using Villuminator modules for illumination because of the exceptionally high level of field uniformity provided by the Villuminator modules. For best results, an appropriate “coupler lens” should be selected to focus light from the Villuminator module into the microscope objective. This white paper explains how to make such selections.

1. INTRODUCTION

Villuminator™ modules are compact optical systems that produce output beams with exceptionally uniform light distributions [1–4]. When applied to fluorescence and/or bright-field microscopy, a coupler lens must be used to focus the output from a Villuminator module into the objective lens, as illustrated in figure 1.

Generally, for standard Villuminator modules [5], any commercial-off-the-shelf (COTS) AR-coated achromatic doublet lens having an effective focal length (EFL) of 150 mm (and diameter in the range between 40 to 50 mm) can be used as a coupler lens (e.g., the lens with part number 49-285 from Edmund Optics, or part number AC508-150-A from Thorlabs, can be used as a coupler lens). However, there may be instances where a coupler lens with a different EFL is desirable. In the next section, we provide basic rules for selecting coupler lenses.

2. HOW TO SELECT A COUPLER LENS

2.1 Condition of using objectives from different manufacturers

Figure 2 illustrates the basic rule for coupler lens selection in situations where you may be using objectives from different manufacturers (note that a fluorescence filter cube has not been included in the figure, as it would be a redundant component for the current discussion). The idea is that, at any specified






Figure 1. Use of a coupler lens to focus light from a Villuminator module into an objective lens.

Selecting Commercial Off-The-Shelf Coupler Lenses for Villuminator Modules Used in Fluorescence and Bright-Field Microscopy
©2024 V-BMB | <https://v-bmb.com>

1

Article in Biophotonics (Jan '2023)



Top Hat Illumination Provides Even Light Distribution Across Samples

LEDs and geometrical optics principles combine to provide effective lighting profiles for microscopy, patterned projection, polymerase chain reaction, and other forms of analysis.

BY ROMAN SEW AND LHAO TAN, ADVANCED PRODUCTS CORP. PTE. LTD.

Even illumination patterns are of great importance when applied to biomedical imaging because they prevent lighting artifacts that could distort the picture of the sample being analyzed. Therefore, suppliers of microscopy solutions have made the integration and replication of these even patterns a priority in their commercial products^{1,2}. With this goal in mind, any bright or dark features of a specimen on a microscope slide must be accurately reproduced for purposes of examination, whether by the human eye or by capturing the specimen's image via a camera sensor contained within a digital microscope. Clarity is particularly important when analyzing conditions in cells, tissues, embryos, or other specimens for research or clinical diagnostics.

Two common illumination techniques used in microscopy are known as critical illumination and Köhler illumination, and each has advantages and disadvantages. Other lighting techniques are also avail-

able for a variety of biophotonics applications. A new lighting method that uses lenses with LEDs to produce flat, top hat illumination has emerged, yielding even distribution of illumination across any plane of interest. This approach greatly benefits life scientists who examine dynamic and intricate samples.

Classical illumination techniques

To appreciate the nuances of how illumination problems can manifest themselves in biomedical imaging, it is important to understand what uniform or even illumination represents (Figure 1) in this context. In optical system design, the technical term used for describing illumination profiles is the irradiance distribution, which refers to the flux per unit area across the plane of illumination.

A dome-shaped irradiance distribution may be considered globally nonuniform across the sample, but it may be locally uniform in specific areas of the sample.

And yet, globally uniform irradiance profiles can still have hot and cold spots — small regions of the illumination profile that have either excessive or insufficient brightness — which may be considered to be localized lighting artifacts in an image. These artifacts may arise from structural properties of the source, such as the coils of tungsten filaments, or even from wires and mechanical joints that run across LED surfaces.

Figure 2 illustrates the basic optical layout of five well-known classical techniques of illumination. For comparison, each system is depicted to project equal flux across equal illumination areas at a screen. Typically, the technique of critical illumination enables projecting a flat irradiance distribution, provided that the source is locally uniform. But sources often do not possess such uniformity.

For instance, if the source is a tungsten filament, then the image of the filament is projected onto the plane of illumination.

If an LED is used, then wires and other electrical and mechanical structures that run across the emitter's surface are projected to the illumination plane. Examples of images displaying these obstructions under a microscope, such as in fluorescence images of cells, are discussed in several technical white papers³.

In critical illumination, sometimes a simple defocus of the flux collector lens (or adding an engineered diffuser) is sufficient to homogenize the local nonuniformities, but this depends on the spatial dimensions and gaps between localized structures of the source. Köhler illumination — which is the most common technique used in bright-field and fluorescence microscopy — eliminates such localized nonuniformities. But if the source is small, such as an LED with a small emitter, then the result is often a globally nonuniform profile, even if edge vignetting is not present. However, in microscopy, other sources of nonuniform patterns in Köhler illumination can also appear in an image⁴.

Collimated illumination typically does not provide globally uniform irradiance distributions, due to the cosine fourth law as it applies to a Lambertian LED source — that is, an LED whose datasheet indicates that its radiant intensity drops to roughly 50% of its peak at an output ray angle of 60° normal to the LED's surface. This principle indicates that when such an LED is located at the back focal plane of a lens, the irradiance across a distant screen in front of the lens can take on a globally nonuniform profile. But if the LED is sufficiently large and locally uniform, then it is possible to produce a uniform illumination profile, provided that the plane of illumination is extended farther away.

This is known as the "searchlight" condition.

Finally, the last two techniques shown in Figure 2 can produce uniform irradiance distributions on a plane, but at the expense of using exotic components, such as fly's eye arrays and light pipes, which can result in expanded system footprints and extended lengths or add to design complexity and cost.

For example, Taro Ichimura and his colleagues have explored the use of fly's eye lenses to produce uniform illuminations at the specimen plane of a so-called trans-scale (multicellular) imaging microscope possessing a large field of view. But their system layout appears rather long⁵.




Figure 1. Global and local nonuniformity of an illumination profile across a plane.




Figure 2. Five well-known classical techniques of illumination.

By distilling the effects of system length into a mathematical formula, it can be shown that a purely geometrical first-order equation governs the condition of equal flux projected onto equal areas of illumination for the critical and Köhler illumination techniques:

$$\frac{1}{R} = \frac{1}{r} + \frac{1}{L}$$

Here, R is the semidiameter of a projection lens (shown in the second system from the top in Figure 2), r is the semidiameter of the flux collector lens, and L is the radial height of the illumination profile measured from the optic axis. For a specified source size, flux collector size, and size of the illuminated area, this

formula determines the size of the projection lens in Köhler illumination such that the total flux and size of the illumination area are the same in both the critical and Köhler techniques.

Since the projection lens is conjugate to the source — that is, the projection lens's plane is the location of the image of the source, which is formed by the flux collector — the formula also dictates the focal length of the flux collector for any choice of distance between the source and the flux collector. Finally, since the flux collector is conjugate to the area of illumination — that is, the plane of illumination is the location of the image of the flux collector, which is formed by the projection lens — this determines the system length.

Courtesy of Laurin Publishing

More info at the V-BMB website



The screenshot shows the V-BMB website interface. At the top, a navigation bar includes the V-BMB logo with the tagline "Precision Control Solutions Made Simple", and menu items for HOME, PRODUCTS, RESOURCES, SHOP, and ABOUT. A callout menu is open under the RESOURCES link, listing: TEC Controller, TEC Assembly, Villuminator™, Motorized Filter Wheel, Motion Actuator, Syringe Pump, Whitepapers, and Videos. A white arrow points from a box labeled "RESOURCES MENU" to the callout menu. Below the navigation is a large image of a multi-channel liquid handling workstation. To the right of the image, the text "INNOVATIVE CONTROL SOLUTIONS" is partially visible. Below this, a message reads "Have questions? We have answers for you!" with two buttons: "LEARN MORE" and "BUY NOW". At the bottom, the email address "sales@v-bmb.com" and a chat icon are visible.



Visit us!
Hall A3
Booth 515A



analytica

APRIL 9-12 | 2024

analytica.de/en

Munich, Germany

Thank You!



<https://v-bmb.com>