

To export STL files from DEXIS IS systems, follow these general steps:

- **Open the ScanFlow Software:** Launch the DEXIS IS ScanFlow software where your intraoral scans are stored.
- **Select the Scan:** Choose the specific scan or case you want to export.
- **Export as STL:**
 - Look for the **Export** or **Save As** option in the software menu.
 - Select the **STL format** as the file type for export. DEXIS IS systems support open STL files, making them compatible with third-party CAD/CAM or 3D printing systems.
- **Save the File:**
 - Choose the destination folder on your computer or external storage.
 - Confirm the export, and the STL file will be saved for further use.

If you encounter any issues, consult the DEXIS IS ScanFlow user manual or contact DEXIS support for detailed guidance.