

AMAZING SKIN AESTHETICS

Phone: (714) 848-2232 | Huntington Beach, CA

CLINICAL FAQ & COMPREHENSIVE MARKET COMPARISON GUIDE

DENSITY™ HIGH-FREQUENCY RADIOFREQUENCY BY JEISYS

Welcome to your streamlined guide to **Jeisys DENSITY™ High-Frequency Radiofrequency** tissue remodeling. Combining patient experience with technical data, this document breaks down common clinical questions and outlines why DENSITY™ represents the absolute elite standard compared to alternative RF devices on the market.

1. GENERAL & PATIENT TREATMENT FAQS

Q: How exactly does the DENSITY™ system work?

A: DENSITY™ delivers targeted high-frequency 6.78MHz radiofrequency energy deep into the dermal skin layers. This energy safely heats localized tissue to induce controlled thermal coagulation, shrinking loose collagen chains and kicking off a progressive biological mechanism where your body builds fresh, structural collagen and elastin fibers for deep lifting and firming.

Q: What should I expect immediately after my treatment? Is there downtime?

A: DENSITY™ is non-invasive with zero true downtime. Immediately following your session, your skin will feel deeply warm (resembling a moderate sunburn) with diffuse redness (erythema) and mild swelling (edema). These are healthy clinical endpoints indicating successful tissue heating and typically resolve entirely within 24 to 48 hours.

Q: Is the treatment painful? How is comfort managed?

A: Unlike older legacy devices, DENSITY™ is highly tolerable. It incorporates a patented dynamic cryogenic gas-cooling membrane that continuously cools the epidermal surface before, during, and after each RF pulse. This masks the deep thermal profile, rendering the procedure comfortable without aggressive topical or oral sedation.

Q: Are there any transient side effects I should know about?

A: Yes. Because RF energy deeply stimulates beneath the epidermis, some patients experience mild texture variations or dry, sandpaper-like tightness around Days 3 to 5 as cellular turnover speeds up. Transient folliculitis (small, temporary bumps due to hair follicle thermal stimulation) can also surface and clears naturally within 7 days. Mild asymmetry from uneven swelling or localized bruising may rarely occur but resolves quickly on its own.

Q: Who is the ideal candidate for DENSITY™? Can it treat all skin types?

A: DENSITY™ is safe and effective for all Fitzpatrick skin types (I-VI). It is ideal for individuals experiencing mild to moderate skin laxity, fine lines, wrinkles, jowl formation, or periorbital/abdominal skin sag. Because RF energy bypasses melanin pathways, it carries an exceptionally low risk of post-inflammatory hyperpigmentation (PIH) compared to traditional lasers.

Q: Why do you recommend coming back in 3 months for another treatment?

A: While immediate tightening occurs due to heat-induced collagen contraction, true biological tissue remodeling takes 3 to 6 months to mature. Clinical data shows that stacking a subsequent treatment at the 3-month mark allows us to build a secondary structural matrix of new collagen exactly as your first treatment reaches its peak performance—compounding and maximizing your long-term lifting longevity.

2. MARKET COMPARISON: WHY DENSITY™ IS GLOBALLY ELITE

The aesthetic market contains various radiofrequency (RF) options, but advanced engineering directly determines your clinical safety, patient comfort, and total lifting outcome:

DENSITY™ vs. Thermage® FLX (The Legacy System)

While Thermage® has long been a household name for monopolar RF, DENSITY™ introduces critical advancements. DENSITY™ features a patented **Sequential Monopolar and Bipolar Energy Delivery System** within a single tip. Thermage® only delivers monopolar energy, which treats deeply but can feel intensely hot and cause uncomfortable spikes. DENSITY™ alternates monopolar and bipolar energy dynamically, treating both deep structural layers and superficial skin textures simultaneously for a more comprehensive lift. Furthermore, DENSITY™ features an advanced real-time skin impedance sensor that measures tissue resistance on *every single shot*, adjusting energy instantaneously to maximize comfort and reduce painful heat spots common with Thermage®.

DENSITY™ vs. Xerf® (The Competitor System)

Xerf® attempts to replicate high-frequency skin tightening but lacks the safety profile and custom configurations of the Jeisys system. DENSITY™ integrates a robust, patented eco-friendly gas-cooling membrane that continuously monitors epidermal temperature during delivery. This protects the outer skin surface from burns while pushing larger bulks of thermal energy deeper into the skin. Xerf's cooling mechanism is less dynamic, limiting total energy delivery and requiring more sessions to match DENSITY's single-session results. Additionally, DENSITY™ offers larger body tips (such as the DENSITY 8.0) explicitly designed to adapt to abdominal and large-surface laxity seamlessly.

DENSITY™ vs. Standard Microneedling RF & Fractional RF Devices

Standard microneedling RF devices require physical needle penetration to deliver energy, which requires topical numbing, carries risks of pinpoint bleeding, and involves a multi-day healing downtime. DENSITY™ is completely non-invasive and needle-free, achieving equivalent deep-dermal heating profiles comfortably through surface contact. Furthermore, standard fractional RF or old fractional systems only treat a percentage of the skin surface, leaving gaps in the treatment grid. DENSITY's patented continuous contact membrane treats 100% of the targeted tissue area uniformly, eliminating untreated gaps for a smoother, tighter result.

3. DIRECT TECHNOLOGY COMPARISON MATRIX

CLINICAL FEATURE	JEISYS DENSITY™	THERMAGE® FLX	XERF® SYSTEM	STANDARD RF / NEEDLING
Energy Options	Sequential Monopolar & Bipolar delivered in a single pulse.	Strictly Monopolar delivery only.	Basic Monopolar or Bipolar configurations.	Basic Fractional or Needle-bound delivery.
Epidermal Cooling	Patented dynamic gas-cooling membrane with real-time epidermal protection.	Standard cryogen cooling bursts. Higher risk of thermal spikes.	Basic contact cooling, limiting max safe thermal bulk delivery.	Often lacks cooling entirely due to needle insertion.
Impedance Feedback	Real-time automatic resistance tracking on 100% of delivered shots.	Periodic calibration adjustments, leading to occasional hot spots.	Standard fixed calibration metrics.	Lacks individual pulse calibration feedback loops.
Skin Coverage	100% uniform surface coverage across the entire treatment block.	100% uniform surface coverage.	Standard surface coverage.	Fractional coverage only (leaves untreated skin gaps).
Downtime & Comfort	Zero downtime. High comfort via advanced rapid gas-cooling cycles.	Zero downtime. Moderate to low comfort (often painful).	Zero downtime. Variable comfort with risk of hot spots.	2-4 days downtime. Requires topical numbing (invasive needles).

Experience the future of aesthetic medicine. Elevate your contours with the unmatched speed, safety, and precision of Jeisys DENSITY™ at Amazing Skin Aesthetics.

Amazing Skin Aesthetics • (714) 848-2232 • Huntington Beach, CA