



# Getting Started with Your RPM Device

## Introduction

Welcome to Aries Solutions! Remote Patient Monitoring (RPM) is designed to help you manage your health effectively from the comfort of your home. This guide will walk you through setting up and using your RPM device to ensure a smooth experience.

### Step 1: Unboxing and Checking Your Device

- Carefully remove your device from the packaging.
- Verify that all components are included (device, charger, user manual, and any accessories).
- Check for any visible damage before use.

### Step 2: Charging and Powering On

- Plug your device into a power source using the provided charger.
- Allow it to charge fully before first use.
- Press and hold the power button to turn it on.

### Step 3: Connecting to a Network

- If your device is cellular-enabled, ensure it has a good signal.
- For Wi-Fi-enabled devices, follow the on-screen instructions to connect to your home network.
- If required, use Bluetooth to sync with your smartphone or other devices.

### Step 4: Setting Up Your Profile

- Follow the on-screen prompts to enter basic health information.
- If your device syncs with an app, download the recommended app and create an account.
- Adjust settings as needed for personalized tracking.

### Step 5: Taking Your First Reading

- Ensure you are in a comfortable position before taking measurements.
- Follow the device instructions carefully (e.g., placing a blood pressure cuff correctly, using a glucose monitor properly).
- Wait for the device to record your data and check the display for results.

### Step 6: Transmitting Data to Your Provider

- Your RPM device will automatically send data to your healthcare provider if properly set up.
- If manual submission is required, follow app instructions to upload readings.
- Contact your healthcare provider if data transmission issues arise.

### Step 7: Troubleshooting and Support

- If your device isn't working correctly, restart it and check connections.
- Refer to the user manual for troubleshooting common issues.
- Contact Aries Solutions customer support for further assistance.

### Conclusion

Your RPM device is a powerful tool for managing your health. Regular use and accurate data submission help your healthcare provider monitor your condition and adjust your care plan accordingly. If you have any questions, reach out to Aries Solutions for guidance!

