Disposable Face Mask

Packings:

Size: 17.5*9.5cm

Quantity: 50/box 2000/Ctn

Company Profile

Established in 2001, it is a company integrating the production, sales and trade of non-woven products. We have more than 200 employees, including 10 R&D personnel, 20 professional technicians and 10 quality management staff.

With two separate factories, our company covers a total area of nearly 60,000 square meters including a Class 100,000 cleanroom and Class 10,000 laboratory.

Our products are widely used in medical protection, personal care, industrial protection, food hygiene and other fields. We have a sound quality management system.

Our main products passed the CE certificates & ISO9001. Some products have been registered in FDA. We have also obtained registration & production license for Class II medical devices in China. Our face masks have been tested by Nelson lab in the United States to meet the highest standards of EN14683-2014 and ASTM F2100-11(2018).

We always put the quality in the first place and keep improving production technology and service level. We would work actively and try our best to provide high quality goods and the better services for overseas and domestic customers.

Certifications

一、FDA Approved

Fiscal Year 2020 CERTIFICATION OF REGISTRATION

This certifies that:

has completed the FDA Establishment Registration (as manufacturer , contract manufacturer) and Device Listing with the US Food & Drug Administration, through

U.S. Agent for FDA Communications: **SUNGO TECHNICAL SERVICE INC.**
6050 W EASTWOOD AVE APT 201, CHICAGO,
ILLINOIS 60630, USA
Telephone: +1-855-957-7779 / E-mail: sungo.group@yahoo.com

Registration Number: 3008048818

Device Listing#: See annex

SUNGO Technical Service Inc. will confirm that such registration remains effective upon request and presentation of this certificate until the end of the calendar year stated above, unless said registration is terminated after issuance of this certificate. SUNGO Technical Service Inc. makes no other representations or warranties, nor does this certificate make any representations or warranties to any person or entity other than the named certificate holder, for whose sole benefit it is issued. This certificate does not denote endorsement or approval of the certificate-holder's device or establishment by the U.S. Food and Drug Administration. SUNGO Technical Service Inc. assumes no liability to any person or entity in connection with the foregoing.

Pursuant to 21 CFR 807.39, "Registration of a device establishment or assignment of a registration number does not in any way denote approval of the establishment or its products. Any representation that creates an impression of official approval because of registration or possession of a registration number is misleading and constitutes misbranding." The U.S. Food and Drug Administration does not issue a certificate of registration, nor does the U.S. Food and Drug Administration recognize a certificate of registration, SUNGO Technical Service Inc. is not affiliated with the U.S. Food and Drug Administration.

Executive Director
Issued: Dec. 26 2019
Cert. No.: 2006US819518
Expiration Date: Dec. 31 2020

二、NQA

Certificate of Registration

This is to certify that the Quality Management System of

XIANTAO SANDA INDUSTRIAL CO., LTD.
 Unified Social Credit Code : 9142900473088731XN
 Operation Address : No.46, East Section of Huangjin Avenue, Xiantao City, Hubei Province, China
 Registered Address : No.46, East Section of Huangjin Avenue, Xiantao City, Hubei Province, China

applicable to

Manufacture and sales of PP non-woven series products (caps, shoes cover, sleeves, working clothing and face masks)

has been assessed and registered by NQA against the provisions of

ISO 9001:2015

This registration is subject to the company maintaining a quality management system, to the above standard, which will be monitored by NQA.
 Certified Clients shall accept regular surveillance assessments, the validity of certificates shall be maintained for the positive result of audit.
 The information of this certificate can be checked on CNCA's website (www.cnca.gov.cn) SNQA's website : www.snqa.com.cn

Mingyu
 Managing Director

Certificate Number: **45960**

Date: 24 January 2019
 Valid Until: 24 January 2022
 EAC Code: 04

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Mingyu
 Managing Director

Certificate Number: **45960**

Date: 24 January 2019
 Valid Until: 24 January 2022
 EAC Code: 04

The use of the UKAS Accreditation Mark indicates accreditation in respect of those activities covered by the accreditation certificate number 015 held by NQA.
 NQA is a trading name of NQA Certification Limited, Registration No: 02021708. Registered Office: Watwick House, Houghton Hall Park, Houghton Regis, Buckingham, LU5 5ZK, UK.
 This certificate is the property of NQA and must be returned on request.

三、CE

Verification of Conformity

Applicant: XIANTAO SANDA INDUSTRIAL CO.,LTD
Address: NO.46 GOLDEN AVENUE,XIANTAO,HUBEI,CHINA
Product(s): Gownall, Lab Coat, Shoe Cover, Pillow Case, Sleeve Cover, Apron, Face Mask, Cap, Surgical Gown, Patient Gown, Isolation Gown, Bed Sheet/Bed Cover
Type(s): See annex
Product Classification: Class I

The submitted technical files including test report of the above products have been reviewed against the self declaration requirements of conformity for CE marking according to Medical Device Directive (93/42/EEC).

Standard(s) used for showing compliance with the essential requirements in the specified directives:
Standard(s): EN ISO 14971:2012; EN ISO 15223-1:2016;
 EN 1041:2013; EN ISO 10993-1:2009/AC:2010;
 EN ISO 10993-5:2009; EN ISO 10993-10:2013

The review result of the technical files and test report support the self declaration for the devices listed above. Where the manufacturer affix's the CE marking to the product listed they must ensure that all the requirements of the appropriate EU directive(s) have and continue to be met.

SUNGO Cert GmbH

Executive Director
 Issued: Dec. 10 2019
 Cert. No.: EU156518
 Expiration Date: Dec. 9 2024

四、Nelson Test Report

Nelson Labs.
A Sotera Health company

Sponsor:
Mingmin Hu
Xiantao Sanda Industrial Co., Ltd.
No. 46 Huangjin Ave. (East Section)
Xiantao, Hubei Province, 433000
CHINA

Study Number: 1033510-S01
Latex Particle Challenge Final Report

Test Article: SD20180308-1
SD20180308-2
SD20180308-3
SD20180308-4
SD20180308-5

Study Number: 1033510-S01
Study Received Date: 23 Mar 2018
Testing Facility: Nelson Laboratories, LLC
6280 S. Redwood Rd.
Salt Lake City, UT 84123 U.S.A.

Test Procedure(s): Standard Test Protocol (STP) Number: STP0005 Rev 05
Deviation(s): None

Summary: This procedure was performed to evaluate the non-viable particle filtration efficiency (PFE) of the test article. Monodispersed polystyrene latex spheres (PSL) were nebulized, dried, and passed through the test article. The particles that passed through the test article were enumerated using a laser particle counter.

Three one-minute counts were performed, with the test article in the system, and the results averaged. Three one-minute control counts were performed, without a test article in the system, before and after each test article and the counts were averaged. Control counts were performed to determine the average number of particles delivered to the test article. The filtration efficiency was calculated using the average number of particles penetrating the test article compared to the average of the control values.

The procedure employed the basic particle filtration method described in ASTM F2299, with some exceptions, notably the procedure incorporated a non-neutralized challenge. In test use, particles carry a charge, thus this challenge represents a more natural state. The non-neutralized aerosol is also specified in the FDA guidance document on surgical face masks. All test method acceptance criteria were met. Testing was performed in compliance with US FDA good manufacturing practice (GMP) regulations 21 CFR Parts 210, 211 and 820.

Test Side: Inside
Area Tested: 81.5 cm²
Particle Size: 0.1 µm
Laboratory Conditions: 21°C, 23% relative humidity (RH) at 1040, 20°C, 25% RH at 1403
Average Filtration Efficiency: 98.9%
Standard Deviation: 0.09

Study Director: Brandon L. Williams
Study Completion Date: CHA PE

801-290-7500 | nelsonlabs.com | sales@nelsonlabs.com | info@nelsonlabs.com | 801-290-7500 | nelsonlabs.com | sales@nelsonlabs.com | info@nelsonlabs.com
Page 1 of 2

Nelson Labs.
A Sotera Health company

Study Number 1033510-S01
Latex Particle Challenge Final Report

Results:

Test Article Number	Average Test Article Counts	Average Control Counts	Filtration Efficiency (%)
1	172	14,013	96.8
2	131	11,853	96.9
3	140	12,677	96.9
4	175	13,246	96.7
5	159	13,458	96.8

801-290-7500 | nelsonlabs.com | sales@nelsonlabs.com | info@nelsonlabs.com | 801-290-7500 | nelsonlabs.com | sales@nelsonlabs.com | info@nelsonlabs.com
Page 2 of 2

Nelson Labs.
A Sotera Health company

Sponsor:
Mingmin Hu
Xiantao Sanda Industrial Co., Ltd.
No. 46 Huangjin Ave. (East Section)
Xiantao, Hubei Province, 433000
CHINA

Study Number: 1033509-S01
Bacterial Filtration Efficiency (BFE)
and Differential Pressure (Delta P) Final Report

Test Article: SD20180308-1
SD20180308-2
SD20180308-3
SD20180308-4
SD20180308-5

Study Number: 1033509-S01
Study Received Date: 23 Mar 2018
Testing Facility: Nelson Laboratories, LLC
6280 S. Redwood Rd.
Salt Lake City, UT 84123 U.S.A.

Test Procedure(s): Standard Test Protocol (STP) Number: STP0004 Rev 15
Deviation(s): None

Summary: The BFE test is performed to determine the filtration efficiency of test articles by comparing the bacterial control counts upstream of the test article to the bacterial counts downstream. A suspension of *Staphylococcus aureus* was aerosolized using a nebulizer and delivered to the test article at a constant flow rate and fixed air pressure. The challenge delivery was maintained at 1.7 - 2.7 x 10⁷ colony forming units (CFU) with a mean particle size (MPS) of 3.0 ± 0.3 µm. The aerosols were drawn through a six-stage, viable particle, Andersen sampler for collection. This test method complies with ASTM F2101-14, EN 14683:2014, Annex B, and AS4381:2015.

The Delta P test is performed to determine the breathability of test articles by measuring the differential air pressure on either side of the test article using a manometer, at a constant flow rate. The Delta P test was designed to comply with MIL-M-39594C, Section 4.4.1.2 and complies with EN 14683:2014, Annex C and AS4381:2015.

All test method acceptance criteria were met. Testing was performed in compliance with US FDA good manufacturing practice (GMP) regulations 21 CFR Parts 210, 211 and 820.

Test Side: Inside
BFE Test Area: ~40 cm²
BFE Flow Rate: 28.3 Liters per minute (L/min)
Delta P Flow Rate: 8 L/min
Conditioning Parameters: 85 ± 5% relative humidity (RH) and 21 ± 5°C for a minimum of 4 hours
Test Article Dimensions: ~160 mm x ~156 mm
Positive Control Average: 2.5 x 10⁷ CFU
Negative Monitor Count: <1 CFU
MPS: 2.9 µm

Study Director: Janelle R. Bentz, M.S.
Study Completion Date: CHA PE

801-290-7500 | nelsonlabs.com | sales@nelsonlabs.com | info@nelsonlabs.com | 801-290-7500 | nelsonlabs.com | sales@nelsonlabs.com | info@nelsonlabs.com
Page 1 of 2

Nelson Labs.
A Sotera Health company

Study Number 1033509-S01
Bacterial Filtration Efficiency (BFE)
and Differential Pressure (Delta P) Final Report

Results:

Test Article Number	Percent BFE (%)	Delta P (mm H ₂ O/cm ²)	Delta P (Pa/cm ²)
1	99.3	3.1	30.4
2	99.4	3.0	29.5
3	99.4	3.3	32.4
4	99.3	3.0	29.6
5	99.7	3.2	31.6

The filtration efficiency percentages were calculated using the following equation:
% BFE = $\frac{C - T}{C} \times 100$
C = Positive control average
T = Plate count total recovered downstream of the test article
Note: The plate count total is available upon request.

801-290-7500 | nelsonlabs.com | sales@nelsonlabs.com | info@nelsonlabs.com | 801-290-7500 | nelsonlabs.com | sales@nelsonlabs.com | info@nelsonlabs.com
Page 2 of 2

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Performed in the microbiological lab under the responsibility of Yvette Register

Analysis Report 18.02617.01
Date 26-06-2018
Page 6/8

Reference: T1810792 - Type II Surgical Mask
Lot NR. 181629

Microbial cleanliness on masks

Date of ending the test: 05-06-2018
Standard used: EN 14683 - §5.2.5 (2014)
Product standard: EN 14683 (2014)
Number of tested masks: 5
Extraction liquid: Peptone 1 g/l, NaCl 5g/l & Tween 20 2g/l
Extraction volume: 300 ml
Extraction time: 5 min.
Counting technique: Membrane filtration
Filtration volume: 100 ml
Culture media: TSA (Tryptic Soy Agar)
Incubation conditions: SDA (Sabouraud Dextrose Agar with chloramphenicol)
3 days at 30°C (TSA)
7 days at 20-25°C (SDA)

Results

# Mask	Mask weight (g)	CFU ^a /mask		Microbial cleanliness	
		Aerobic microbial count (bacteria)	Fungi count (SDA)	Σ CFU/mask	Σ CFU/g
1	2.64	6	<3	<9	<4
2	2.60	9	<3	<12	<5
3	2.63	6	<3	<9	<4
4	2.58	18	9	27	11
5	2.63	12	3	15	6

* CFU : Colony Forming Unit

Performed under accreditation in the microbiological lab under the responsibility of Yvette Register

EN14683 Type IIR

Analysis Report 18.02617.02
Date 26-06-2018
Page 1/9

Your notice of 03-05-2018 Your reference Date 26-06-2018

Analysis Report 18.02617.02

Required tests:
EN 14683 (2014) EN 14683 - annex B (2014) Bacterial filtration efficiency
EN 14683 (2014) EN 14683 - annex C (2014) Surgical masks - Breathability (differential pressure)
EN 14683 (2014) EN 14683 - §5.2.5 (2014) Microbial cleanliness on masks
EN 14683 (2014) ISO 22669 (2004) Surgical masks - Splash Test

Identification number	Information given by the client	Date of receipt
T1810793	Type II R Surgical Mask Lot NR. 181629	02-05-2018

Yvette Register

Order responsible

This report may be reproduced, as long as it is presented in its entire form, without written permission of CENBEL.
The results of the analysis cover the received samples. CENBEL is not responsible for the representativeness of the samples.
In assessing compliance with the specifications, we did not take into account the uncertainty on the test results.

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Analysis Report 18.02617.02
Date 26-06-2018
Page 2/9

Reference: T1810793 - Type II R Surgical Mask
Lot NR. 181629

Bacterial filtration efficiency

Date of ending the test: 31-05-2018
Standard used: EN 14683 - annex B (2014)
Product standard: EN 14683 (2014)
Mask description: Non woven masks - 3 ply (blue outside/white inside)
Number of tested masks: 5
BFE Area tested: ± 46 cm²
Masks conditioning: 21 ± 5°C and 85 ± 5% RH
Side of the mask in contact with the bacterial challenge: Inner side
Challenge bacterial strain used: *Staphylococcus aureus* ATCC6538
Bacterial challenge per test: 2200 ± 500 CFU
Total test time: 1 min. delivering challenge + 1 min. without challenge (air flow continuing)
Flow rate: 28.3 l/min.
Positive control: Tests performed with no filter material in the air stream
Negative control: Test performed without challenge

Performed in the microbiological lab under the responsibility of Yvette Register

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Results

B = Bacterial filtration efficiency (%)

$$B = \frac{(C - T)}{C} \times 100$$

With C = mean of the total plate counts for the positive control runs
T = total count for the tested mask

# Mask	B (%)
1	> 99,9
2	99,9
3	99,9
4	99,6
5	99,9

Mean particle size of the bacterial challenge aerosol: 3,0 µm

Controls

Mean positive controls 2395 CFU
Negative control < 1

Performed in the microbiological lab under the responsibility of Yvette Rogister


Reference: T1810793 - Type II R Surgical Mask
Lot NR. 181629

Surgical masks - Breathability (differential pressure)

Date of ending the test 30-05-2018
Standard used EN 14683 - annex C (2014)
Product standard EN 14683 (2014)

Mask description Non woven masks – 3 ply (blue outside/white inside)
Number of tested masks : 5
Number of areas per mask 5 (see figure)
Dimension of the areas : Disc whose diameter is 2.5 cm
Surface areas : 4,9 cm²
Flow rate : 8 l/min.
Direction of the air flow : From the inside of the mask to the outside
Masks conditioning : 21 ± 5°C and 85 ± 5% RH

Figure : Distribution of the areas in the mask

Results Δ P

Performed in the microbiological lab under the responsibility of Yvette Rogister

	Mask 1	Mask 2	Mask 3	Mask 4	Mask 5
Area A	33,0	39,1	34,0	34,2	37,7
Area B	31,6	34,4	35,4	34,0	36,7
Area C	34,2	33,2	35,2	30,4	34,6
Area D	32,2	33,6	32,2	32,6	30,6
Area E	31,6	32,6	31,6	32,4	30,6
Average Δ P (Pa/cm ²)	32,5	34,6	33,7	32,7	34,0

Performed in the microbiological lab under the responsibility of Yvette Rogister

Reference: T1810793 - Type II R Surgical Mask
Lot NR. 181629

Microbial cleanliness on masks

Date of ending the test 01-06-2018
Standard used EN 14683 - §5.2.5 (2014)
Product standard EN 14683 (2014)

Number of tested masks 5
Extraction liquid Peptone 1g/l, NaCl 5g/l & Tween 20 2g/l
Extraction volume 300 ml
Extraction time 5 min.
Counting technique Membrane filtration
Filtration volume 100 ml
Culture media TSA (Tryptic Soy Agar)
Incubation conditions SDA (Sabouraud Dextrose Agar with chloramphenicol)
3 days at 30°C (TSA)
7 days at 20-25°C (SDA)


Results

# Mask	Mask weight (g)	CFU*/mask		Microbial cleanliness	
		Aerobic microbial count (bacteria)	Fungi count (SDA)	Σ CFU/mask	Σ CFU/g
1	3.04	3	<3	<6	<2
2	2.99	6	3	9	4
3	3.00	<3	<3	<6	<2
4	2.97	3	<3	<6	<3
5	3.02	6	<3	<9	<3

* CFU : Colony Forming Unit

Performed under aeration in the microbiological lab under the responsibility of Yvette Rogister

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NOTE:
The physical properties mentioned above represent this lot of non-woven fabric. These values should not be considered as specification values of this type of material. Shelf time: 1 year (From date of manufacture).

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