

THE REGENETRON PEMF SYSTEM



The Regenatron is a Solid-State Digital unit. Easy to use, 5 separate settings or use 2 at a time. Light weight for moving around or travel. It operates at 2544 Gauss.

It is ready for 110v USA or 220v International power requirements. Weight 3 lbs. Safety tested and ready for use, 3-year warranty.

This unit has two complete separate machines inside. One of the machines delivers the primary treatment frequency and the other delivers a sweep of frequencies. This technology provides a much more effective treatment than every other competing PEMF machine. The device is comprised of two independent therapy devices that can be used together or separately.

The device provides Pulsed Electro Magnetic Fields (PEMF) stimulation and healing and offers seven different treatment cycles that emulate different brain wave frequency ranges and one treatment cycle for PEMF subtle energy therapy.

PEMF THERAPY SOLUTIONS

424-272-6011

GETTING STARTED:

Connect the power cube into the MAINS and into the MAINS socket on the back panel of the device.

Always plug a treatment coil connector into the connector socket on the back panel of the device before selecting a treatment.

INSERTING/CHANGING A TREATMENT

Do not plug or unplug a treatment coil during the treatment cycle.

Insert a coil connector by pushing it into the connector socket and twisting clockwise to lock it in place. Remove the treatment coil by twisting the connector counterclockwise and pulling it out.

BEGINNING YOUR TREATMENT

Users can remain fully clothed and no direct contact between the coil and the skin is necessary for the treatment to be effective.

Place the treatment coil on or around the desired treatment area. The closer the coil is to the treatment area, the more effective the treatment will be.

Press the desired push button to begin one of seven preset brainwave treatment cycles along with the Subtle Energy Therapy treatment. The corresponding blue light comes on. For one of the brainwave treatments alone, press the Subtle Energy pushbutton to turn it off.

The Subtle Energy treatment can be used alone by pressing its corresponding pushbutton. You may switch back and forth between treatments by pressing their respective pushbuttons. Each treatment has a corresponding blue light blinking pattern.

PAUSING/RESETTING/ENDING TREATMENT

Pressing a lit pushbutton after treatment has begun places the device on pause. Pressing the lit pushbutton once more resumes the treatment.

Holding a lit pushbutton down until the device beeps, stops the treatment cycle.

Extensive literature and research correlate the different brain wave frequencies to different states of brain activity and states of mind:

- Delta brain waves are associated with deep sleep (1 to 3 pulses per second)
- Theta brain waves are associated with light meditation, the transition before sleep and just after waking (4 to 7 pulses per second)
- Alpha brain waves are associated with deep relaxation and creativity (8 to 12 pulses per second)
- Beta brain waves are associated with a heightened state of alertness, focused concentration and critical reasoning (13 to 30 pulses per second)
- Active Alert (13 to 100 pulses per second)
- Gamma waves are associated with bursts of insight and high-level information processing (30 to 100 pulses per second)

The different programs provided by the device create an environment for the brain to align with a range of PEMF with frequencies that emulate the different brain wave states.

BEEP CODES SEQUENCES:

- Single Beep = End of Session. Press any treatment pushbutton to begin a new treatment.
- Three Short beeps = Unit is overheating. Shut down and let the device cool down for ten minutes. Restart by pressing any treatment pushbutton.

INDICATIONS FOR USE

- Reduce pain
- Reduce inflammation
- Heal faster
- To be used with other protocols to open the cells for deeper healing and penetration
- Improve your energy
- Enhance your brain wave states

This device has three components, a control unit that generates pulsed electromagnetic fields of different intensities, a treatment coil that transfers the therapeutic pulsed electromagnetic fields to the body and a power supply that is plugged into the mains and into the back of the control unit to provide electrical power. It is designed to allow users to safely treat inflamed and painful areas of the body when used according to the specified conditions of use. It is important to read this manual before using it. No undesirable side effects have been reported with use the use of this device. Pre-clinical testing has demonstrated that this device is safe for use when used according to the instructions in this manual.

PUSH the “Subtle Energy Therapy” button (the one with a little play sign), “push play” and that is your 10 pulses per second. It is great for inflammation, for cell charge, and for lymphatic drainage or the chelation of toxins or heavy metals. This mode has the benefits of bone healing and reducing inflammation – perfect for healing after surgery or strengthening the body systems.

This unit transmits deep healing to all cells, to repair, replace and detox so these cells can better absorb nutrients and push out toxins.

USING THE BODY MAT AND SLEEP PADDLE

The mat and the paddle are the perfect combination for sleep, a quick recharge or a focused brain balance. You can start and stop in case you need a break during your session. Just push once to start, and when you push the button again it will start to blink, and the pulse will stop. When you're ready to continue just push it again and start it back up. If you want to restart from the beginning just push and hold and the machine will beep and go back to a green screen to begin again.

SLEEP AND HEAL

Let's talk brain waves...and using the paddle and the brain waves buttons down here at the bottom of the unit. At the side of the machine are the Delta, Theta and Sleep Sweep buttons. These are particularly for people with trouble sleeping or waking up exhausted. You can use the 8 hours of sleep with either the Delta or Theta pulse rates or use the Sleep Sweep, which combines the two.

You might want to try: one night using Delta, one night using Theta, and one night using the Sleep Sweep and see how you sleep by placing the paddle under your pillow (it is going to pulse the two waveforms through your brain). The paddle may get a little bit warm under your pillow, but this is normal.

MEDITATIVE AND CALMING

People with a lot of stress or needing to ground themselves a little bit better would use the Alpha mode. Or if you want to meditate, and slow down your mind, use either the mat or the paddle.

When the “Alpha button,” is pushed, there is a 7.8 pulse rate which is designed to be very similar to the Schumann resonance 7.83, and basically this resonance is going to allow the brain to sleep into waveforms that are more calm, relaxing, meditative and grounding.

ACTIVE ALERT

The next part of the machine is “Active Alert”. The Theta and Gamma are faster brain waves so you can use them together or try them individually each morning. If you want to Theta by itself, see how you feel. Gamma by itself, see how you feel. Or just do the active alert for both. Have fun with this setting and it might be a good idea to keep a journal to document what and how you feel after using each setting.

BE MORE ALERT IN THE AM

This is a super-fast pulse rate and it accesses a higher brain wave form to make you more alert and more clear-minded. For students or older individuals, this higher brain wave makes processing a little faster, and helps you process information and focus better. About 30 minutes in the morning and you will notice a nice, clear focused mind throughout the day.

TRAIN YOUR BRAIN

PEMF Therapy when applied using applicators such as pads or pillows stimulates the brain using pulsating magnetic fields, which improves circulation and neuro-chemical imbalances that cause depression.

By enhancing the cells and nutrient intake, this unit also balances hormonal functions including serotonin (which helps you feel calm, serene, optimistic and self-confident), dopamine (responsible for making you feel excited, motivated, and energized) and the stress hormone, cortisol (which revs you up into a high gear when you need it).

Regular therapy at home is easier to fit in the schedule and benefits can be enjoyed long-term, providing mental and cognitive performance enhancement, it's much more than a little boost of serenity. PEMF therapy also has anti-inflammatory effects on the brain tissues too just like any other part of the body. Brain inflammation can also be responsible for brain fog or decline in cognitive function, depression, stress and anxiety. Research also suggests that you can reverse these trends and take charge of your brain health in as little as two weeks if you give your brain the frequencies it needs to recharge itself to operate effectively. Research also proves pulsed electromagnetic field therapies (PEMF) have beneficial effects upon Alzheimer's disease and its underlying neurophysiological abnormality.

In-depth research using magnetic field therapy for the entire body, or on all the cells of the body, show that cellular functions are also affected. For example, during deep sleep, we tend to regenerate more. When sleeping all night in magnetic field tuned to Schumann resonance amazing health benefits have been reported. Usually reported at 7.8 Hz, the Schumann resonance ensures life on earth.

Schumann – Alpha brain waves are present in deep relaxation and usually when the eyes are closed, when you're slipping into a lovely daydream or during light meditation. It is an optimal time to program the mind for success and it also heightens your imagination, visualization, recall, learning and concentration.

The Schumann resonance is known to be the earth's frequency. Using the Alpha setting is a great way to synchronize with the earth and ground yourself. Try this cycle to when you are feeling Anxious or stressed.

Helmholtz Field. Helmholtz field is created when you bring two magnetic fields together. Sandwich your body between the mat and the paddle to charge up those stubborn low cell membranes between the attachments to recharge and reset.

PLAY AROUND WITH IT

Experiment with it, play with it, see how this machine can help your brain and your body with pain and inflammation and regeneration. Use the Therapy Mat to detox and enhance the absorption of nutrients. Improve sleep, use all brainwaves (delta and theta) when using the Alpha and Active Alert settings. Use both the paddle and the mat to create a healing field of energy and promote healing your mind and body all at one.

You will notice relief from pain, inflammation and exhaustion, and the all-important benefit of healing while you get Regenerative sleep.

CONTRAINDICATIONS

DO NOT USE within 25cm of pacemakers, defibrillators or any other implanted electronic devices. **DO NOT USE** soon after taking any medication.

Consult your health care provider to discuss using the device if you think you are at risk for any of these contraindications.

PRECAUTIONS

Caution: Users with heart problems should consult their physician before using the device.

Caution: Users with suspected or diagnosed epilepsy should consult their physician before using the device.

Caution: Users taking pain, anxiety, depression or any other medication should be carefully monitored when using the device as medication effectiveness may be intensified.

USER CONTROLS AND SYMBOLS

The device has eight different push button switches that each begins a treatment for a preset amount of time in a specific range of frequencies.

There are no known or reported adverse or allergic reactions with the use of this device. Detoxification may occur after treatment. Drink plenty of water after treatment.

In case of any adverse reaction, stop using the device and consult a physician.

WARNINGS

DO NOT use the device near credit cards, security access cards, car keys, hearing aids, watches, jewelry and metal chains, cell phones, iPODS, laptops, remote controls, or any other electronic media. The electromagnetic fields may disrupt their functioning and/or demagnetize them.

DO NOT USE in wet environments. Do not immerse any part or pour any liquids on the device. Keep away from sources of heat and moisture.

Keep the device out of the reach of children. Children may be at risk of strangulation with the cord and treatment coil pigtails or of asphyxiation with the packing materials.

TROUBLESHOOTING

With the PEMF power cube connected to the MAINS and to the device, press any of the treatment pushbuttons on the top of the device control unit. Each treatment has a corresponding blue light blinking pattern.

If the light does not come on, does the control unit make any beeping sounds?

If the control unit emits a beep sound, proceed to identify the beep sequence to remediate.

MAINTENANCE Service

This device has an expected service life of five years and no components degrade abnormally over time. The device has no user serviceable parts and must be returned to the manufacturer for servicing. Additional treatment coils and power cords are available from the manufacturer.

REMOVING THE COVER OR TAMPERING VOIDS THE WARRANTY.

Contact the manufacturer if something unexpected happens or for assistance in setting up, using or maintaining the device.

This device is an electronic device. Electronics should never be disposed of with regular trash. Take non-working electronics to an electronics recycling center.

CLEANING

There is no mandatory or scheduled cleaning, maintenance or sterilizing necessary. If you choose to clean the device:

DO NOT immerse any part or pour any liquids on the device, | Disconnect the device from the MAINS before cleaning.

STORAGE & TRANSPORTATION

Keep the device dry and store in a dry place. Storage in a damp place may cause corrosion. Store and transport at temperatures between -25 degrees centigrade to +70 degrees centigrade at a relative humidity of up to 93% non-condensing.

Remove the PEMF power cube and the detachable treatment coil from the control unit and pack all parts of the device securely before transporting.

WARRANTY

The Warranty Registration Card must be filled out and returned to the manufacturer within 30 days of the date of purchase to activate the warranty.

The manufacturer warrants this device to operate properly for a period of three years from the date of the original purchase invoice. In the event of a malfunction during the warranty period, the manufacturer will, at its discretion, replace or repair the device to its original operating condition. Freight and insurance to and from the manufacturer's repair facility are not included. Freight and Insurance are the responsibility of the registered owner.

A RMA number must be obtained from the manufacturer prior to returning any device or accessory for service. The device must be delivered freight paid with the RMA, the registered owner's name and address, and a brief description of the difficulties encountered. The device is to be shipped to the address designated by the manufacturer with the RMA number.

Such service, repair or adjustment of the device is guaranteed to the original purchaser provided the device has not been tampered with, does not have any physically broken parts and the control unit was not opened, altered or damaged as a result of misuse, accident, water, grit, impact, or lack of proper care.

CONTACT INFORMATION:

Regeneration Devices Ltd. +44 1172 050411 info@regenerationdevices.com

Vesper Manufacturing 14556 Weddington Street Sherman Oaks, CA 91411 U.S.A

(818) 267-0095

<p style="text-align: center;">BRAIN WAVES</p> <p>Sleep Sweep</p>	<p>This setting will steer your brain to the Delta and Theta Brainwaves to fall asleep, stay asleep and wake up refreshed. Most people put the paddle under their pillow for the 8-hour cycle.</p>
<p>Delta</p>	<p>Delta is the slowest band of brainwaves. When your dominant brainwave is delta, your body is healing itself and “resetting” its internal clocks. You do not dream in this state and are completely unconscious</p>
<p>Theta</p>	<p>Theta Reduced consciousness, deep meditation, dreams, light sleep, REM sleep and extreme relaxation. Hypnosis and self-programming. Increased creativity and retention of learned material</p>
<p>Beta</p>	<p>Beta brain waves are associated with normal waking consciousness and a heightened state of alertness, logic and critical reasoning. While Beta brain waves are important for effective functioning throughout the day.</p>
<p>Gamma</p>	<p>Gamma is associated with the formation of ideas, language, and recall processing and various types of learning.</p>

Alpha	Alpha brain waves are present in deep relaxation and usually when the eyes are closed, when you're slipping into a lovely daydream or during light meditation. It is an optimal time to program the mind for success and it also heightens your imagination, visualization, recall, learning and concentration. The Schumann resonance is known to be the earth's frequency. Using the Alpha setting is a great way to synchronize with the earth and ground yourself. Try this cycle to when you are feeling Anxious or stressed.
Helmholtz Field	A Helmholtz field is created when you bring two magnetic fields together. Sandwich your body between the mat and the paddle to charge up those stubborn low cell membranes between the attachments to recharge and reset
Active Alert	This setting will steer your brainwaves to the Beta and Gamma Brainwaves to be more alert and attain better focus and concentration. Most people use the paddle under their pillow and start their day with this setting.



Great for all ages...

PEMF THERAPY SOLUTIONS
424-272-6011