

HOW TO USE YOUR TELEHEALTH REMOTE STETHOSCOPE PEDESTAL



THE FIRST STEP IN ALL TELEMEDICINE CONSULTS IS TO TELEPHONICALLY CONTACT THE MEDICAL PROVIDER

Make sure the Devices are charged, and the cables are plugged fully into the stethoscope & interface. To get to the Home screen place a small pen into the hole on the left side of the display and tap lightly until the home screen (As seen above) appears.

Alexa Communicate & Drop-In Screen

- To enter the Alexa app, tap the Alexa Icon twice
- Wait a few seconds and the Alexa app should appear
- Tap on the communicate Icon at the bottom of the screen
- The Drop-In Icon will then appear at the top of the screen
- Tap on the Drop-In Icon



Selection of Telehealth Provider to connect with

- With the device in the Drop-In mode, you will see several devices listed such as Dr. Jones Stethoscope, Nurse Smith Stethoscope, Etc. (please use the devices with the word stethoscope)
- Select the person you wish to conduct the telemedicine consultation with and tap on that name.

Communication between the Patient and Doctor

- When the connection between the Patient & Doctor starts the screen may be blurry for a few seconds, this allows the other person to close the camera lens if needed.
- The Screen will clear and both parties should be able to see each other. Because the audio feed is used by the stethoscope, the Nurse or Doctor will not hear your voice through the pedestal a second device such as a phone or Telemedicine Tablet will be needed in conjunction with the Pedestal. Conduct the Telemedicine Consult as long as required

Disconnecting the Conference

- Once the telemedicine consultation is complete, tap the display screen and then tap the red Phone Icon to disconnect.

See [TELEHEALTHFL.COM](https://www.telehealthfl.com) for more information