



INFRARED THERMOMETER











INSTRUCTIONS

Index

1. Symbols used
2. Introduction
3. Caution
4. Components
5. Display screen
6. Features
7. Temperature ranges
8. Operation
9. Setting
10. Warning
11. Care of the product
12. Troubleshoot
13. Compliance

1.Symbols used

	DC 3V2X15V AAA
	Type BF equipment
	Disposal should be done in accordance with the rules/legislation applicable in India specified for the disposal of electronic waste. Do not dispose the unit in general household waste
RoHS	The device is in accordance with ROHS 2011/65/EU
	The material is PVC
	CHN Green Point Mark
	Keep the unit in a dry area away from moisture or excessive humidity
	Mandatory reference manual
	This symbol indicates instructions or directions that must be followed in order to avoid damages to user or other people. to avoid damages to user or other people

2 Introduction

The product is a professional thermometer for measuring the forehead temperature of adults and children


The temperature display depends on the status of the skin.

3. Caution

- Please read this user manual before switching the unit on.
- Operating temperature is between 10°C and 40°C.
- Please keep the thermometer away from electrical appliances.
- The unit should be stored only between 0°C to 50°C and below 85%RH.
- Keep away from sources with electro-magnetic radiation such as radios, mobiles, etc.
- Please do not expose the unit to sun, fire or water.
- Keep the unit safe from any impact. Do not use the unit in case of damage.
- Please ensure there is no sweat on the forehead and do not wear a hat or scarf as this would influence the measurement values.
- Ensure that the measuring distance is between 1-3 cm.
- In case of rapid changes in ambient temperature, use the unit only after 30 minutes to ensure more accurate measurement values.
- In case of sweat on the forehead or abnormal body temperature, please take the temperature of the earlobe.
- TO clean, use alcohol to lightly wipe the surface.
- This thermometer is not intended to be used as a substitute for consultation with your physician or doctor.
- Please do not attempt to open or repair the product. For repair of the product, please customer care.


4. Components

- A.LCD display area
- B.setting
- C.MODE
- D.detecting area
- E.Down
- F.Up
- G.testing key
- H.Battery cover



5. Display screen

- A.testing result
- B.Memory
- C.Alarm signal
- D.Body surface mode
- E.Internal mode
- F.battery mark
- G.Display mode in $^{\circ}\text{F}$
- H.Display mode in $^{\circ}\text{C}$



6. Features

- Non-contact and with a high accuracy of measurement.
- Two kinds of measurement mode-body temperature or Surface temperature.
- Exchange between C and $^{\circ}\text{F}$
- Sound alert function to notify high body temperature.
- Storage capacity of the last 32 sets of measurement data.
- Function to save the data and to Auto Turn off.
- Function to auto select the range of measurement and LCD display resolution- 0.1°C .(0.1°F)
- Backlight LCD digital display

7. Temperature ranges

a. Measuring position


Measuring position	Temperature in °C	Temperature in °F
Anus	36.6-38.0	97.8-100.4
Oral	35.5-37.5	95.9-99.5
Axilla	34.7-37.3	94.4-99.1
Ear	35.8-38.0	96.4-100.4

b Normal body temperatures for different age groups

Age group(years)	Temperature in °C	Temperature in °F
0-2	36.4-38.0	97.5-100.4
3-10	36.1-37.8	97.0-100.
11-65	35.9-37.6	96.6-99.7
65	35.8-37.5	96.4-99.5

8. Operation

- Insert the battery.
- Wait for 10 minutes for the first time use.
- In case of use after a long period, please wait for 1-2 seconds for the unit to switch on as it detects the environment temperature initially.
- Focus on the forehead keeping a distance of 1-3 cm away and then press the button to start. The temperature will be displayed instantly. Please ensure the forehead is free of hair, sweat, cosmetic or hat.



In case of rapid changes in environment temperature, use the unit for checking the body temperature by focusing behind the ear. Keeping a distance of 1-3 cm away.



9. Setting

9. 1 Temperature setting-f1

- Keep pressing the MODE button for 3 seconds till the screen display shows F1. Choose the ▼ arrow for °C and the ▲ arrow for °F.





9.2 Alarm-F2

- Keep pressing the **MODE** button and continue till the screen display shows F2.
- Press ▲ and ▼ arrow for increasing and decreasing the alarm temperature by 0.1 °C (0.1 °F)


Attention: The default alarm temperature is 38 °C (100.4 °F)

9.3 Buzzer On/Off-F3


- Press **MODE** button for 3 seconds and then continue to press mode button two times more. So that the display shows F3.

- Press  arrow button to close the buzzer. (The sound icon  will disappear)
- Press  arrow button to open the buzzer. (The sound icon  will display on LCD screen)



9.4 Abnormal display

- In case of any erratic display, press the  and switch off and re-start the unit.
(The unit is meant for measuring the temperature of the human body within the range of 32.0 to 42.5°C or 89.0 to 108.0 °F.)

9.5 Measuring mode

- Press the  arrow to change the mode from Body to Surface temperature measurement.
(The temperature range in this mode is 0.0 to 100.0°C or 32.0 to 212.0 °F.)

9.6 Data saving

- The test results are automatically saved and is displayed at the right bottom of the screen
- Press  or  arrow to check the history

9.7 Change the battery

- The battery should be replaced when the display shows "sign Steps.
- Open the battery shell and remove the battery
- Insert fresh batteries in to the compartment with the correct polarity+-, otherwise the unit may be damaged
- Please do not use Rechargeable batteries
(In case the unit is not used for a long period, make sure to remove the battery.)

9.8 Device life

- The unit can be used for taking temperatures for 40,000 times.

9.9 Protection against electrical shock

- AS per Type BF.

10. Warning



- Handle the display screen with utmost care to avoid damage.
- Do not use rechargeable batteries.
- Do not use the unit under direct sunlight.
- Keep away from water.

11. Care of the product

- Use the scan head with care as this is the most important part of the unit.
- TO clean, use cotton dipped in 70% alcohol.
- Keep the unit in a dry place, avoiding contact with water and away from sunlight and dirt.

12. Troubleshoot

- Please follow the below troubleshooting to attend to issues arising during operation of the unit.
- If the problem persists, contact customer care.

Sr.no.	Problem	Display	Reason	Solution
1	No display in LCD screen		<ul style="list-style-type: none"> In Body Temperature mode the temperature is below 32°C or above 42.5°C 	Use within the specified range
2	Display shows "Lo" or "HI"	 	<ul style="list-style-type: none"> Surface covered by sweat or hair Unstable air environment Unit is too far from scan area Rapid and large change from indoor to outdoor temperature 	<ul style="list-style-type: none"> Clear the scan area Wait till the environment temperature stabilizes Use the unit within 1-3cm from the scan area Wait for 20 minutes before taking the body temperature

13.compliance

- Although the product is in compliance with EN 60601-1-2 certification for electro-magnetic radiation, care should be taken to avoid such interference from such radiation while using the unit.
- This product conforms to the following standards.
- EN 980: Graphical symbols for use in the labelling of medical devices
- EN: 1041: Information supplied by the manufacture with medical devices
- EN: 60601-1: Medical electrical equipment Part 1: General requirements for safety (IEC: 6060-1)
- EN 60601-1-2: Medical electrical equipment Part 1-2: General requirements for safety Collateral standard