


Rezzimax® Tuner Pro II Rezzipes



TUNER PRO
II


REZZIMAX®
Tune Out Pain. Tune Into Life.



Introduction

This Rezzipe Book is a step-by-step guide intended to assist you in using the Rezzimax Tuner Pro II. Follow each step as you would a recipe. It will provide you with basic techniques for supporting your body's natural healing. It is divided head-to-toe based on where pain or tension is currently felt. After following these, feel free to develop your own Rezzipes and let us know! For additional support, view our [online training videos](#) or enjoy this [complimentary training course](#).

The Basics:

To incorporate the Vagus Nerve (the nerve in charge of natural healing) there are a few things that pertain to each of the steps within the Rezzipes.

1. Place your tongue between your teeth (this helps to unclench your jaw)
2. Hum gently along with the resonance of the Tuner Pro II (this helps your Vagus Nerve to magnify the positive effects)
3. Press the two circle buttons simultaneously, they will flash once released.
4. There are 4 modes; green/blue/red/white that change to the next by pressing the button between the two arrows. Completion of a mode is signaled by two distinct pulses.
5. If you are sensitive to vibration, please begin with the Sensitive Rezzipe in this guide.

Table of Contents

Head & Neck	<u>4</u>
Jaw & TMJ	<u>6</u>
Loss of Taste and Smell	<u>8</u>
Shoulder & Chest	<u>10</u>
Elbows & Hands	<u>12</u>
Stomach	<u>14</u>
Back	<u>15</u>
Legs	<u>17</u>
Knees	<u>18</u>
Feet & Ankle	<u>20</u>
Whole Body De-stress	<u>22</u>
Sensitive	<u>24</u>
Primitive Reflexes	<u>26</u>
Stress Down-Regulation	<u>33</u>
Appendix	<u>36</u>

Head & Neck

This Rezzipe supports your body's natural healing response for headaches, neck pain, etc.



Step 1:

- Optional: Use Neck Pillow
- Tuner: Accessory Tongs
- Mode/Level: Blue
- Duration: 5 min
- Pressure: 1
- Location: Back of Neck

Attach tongs to Tuner and place tongs on both sides of the neck with Tuner resting against the spine slowly rotating head side to side, keeping pressure against the device.



Step 2:

- Tuner: Accessory Tongs
- Mode/Level: Green
- Duration: 1-2 min
- Pressure: 1
- Location: Under jaw

Place tongs under the jaw holding the Tuner with both hands



Step 3:

- Tuner: Accessory Tongs
- Mode/Level: Choose a Comfortable Level
- Duration: 1-2 min
- Pressure: 1
- Location: Center of Eyebrows



Place tongs at center of eyebrows

Step 4:

- Tuner: Accessory Tongs
- Mode/Level: Desired Level
- Duration: 1-2 min
- Pressure: 2
- Location: Between Eyes

Place tongs close together between the eyes at the top of the nose with light pressure towards your forehead



Step 5:

- Tuner: Accessory Tongs
- Mode/Level: Desired Level
- Duration: 7-8x
- Pressure: 1
- Location: Back and Forth Over Head

Run the tongs over the top of your head back and forth 7-8 times



Step 6:

- Tuner: Accessory Tongs
- Mode/Level: Desired Level
- Duration: 1-2 min
- Pressure: 1
- Location: Cheekbones

Place tongs on your cheekbones



Step 7:

- Tuner: Accessory Tongs
- Mode/Level: Desired Level
- Duration: 7-8x
- Pressure: 2
- Location: In Mouth

Place tongs inside your mouth with the tongs stretched outside the teeth. Angle up and down and open and close your mouth slowly 7-8 times

Jaw & TMJ

This Rezzipe supports your body's natural healing response for jaw pain.



Step 1:

- Optional: Use Neck Pillow
- Tuner: Accessory Tongs
- Mode/Level: Blue
- Duration: 5 min
- Pressure: 1
- Location: Back of neck

Attach tongs to Tuner and place tongs on both sides of the neck with Tuner resting against the spine slowly rotating head side to side, keeping pressure against the device.



Step 2:

- Tuner: Accessory Tongs
- Mode/Level: Desired Level
- Duration: 1-2 min
- Pressure: 1
- Location: Center of Eyebrows

Place tongs at the center of the eyebrows



Step 3:

- Tuner: Accessory Tongs
- Mode/Level: Desired Level
- Duration: 1-2 min
- Pressure: 2
- Location: Between Eyes


Place tongs close together between the eyes at the top of the nose with light pressure towards your forehead



Step 4:

- Tuner: Accessory Tongs
- Mode/Level: Green
- Duration: 1-2 min
- Pressure: 1
- Location: Under Jaw

Place tongs under the jaw holding the Tuner with both hands



Step 5:

- Tuner: Accessory Tongs
- Mode/Level: Desired Level
- Duration: 1-2 min
- Pressure: 1
- Location: Under ears


Place tongs along the jaw line with the ends of the tongs under your ears holding the Tuner with both hands



Step 6:

- Tuner: Accessory Tongs
- Mode/Level: Desired Level
- Duration: 7-8x
- Pressure: 2
- Location: In Mouth

Place tongs inside your mouth with the tongs stretched outside the teeth. Angle up and down and then slowly open/close your mouth 7-8 times



Step 7:

- Tuner: Accessory Clear Tip
- Mode/Level: Higher level
- Duration: 1-2 min per area
- Pressure: 2
- Location: Areas with Remaining Tension

Remove tongs, install clear tip, place on any remaining areas with tension

Loss of Taste and Smell

This Rezzipe supports your body's natural healing response for loss of Taste and Smell.



Step 1:

- Tuner: Accessory Tongs
- Mode/Level: Desired Level
- Duration: 1-2 min
- Pressure: 1
- Location: Center of Eyebrows

Place tongs at the center of the eyebrows



Step 2:

- Tuner: Accessory Tongs
- Mode/Level: Desired Level
- Duration: 1-2 min
- Pressure: 2
- Location: Between Eyes

Place tongs close together between the eyes at the top of the nose with light pressure towards your forehead



Step 3:

- Tuner: Accessory Tongs
- Mode/Level: Desired Level
- Duration: 1-2 min
- Pressure: 1
- Location: On Cheekbones

Place tongs on your cheekbones



Step 4:

- Tuner: Accessory Tongs
- Mode/Level: Green
- Duration: 1-2 min
- Pressure: 1
- Location: Under Jaw

Place tongs under the jaw holding the Tuner with both hands



Step 5:

- Tuner: Accessory Tongs
- Mode/Level: Desired Level
- Duration: 7-8x
- Pressure: 2
- Location: In Mouth

Place a drop of peppermint or lemon or other spice on your tongue, then place tongs inside your mouth with the tongs stretched outside the teeth. Angle up and down and open and close your mouth 7-8 times

Try different flavors, spices, food grade oils 1-2 times daily with step 5

Shoulder & Chest

This Rezzipe supports your body's natural healing response for shoulder and chest pain.



Step 1:

- Optional: Use Neck Pillow
- Tuner: Accessory Tongs
- Mode/Level: Blue
- Duration: 5 min
- Pressure: 1
- Location: Back of neck

Attach tongs to Tuner and place tongs on both sides of the neck with Tuner resting against the spine slowly rotating head side to side keeping pressure against the device.



Step 2:

- Tuner: Wing
- Mode/Level: Higher Level
- Duration: 1-2 min
- Pressure: 3
- Location: Arms

Gently scrape your arm with the wings starting from the fingertips all the way to the shoulder, repeat on other arm



Step 3:

- Tuner: Accessory Tongs
- Mode/Level: Higher Level
- Duration: 1-2 min
- Pressure: 3
- Location: Chest

Place Tuner over chest, with tongs around shoulder, repeat on other side of chest



Step 4:

- Tuner: Accessory Tongs & Armband
- Mode/Level: Higher Level
- Duration: 1-2 min (or much longer)
- Pressure: 3
- Location: Shoulder

Strap on to desired shoulder with armband



Step 5:

- Optional: Use Neck Pillow
- Tuner: Accessory Tongs
- Mode/Level: Higher Level
- Duration: 1-2 min
- Pressure: 1
- Location: Armpit

Place Tuner under armpit and gently hold arm to body



Step 6:

- Tuner: Accessory Clear Tip
- Mode/Level: Higher Level
- Duration: 1-2 min
- Pressure: 2 -3
- Location: Armpit

With hand on top of head, press Tuner Tip against the back side of the armpit for 1-2 minutes followed by underneath the front muscle of the armpit with upwards pressure for 1-2 minutes



Step 7:

- Tuner: Accessory Clear Tip
- Mode/Level: Higher Level
- Duration: 1-2 min
- Pressure: 2-3
- Location: Armpit

With arm at side, press Tuner Tip underneath the front muscle of the armpit with upwards pressure for 1-2 minutes

Elbows & Hands

This Rezzipe supports your body's natural healing response for elbow and hand pain.



Step 1:

- Optional: Use Neck Pillow
- Tuner: Accessory Tongs
- Mode/Level Blue
- Duration 5 min
- Pressure 1
- Location Back of neck

Attach tongs to Tuner and place tongs on both sides of the neck with Tuner resting against the spine slowly rotating head side to side, keeping pressure against the device.



Step 2:

- Tuner: Wing
- Mode/Level Higher Level
- Duration 1-2 min
- Pressure 3
- Location Arms

Gently scrape your arm with the wings starting from the fingertips all the way to the shoulder, repeat on other arm if desired



Step 3:

- Tuner: Accessory Tongs & Armband
- Mode/Level: Blue
- Duration: 5 min
- Pressure: 1
- Location: Forearm

Strap Tuner on forearm with armband and with tongs around elbow, repeat on other forearm if desired



Step 4:

- Tuner: Accessory Tongs & Armband
- Mode/Level: Blue
- Duration: 5 min
- Pressure: 1
- Location: Wrist

Strap Tuner on wrist with armband and with tongs around hand, repeat on other wrist if desired

Stomach

This Rezzipe supports your body's natural healing response for stomach and gastrointestinal pain.




Step 1:

- Tuner: Accessory Blue Brain
- Mode/Level: 1x Red, 4x Green
- Duration: 10 min
- Pressure: 2
- Location: Belly Button

Install blue brain, place on belly button for 1 red cycle and 4 green cycles

Step 2:

- Tuner: Accessory Clear Tip
- Mode/Level: High
- Duration: 15 Seconds - repeat 7 times
- Pressure: 2
- Location: Right Hip - Trace over the large intestine slowly



Remove blue brain, install clear tip, turn on high, begin at right hip go to rib cage, pull slowly as you cross to left side of rib cage. Drop slowly to left top of pelvis, go back toward the right side, stopping for 10 seconds on over the bladder (a couple inches below the belly button)

Back

This Rezzipe supports your body's natural healing response for shoulder and chest pain.



Step 1:

- Optional: Use Neck Pillow
- Tuner: Accessory Tongs
- Mode/Level: Blue
- Duration: 5 min
- Pressure: 1
- Location: Back of Neck

Attach tongs to Tuner and place tongs on both sides of the neck with Tuner resting against the spine slowly rotating head side to side, keeping pressure against the device



Step 2:

- Tuner: Wing
- Mode/Level: Higher Level
- Duration: 1-2 min
- Pressure: 3
- Location: Back

Scrape back starting from low back going up towards your neck



Step 3:

- Tuner: Accessory Not Applicable
- Mode/Level: Higher Level or green pattern
- Duration: 2-10 min
- Pressure: 1
- Location: Small of Back



Place Tuner in small of your back and lay down on it

Step 4:

- Tuner: Runners Pack
- Mode/Level: Higher Level or green
- Duration: 10-30+ min
- Pressure: 2
- Location: Back


Place Tuner inside the Runners pack and place the mesh against the low back



Step 5:

- Tuner: Runners Pack
- Mode/Level: Higher Level or Green
- Duration: 10x
- Pressure: 2
- Location: Sitting

With Tuner on the back straighten the leg, flex and extend foot for 15 seconds – repeat 10 times on each leg.



Step 6:

- Tuner: Accessory Not Applicable
- Mode/Level: Gentle Settings Low Level
- Duration: 10 min
- Pressure: 2
- Location: Sitting on the Tuner

Wrap Tuner in a heavy towel or blanket and sit on the Tuner like sitting on a saddle. Low level, hum to the device.

Legs

This Rezzipe supports your body's natural healing response for leg pain.



Step 1:

- Tuner: Wing
- Mode/Level: Higher Level
- Duration: 1-2 min
- Pressure: 3
- Location: Foot - Toes

Scrape foot and between toes (every bit of skin) with firm pressure with Tuner Wing



Step 2:

- Tuner: Wing
- Mode/Level: High
- Duration: 1-2 min
- Pressure: 3-4
- Location: Leg

Scrape all muscles in leg with firm pressure with Tuner Wing



Step 3:

- Tuner: Accessory Tongs
- Mode/Level: Higher Level
- Duration: 6-7x
- Pressure: 1
- Location: Hip

Install tongs, place on hip, with leg extended flex and relax leg/foot 6-7 times. Repeat steps (1-3) on other leg if needed



Step 4:

- Tuner: Accessory Tongs
- Mode/Level: Green
- Duration: 2 min
- Pressure: 1
- Location: Feet

Place Tuner on the ground with instep of feet on tongs and big toe resting on wings

Knees

This Rezzipe supports your body's natural healing response for knee pain.



Step 1:

- Tuner: Wing
- Mode/Level: High
- Duration: 1-2 min
- Pressure: 3
- Location: Legs

Scrape all muscles in lower/upper leg beginning at ankle and go up towards your hip with firm pressure as tolerated using the Tuner Wing, repeat on other leg if desired



Step 2:

- Tuner: Accessory Tongs
- Mode/Level: Higher Level
- Duration: 6-7x
- Pressure: 2
- Location: Hip

Install tongs, place on hip, with leg extended flex foot 6-7 times, repeat on other hip if desired



Step 3:

- Tuner: Accessory Tongs & Pillow
- Mode/Level: High
- Duration: 6-7x
- Pressure: 1
- Location: Under Kneet


Spread tongs about 1 inch apart, place under knee, flex foot 6-7x, repeat on other knee if desired



Step 4:

- Tuner: Accessory Tongs & Pillow
- Mode/Level: High
- Duration: 6-7x
- Pressure: 1
- Location: Under Knee

Completely spread tongs apart to contact sides of knee, place on pillow under knee, flex foot 6-7x, repeat on other knee



Step 5:

- Tuner: Accessory Tongs
- Mode/Level: High
- Duration: 6-7x
- Pressure: 1
- Location: On Knee

Place Tuner on top of knee, flex foot 6-7x, repeat on other knee if needed

Feet & Ankle

This Rezzipe supports your body's natural healing response for feet and ankle pain.



Step 1:

- Tuner: Wing
- Mode/Level: Higher Level
- Duration: 1-2 min
- Pressure: 3
- Location: Foot


Scrape foot and between toes with firm pressure with Tuner Wing



Step 2:

- Tuner: Wing
- Mode/Level: High
- Duration: 1-2 min
- Pressure: 3
- Location: Leg arm

Scrape beginning at ankle and go up towards your hip with firm pressure using the Tuner Wing



Step 3:

- Tuner: Accessory Tongs
- Mode/Level: Higher Level
- Duration: 6-7x
- Pressure: 1
- Location: Hip

Install tongs, place on hip, with leg extended flex foot 6-7 times, repeat on other hip



Step 4:

- Tuner: Accessory Tongs & Strap
- Mode/Level: Green
- Duration: 5-10 min
- Pressure: 1
- Location: Calf - Walking

Strap Tuner to calf with black velcro strap and with tongs around ankle, repeat on other calf



Step 5:

- Tuner: Accessory Tongs
- Mode/Level: Green
- Duration: 2 min
- Pressure: 1
- Location: Feet

Place Tuner on the ground with instep of feet on tongs and big toe resting on wings

Whole Body De-stress

This Rezzipe supports your body's natural healing response for general stress throughout the body.



Step 1:

- Tuner: Accessory Tongs
- Mode/Level: Blue
- Duration: 5 min
- Pressure: 1
- Location: Back of neck

Attach tongs to Tuner and place tongs on both sides of the neck with Tuner resting against the spine slowly rotating head side to side, keeping pressure against the device.



Step 2:

- Tuner: Accessory Tongs
- Mode/Level: Green
- Duration: 1-2 min
- Pressure: 1
- Location: Center of Eyebrows

Place tongs at the center of the eyebrows



Step 3:

- Tuner: Accessory Tongs
- Mode/Level: Desired Level
- Duration: 1-2 min
- Pressure: 2
- Location: Between Eyes

Place tongs close together between the eyes at the top of the nose with light pressure towards your forehead



Step 4:

- Tuner: Accessory Tongs
- Mode/Level: Desired Level
- Duration: 7-8x
- Pressure: 1
- Location: Back and Forth Over Head

Run the tongs over the top of your head back and forth 7-8 times



Step 5:

- Tuner: Accessory Tongs
- Mode/Level: Desired Level
- Duration: 7-8x
- Pressure: 2
- Location: In Mouth

Place tongs inside your mouth with the tongs stretched outside the teeth. Angle up and down repeatedly and then open and close your mouth 7-8 times



Step 6:

- Tuner: Accessory Tongs
- Mode/Level: Desired Level or Blue
- Duration: 1-2 min
- Pressure: 1
- Location: Stomach

Place Tuner on stomach with tongs at the base of your stomach



Step 7:

- Tuner: Accessory Tongs
- Mode/Level: Blue
- Duration: 5-10 min
- Pressure: 1
- Location: Feet

Place Tuner on the ground with balls of feet on tongs and big toe resting on wings. While humming, think of things for which you are profoundly grateful for.

For Sensitive Individuals

CRPS - POTS - Lyme - Fibromyalgia - MCAS - Trigeminal Neuralgia + others...

This Rezzipe supports your body's natural healing response by helping your body slowly get used to the vibrations from the Tuner.

Note that this process may take multiple sessions depending on your current conditions. As each step becomes comfortable, continue to next step.

Avoid areas of CRPS (Complex Regional Pain Syndrome) or excessive pain. Stick to areas that are not too sensitive until the nervous system calms down, and that can take weeks for sensitive people.



Step 1:

- Tuner Attachment: None
- Mode/Level: Lowest Level in Gentle Settings
- Duration: 1 min
- Pressure: Not Applicable
- Location: As Near as Tolerable (a level that will not cause an increase in pain) - wrap in a towel or blanket if needed.

Place Tuner on the ground with instep of feet on tongs and big toe on wings. Once you are able to tolerate the vibration near you, continue to the next step (If sensitive on feet, skip to step 3)



Step 2:

- Tuner: Accessory Tongs
- Mode/Level: Lowest level in Gentle Settings
- Duration: 1 min
- Pressure: Minimal
- Location: Feet or Foot (Depending on Sensitivity)

Place Tuner on the ground with instep of feet on tongs and big toe on wings, once comfortable, incrementally increase the level (if sensitive on neck, skip to step 4)



Step 3:

- Optional: use neck pillow
- Tuner: Accessory Tongs
- Mode/Level: Lowest
- Duration: 1 min
- Pressure: 1
- Location: Neck

Place Tuner on the back of your neck, once comfortable, slowly increase the level - perhaps a few days later...



Step 4:

- Tuner: Accessory Tongs
- Mode/Level: Lowest
- Duration: 1 min
- Pressure: 1
- Location: Center of Eyebrows

On the lowest level, place on the center of your eyebrows, incrementally increase the level to a comfortable level.



Step 5:

- Tuner: Accessory Tongs
- Mode/Level: Lowest
- Duration: 7-8x
- Pressure: 1
- Location: In Mouth

On the lowest level, place inside mouth, once comfortable, angle up and down and open and close mouth 7-8 times

Step 6:

- Tuner: Accessory Tongs
- Mode/Level: Lowest
- Duration: 1 min
- Pressure: 1
- Location: Skipped Area

Once your body is comfortable at all other areas, you may begin the step (2-4) that was skipped when you were too sensitive previously

Primitive Reflex Integration using a Tuner

1: Moro (startle) Reflex - Sudden sound or movement



A. Begins to develop 9 weeks in Utero - inhibited by 2-4 months of life - by about 1 year the vagal system develops and inhibits Moro

B. Retained? You may see Hypersensitive/reactive, Poor impulse control, motion sickness, poor coordination, visual perception difficulties, light/sound sensitive, difficulty with change

C. Rezzi Solution: Starfish exercise, Jumping jacks

D. Tuner 10 x scraping tailbone to head (best on all fours)

- Medium-High intensity
- Pressure 2
- Duration 1 minute
- Place Tuner on stomach for 1-2 minutes







2: Palmar Reflex – Pressure in the palm causes fingers to close

- A.** Begins 11 weeks in utero – inhibited by 2-3 months
- B.** Retained? Poor dexterity (Babkin response) – poor manual (hand) skills – speech difficulties

C. Rezzi Solution:

- Scrape from fingertips to neck.
- Use clear tip or wing corner on Tuner to make “x” patterns across palm; on high setting, 10x per hand



3: Babinski Reflex – Toes curl = plantar response – Big toe extends and toes flare = classic Babinski

- A.** Present at birth, fading by 1 year
- B.** Retained? Toe walking, poor balance, tripping/falling
- C.** Rezzi Solution:

- Scraping the body with Tuner on as high as tolerated from toes to head, 66% of time on side of weakness, 33% of time on strong side. 20 minutes daily or 3x/ week at a minimum
- Use Tuner wing tip, device as strong a level as tolerated, to run from heel along outside of foot along the small toes and then toward big toe as shown in test position. Repeat 12 times per side.



4: Spinal Galant Reflex

- A.** Begins 20 weeks in utero - inhibited by 3-9 months
- B.** Retained? Fidgeting - bedwetting - poor concentration - poor short-term memory - 1 or 2 sided postural issues



C. Rezzi Solution:

- Scraping the spine and sides of the spine with Tuner on as high as tolerated from tailbone to top of the neck, 66% of time on side of weakness, 33% of time on strong side. 5 minutes daily




5: Rooting Reflex - Searching, Sucking, & Swallowing Reflex

- A.** Begins 24-28 weeks in utero - inhibited by 4-6 month
- B.** Retained? Hypersensitive lips or mouth - tongue thrust - speech difficulties - poor dexterity (Babkin response)



- C.** Rezzi Solution: Run Tuner on medium to high as tolerated, from ear to corner of mouth & across the lips 10 times per side
- D.** Intraoral Tuner on tongue



6: Symmetrical Tonic Neck Reflex (STNR)

- A.** Begins 6-9 months of life Head goes up, arms extend.
Inhibited by 9-11 months
- B.** Retained? Poor posture, slump sitting, ape walk, toe walking, poor coordination, difficulty sitting still, difficulty concentrating



- C.** Rezzi Solution: Scraping the body with Tuner on as high as tolerated from toes to head, 66% of time on side of weakness, 33% of time on strong side. 20 minutes daily or 3x/week at a minimum



7: Asymmetrical Tonic Neck Reflex (ATNR)

- A.** Begins in utero - inhibited by 3-4 months of age
- B.** Retained? Poor eye-hand coordination, poor handwriting, difficulty reading and tracking



- C.** Rezzi Solution: Scraping the body with Tuner on as high as tolerated from toes to head, 66% side of weakness, 33% of time on strong side. 20 min 3x/week

8: Tonic Labyrinthine Reflex (TLR) 2 parts - flex neck; arms together - Extend neck; arms apart

A. Flexion begins in utero - extension begins at birth

B. Retained? Poor posture, weak muscle tone, poor visual/spatial perception, poor sequencing, dislike of sports, poor organization, poor sense of time/rhythm, fatigue easily



C. Rezzi Solution: Scraping the body with Tuner on as high as tolerated from toes to head, 66% of time on side of weakness, 33% of time on strong side. 20 minutes daily or 3x/week at a minimum

D. Use Tuner in Runners Pack for 30-60 minutes each day on blue or green algorithm with pack oriented toward side of weakness or on stomach

9: Landau Reflex (Superman) Inhibits TLR - engages extensors

A. Emerges 4-5 months of life, inhibited ~ 3.5 years

B. Retained? Affects balance and tone through movements - awkward or stiff running - difficulty hopping, skipping, jumping



C. Rezzi Solution: Scraping the body with Tuner on as high as tolerated from toes to head (avoiding sensitive areas), 66% of time on side of weakness, 33% of time on strong side. 20 minutes daily or 3x/week at a minimum

- Use Tuner in Runners Pack for 30-60 minutes each day on blue or green algorithm with pack oriented toward side of weakness or on stomach



Principles of the Brain and the Stress response

The brain is wired to create “superhighways” which create patterns of behavior and automatic reaction.

- Good and bad experiences make us more efficient.
- Strong beliefs and patterns of behavior - wired neural pathways are wrapped with myelin sheaths to block out interference, creating “super-highways.”

Our Stress Down-Regulation program “bombs” these old super-highways, putting ruts and holes in them so it will be easier to create new and powerful neural pathways free from the intrusive trauma patterns. Examples:

- Porn addict who struggles with the chemicals and wiring in his/her brain and wants to view the porn (the old highway) needs to be able to more easily access the new highway of choosing a better activity and avoiding the old patterns of choice.
- Alcohol addict who recalls bad experiences on a regular basis that wants to alleviate the tension caused by bad memories. The brain seeking relief will find something to fill the void...
- We wonder if it is possible to create new “brain muscle memory” to heal and down-regulate the nervous system.

Preparing to undergo the Stress Down-Regulation process.

- Identify what you want to work on – a feeling, a memory, an experience, a fear, etc
- Journaling can help

The Stress Down-Regulation Technique

1. Start by thinking of stress or worry of the past, present or future. Note current pain/stress level from Low (1) to High (10).

2. Strap the Rezzimax Tuner to the non-dominant lower extremity. Turn on an algorithm pattern.

- Billions of bits of carefully calibrated information will be sent to the brain at the same time you are thinking about a Stressful event or a Worry or Fear you may have. This will help break up painful patterns. The resonance helps by sending the brain data from hundreds of thousands of sensors in joints, ligaments, muscles, and cells throughout the body. The balanced resonant vibration simply accelerates the process.

3. How to do it



1. Place the knot in the middle of a strap or rope in the space between the door and the door frame near the hinges and shut the door tightly. Place it at the middle hinge level to anchor the strap tightly. Alternatively hug an exercise ball, pillow, chair or your legs
2. Put your dominant (dominant side is generally the side that you use to write) foot back, the non-dominant foot in front, so you are balanced, and then squat down slightly and pull firmly against the strap to engage the core muscles. Keep this tension against the strap while performing the exercise.



3. Turn your head to the dominant side. Turning your head to the dominant side keeps the body from turning on the stress chemicals. (Imagine the face of someone experiencing some type of stressful event, their jaws are clenched, neck muscles are tense and tight, and they may even have bulging veins in their temples)! Turning the head to the side weakens the strength of one of those neck muscles (the sternocleidomastoid on the dominant side). Turning our heads to the dominant side weakens the physiological stress in the neck area to allow the calming chemicals to flow more easily.

4. Stress Mode

- Breathe normally for 20-30 seconds while thinking about all the details you can recall about a particular stressful event, and become aware of any senses you can recall, such as sights, sounds, smells, feeling, time of day and year, who you were with, etc. The goal is to briefly bring up the stress chemicals from the sympathetic nervous system caused during the particular events.

Take a big breath of air and then let the air out in 3 distinct chunks (out, hold, out, hold, out the rest of the way.) The hold is brief and then continue breathing out. This causes breakup of trauma chemicals.

- When the breath is fully out of the lungs, hold at this state for as long as you can (while pulling for all you can!)
- It is at this stage, when the breath is all out of your lungs, that you need to tell yourself to stop thinking about the stressful event and begin thinking about the things that are most important to you. Thinking about the things for which you are profoundly grateful will provide your brain with material with which it can build new connections as you break up the chemicals of prior stressful events. Think of as many “gratitude and happy thoughts” as you can until you cannot wait any longer to come up for air. When you must finally breathe in, relax completely and breathe deeply for several seconds before performing the exercise again.
- Performing this exercise creates a vacuum state in the brain where you can shut off the chemicals of the stressful event (and anxiety and depression and allergy and...) chemicals and turn on the growth and healing chemicals. This heals wounds created by past experiences and future fears. You will be confusing the chemicals associated with the stressful event and breaking them up! Remember, wait until you absolutely need to, to come up for air!
- You can close your eyes during this whole process if you wish.

5. Repeat the Process 7 Times; perform daily for 6 weeks for best results. Note pain/stress level on a scale of 1-10 at the end of the exercise.

6. Triggering the Healing Chemicals in a trained brain

Don't let the fear of past, present or future stressful events, anxiety, stress, or allergies keep you from being who you need to be!

Appendix

Note: The techniques included may be helpful for the following conditions.

Page numbers are where the Rezzipe starts. Complete the whole Rezzipe for best results.

Addictions	32	Fibromyalgia	22
Adhesive Capsulitis	15	See individual body part sections as needed	
Allergies	34	Focus	16, 22
Anxiety	34	Frozen Shoulder	10
Arthritis	22	Gastrointestinal Issues	14
See individual body part sections as needed		• Gastrointestinal Upset	
Attention-Deficit Disorder (ADD)	22	• Gastroparesis	
Attention-Deficit/Hyperactivity Disorder (ADHD)	22	• Internal Organs	
Autism	22, 26-30	• Irritable Bowel Syndrome (IBS)	
Autonomic Nervous System Dysfunction (ANS)	32, 34	Hammer Toe	20
Balance	27, 31	Headache	4
Bedwetting	28	Heart Issues	10, 22
Blood Flow	12, 22	Hypotonia	22
Brain Fog	4	See individual body part sections as needed	
Brain Injuries	4, 22	Loss of Senses	
Broken Bones	22	• Hearing	4
See individual body part sections as needed		• Smell	8
Cancer Pain	22, 24	• Taste	8
Carpal Tunnel	12	Lymphatic Flow/Drainage	10
Chronic Fatigue Syndrome	4, 31	Migraine	4
Complex Regional Pain Syndrome (CRPS)	24	Muscle Strain	
Concussion	4	See individual body part sections as needed	
Core Strengthening	14, 16	Muscle Tension	22
Crohn's Disease	14	See individual body part sections as needed	
Craniofacial	4, 6	Neuralgia	24
Dead Arm	10, 12	See individual body part sections as needed	
Dental	5, 7, 9	Neuropathy	24
Depression	34	See individual body part sections as needed	
Diaphragmatic Training	10, 14	Numbness/Tingling	22, 24
Diverticulitis	14	See individual body part sections as needed	
Earache	4	Pain (Chronic and Acute)	24
Fascia Release	26	See individual body part sections as needed	
See individual body part sections as needed		Paralysis	22
		See individual body part sections as needed	

Parkinson's Disease	22
Pectoralis Minor	10
Pelvic Unlocking	14 , 16
Plantar Fasciitis	20
Post-Traumatic Stress Disorder (PTSD)	33
Postural Orthostatic Tachycardia Syndrome (POTS)	22
Pre- and Post-Surgical Rehab	22
See individual body part sections as needed	
Primitive Reflexes	
• Asymmetrical Tonic Neck Reflex (ATNR)	29
• Babinski	27
• Landau	31
• Moro (Startle)	26
• Palmar	27
• Rooting	28
• Spinal Galant	28
• Symmetrical Tonic Neck Reflex (STNR)	29
• Tonic Labyrinthine Reflex (TLR)	30
Proprioception	22
Reading Difficulties	22
Restless Leg Syndrome	20
Rotator Cuff Tear	10
Scar Tissue	
See individual body part sections as needed	
Sciatica	15 , 17
Seizures	4 , 22
Sensitive	24
Sensory Processing Disorder	22
Sinus Congestion	4
Sleep Disorders	22
Soreness/Stiffness	
See individual body part sections as needed	
Sports Performance Enhancement	30
Stress	22 , 33
Subscapularis	10
Swelling	22
See individual body part sections as needed	
Teething	4 , 6
Temporomandibular Disorders (TMD)	6
Temporomandibular Joint Dysfunction (TMJ)	6
Tennis Elbow	12
Tinnitus	4
Trauma	32 , 34
See individual body part sections as needed	
Tremors	22
Trigger Finger	12
Vocal Chord Tone	4
Whiplash	
See individual body part sections as needed	

