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Lumières continues

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Continuous Lights

When you press:

1 time: 650 and 850 nanometers (nm), red light and infrared

2 times: 660nm red light

3 times: 850nm or 940nm infrared, will be almost invisible to the eye

The green button indicates that your product is on

Pulse Mode:

When you press the button for 3 seconds, it turns blue. The same powers as those of the continuous light mode will be available, but the treatment will be in pulsation.

Ideal for blood circulation activation, lymphatic drainage, reduction of tissue warming, reduction of inflammation, and much more.

The blue button indicates that your product is on

Recharging your Goldride product:

When the light flashes red, your product is charging. The button will stop flashing and remain red continuously when your product is full.

The recharge time may vary, depending on the product. You can still use it even if it is not 100% full.

*Use the provided charging cable. *

Usage:

Use your Goldride product 1 to 2 times per day for at least one week, 1 treatment every 1 hour. *Built-in timer*

Reduce treatments to once a day if desired or maintain 2 treatments per day.

Do not do more than 1 treatment in a row in the same location. Allow 60 minutes before doing a 2nd treatment.

Benefits:

Temporary relief of muscle and joint pain, muscle relaxation, relaxation of muscles and tissues, pain relief treatment for arthritis and osteoarthritis, natural anti-inflammatory, blood activation, production of endorphins, and much more!

Wavelengths in the spectrum of 630 and 940 nm provide optimal biotic response. Super-bright LEDs (light-emitting diodes) are concentrated, and safe light sources used to produce specific light powers that provide the therapeutic benefits of light therapy without undesirable side effects. Red light therapy offers athletes significant benefits in rehabilitation after injury and performance enhancement in all sports.

Cleaning your Goldride product:

To clean your product, please:

1) Use a damp cloth (baby wipes, Lysol wipes, etc.)

2) Clean the desired area.

3) Let it dry before storing the product.

*If used on multiple animals, please disinfect the product each time. *