



2025/26




schuchmann®





**„My momo tricycle. and I have long been
an unbeatable team.“**

Sophia has grown up with her **momo tricycle**.



10
–
33

Sitting + Positioning.

34
–
49

Transport + Travelling.

50
–
67

Walking + Cycling.

68
–
85

Standing + Transfer.

86
–
99

Bathroom + Toilet.

100
–
107

General guidelines.

Whether sitting, driving, walking or standing – our aids accompany and support you in everyday life. They help you get out to meet your friends and promote your development and health. Every aid can be individually adapted to your body, and also to suit your tastes.

Contents.





**Hello!**

We are Miriam and Torsten Schuchmann. Our family-run company is now in its second generation, our company has over 40 years of experience in the design, manufacture and provision of childrens' rehabilitation equipment.

It is a matter close to our hearts to support children with individual requirements as optimally as possible – using products which are just as individual as you are.

What inspires us.

The reason we have focused our efforts on childrens' rehabilitation equipment has its origins in our own family. We have first-hand experience and completely understand the challenges associated with a severely disabled person with multiple disabilities.

Miriam's sister was born with severe cerebral palsy. When she was young the development of rehabilitation equipment was at a very early stage. Her parents recognised what could be done to them to improve her participation in everyday life, and a business concept was born.

What we inspire.

Our aim is to support children with disabilities so that they can become more mobile and independent. At the same time, our task has more far-reaching implications. Our 180-strong team work closely with specialist dealers, doctors and therapists so that we can truly understand the requirements of the children and their families. This is an extremely moving process for all involved.

Our products represent more than just products, they become companions in life. We want to get the market moving – with real innovations in rehabilitation equipment which fulfil the maximum demands on functionality, design and safety.

If you have any questions or suggestions on general topics or on our products, we – and of course the entire Schuchmann Team – are happy to help at any time.

Sitting + Positioning.



smilla.
14–17



madita-fun mini.
18–21



madita-fun.
22–25



new

maxi.
26–29



new

tina.
30–33

Transport + Travelling.



carrot 3.
38–41



new

marie integral.
42–45



new

charlie integral.
46–49

Walking + Cycling.



malte.
54–57



malte outdoor.
58–59



marcy.
60–63



new

momo tricycle.
64–67

Standing + Transfer.



till.
72–77



todd.
78–81



tim.
82–85

Bathroom + Toilet.



boris.
90–93



basti.
94–97



bibi.
98–99

General guidelines.



**Cleaning +
disinfection.**
102–103



**Fabric collection +
Name embroidery.**
104–107

Contents.



Key.

Briefly and concisely explained.



Indoor



Outdoor



Wet area



Car



Max. load

Angular adjustment
of central column

Back angle



Tilt in space



Seat height



Seat angle



Seat depth



Seat width



Elbow height



Handle height

Inside leg
length

Body size



Access point

GMFCS
classification

You can choose your own fabric for your aid from our extensive and colourful range
see pages 94 to 95.



Harness possibilities
We offer many different harness possibilities to make you feel safe, and to allow you to use your aid whilst positioned straight and upright (see price list).



Your seating aid belongs to you, and everyone should see that.
We can embroider your name onto the back unit. You can find more information on page 95.



Diverse accessories program for optimum support & individual equipment. You can find more information in the price list.

Our icons quickly show you how and where you can use our aids.

You can see at first glance if our products are suitable for inside or outside use, and how they can be adapted for individual body measurements and requirements.

You have the choice:

of many great fabrics, straps, accessories and embroideries. So that your aid is as individual as you are!

GMFCS.

Using this, finding your daily companion is child's play.

GMFCS stands for a classification of the gross motor skills of children and young adults. Doctors and therapists use this system to classify the degree of motor impairments in children suffering from cerebral palsy. Here five stages have been defined from Level 1 to Level 5. Each level defines how and whether children are able to move independently or whether they are dependent on aids and mobility aids.

GMFCS between the 6th and 12th birthday:

GMFCS Level I



The children can walk freely at home, in school, outdoors and in the surrounding area. The children can climb over and get down from kerbstones without additional physical aids, and climb stairs without holding on to railings. Running and hopping is possible, but speed, balance and coordination are limited. The children can participate in physical activities and sports depending on their personal requirements and the environmental conditions.

GMFCS Level II



In most situations, the children can walk freely, but experience difficulties over long distances and on uneven ground, in case of slopes, in crowds of people, on narrow paths or when carrying objects. The children can climb up and down stairs by holding onto railings, or can achieve this with physical assistance if no railings are available. Outdoors and in their surroundings, they are able to walk with physical support or using a walking aid, or they use a wheelchair over longer distances. Running and jumping are at best restricted. The restrictions in the execution of physical-motor tasks can make adaptations necessary which permit children to participate in physical activities and sports.

GMFCS Level III



The children are usually able to walk indoors with a walking aid. When the children are seated, they may require a seat strap to straighten the pelvis and for balance. When transitioning from sitting to standing and from the floor to standing, physical assistance is required from a person or a supporting surface. Over longer distances, the children use a form of rolling mobility aid. Children may be able to go up and down stairs by holding onto railings under supervision or with physical support. The limitations of walking may make adaptations necessary which first permit participation in physical activities and sports, including an active wheelchair or an electric wheelchair.

GMFCS

We have assigned the appropriate GMFCS level to each of our products. This should assist in selecting the right aid for your child's requirements. The classification serves merely as a guide, and may deviate from child to child.

GMFCS Level IV



The children use methods of movement which usually involve physical support or electrical mobility aids. The children require adapted seat devices for torso and pelvic control; physical support is also required for most transfers. At home, the children move along the floor (rolling, creeping or crawling), walk for short distances with physical assistance or use electrical mobility aids. When the children are brought into a standing position, they can use a walker at home or in school. At school, outdoors and in the surrounding area, children are pushed in an active wheelchair or they use an electric wheelchair. Due to their restrictions in movement, adaptations are required which first permit physical activities or participation in sports, including physical support and/or electrical mobility.

GMFCS Level V



The children are transported in all situations in a wheelchair. The children are restricted in the possibilities for holding their head and torso upright against the force of gravity and in their arm and leg movements. Supporting technologies are used in order to improve straightening of the head, sitting, standing and/or mobility, but the existing limitations cannot be compensated for by the additional aids used. Transfers require the physical support of an adult. At home, the children are able to move short distances along the floor or are carried by an adult. The children may be able to move independently with an electric wheelchair through comprehensive adaptation of the sitting position and the controls. The restrictions in mobility make adaptations necessary for the children to participate in physical activities and sports, including physical support and the use of electrical mobility aids.


You can find out more information on the GMFCS classification between children's 12th and 18th birthdays on www.schuchmann.co.uk.



Sitting + Positioning.



Cooking together – tastes delicious!



Leah loves cooking with her dad.
Well positioned in her **madita-fun.**, the fresh fruit
and veg almost cuts itself. With a little help from
Leah! It's a good thing that **madita-fun.** can be
easily adjusted to a suitable working height.

Sitting + Positioning.



After school, Leah likes to help with the cooking. She is also well positioned in her **madita-fun** therapy chair for this. Her therapy chair helps her to adopt the best possible posture. Sitting in the right position is important so that Leah's bones and muscles can develop well.

Sitting is a basic posture for humans, one which holds even more significance for children with disabilities. A correct sitting position is hugely important in order to prevent deformations. Even for the youngest children we should make sure that the spine, pelvis and hips are not damaged by poor posture.

The correct seating aid, which can be individually tailored not only prevents misalignment it can also make sitting enjoyable for children. The optimum sitting position has positive effects on your back, digestion and breathing. With the correct support children can feel less fatigued and can play with friends and participate in everyday life. The correct seating aid can also relieve strain on parents, carers and teachers too.

Discover the world from a sitting position

In a sitting position even the littlest ones can start to explore their environment. You don't have to worry about them maintaining their posture as every seating aid we provide will give them just the right amount of support. Maximum freedom of movement can be achieved with **madita-fun**'s removable back unit. The freedom provided to the child in the upper back motivates them to sit up more and provides more space for activity and to explore what's around them. However, children should never remain sitting in the same position for prolonged periods of time for this reason our seating aids offer many options to change and vary the users posture during use.

It is important to us to tailor all seating aids to the requirements of children, parents and therapists. Therefore they are variable in use, easy to handle and offer very different functions to ensure that children feel comfortable. Not only in a sitting position, but also when reclining. For this reason, we also have a cradle system **tina**. for children who suffer from extreme contractions and sensitivity to pressure. In this mobile cradle, little ones feel comfortable, can relax and can also be brought into a sitting position for brief periods of time.

Sitting + Positioning.



smilla.
14-17



madita-fun mini.
22-25



madita-fun.
18-21



maxi.
26-29



tina.
30-33



GMFCS
Level I-III

kg
max. 60



(-) 3°-17,5°

(-) 5°-15°

16-40 cm

17-40 cm



I am your seating aid for everyday use.

I promote active sitting and support your development. The great thing is whatever size you are, I will grow with you! You can simply adjust my seat and footplate up or down so that you can sit comfortably at the table with your big brother or sister. I have been manufactured with loving attention to detail using wood, so I suit the other furniture in your home perfectly! With my many accessories, I can easily be adapted to suit your individual requirements.



You decide how I should look like.
Combine my frame in anthracite, natural beech or white stain with a selection of great colours for the covers!



Height adjustment made easy!

Simply press the button and the seat and footplate can be moved up or down.



Toolless angle adjustment!

The seat angle can be adjusted by turning the twist grip.

Push handle

for easy transfer.

Backrest

can be adjusted in height and angle.

Armrests

can be adjusted in height.

Seat plate

can be adjusted in height, depth and angle.

Base foot

with plastic sliders.

Headrest
supports the head control.

Removable pad cover
made from breathable material.

Thorax pelotte pads
facilitate a stable and upright sitting position.

Thigh guide
provides the necessary support.

Wooden frame
with integrated aluminium profiles.

Foot plate
can be adjusted in height, depth and angle.

Fast height adjustment

via the aluminium profiles
integrated in the frame.

**Frame colours**

- Beech natural
- Beech white
- Beech slate grey



		Size 0 nursery	Size 0	Size 1 nursery	Size 1	Size 2
Seat width		17 - 28 cm		20 - 35 cm		30 - 40 cm
Seat depth		16 - 26 cm		25 - 35 cm		30 - 40 cm
Max. seat height		42 cm	62 cm	42 cm	62 cm	
Back height*		25 cm		29 cm		34 cm
Back angle		(-) 3° - 17.5°				
Seat angle		(-) 5° - 15°				
Footrest quick height adjustment min.**	Standard	12 cm				
	with 50 mm dual castors	17 cm				
	with 75 mm dual castors	20 cm				
Height of the footrest at floor level min.***	Standard	2.5 cm				
	with 50 mm dual castors	7.5 cm				
	with 75 mm dual castors	10.5 cm				
Lower leg height with quick height-adjustable footrest min.		14 cm				
Lower leg height with footrest low to the ground min.		18.5 cm				
Total width		46 cm	47 cm	53 cm	54 cm	59 cm
Length base frame		63 cm	69 cm***	63 cm	69 cm***	69 cm***
Max load.		30 kg		45 kg		60 kg
Weight		9.2 kg	9.7 kg	9.6 kg	10.2 kg	12 kg
* Minimum back height measured from seat surface to upper edge of back						
** adjustable lower leg length max. = selected/required seat height - footrest height min.						
***long stand + 10 cm						



You can find more information on fabrics & embroidery on pages 104 - 107.

You can find more information on accessories & harnesses in the price list.



madita-fun mini.

Small, but impressive.

GMFCS
Level III-V

kg
max. 20



(-) 5° - 25°

(-) 5° - 42°

12 - 19 cm

25 - 62 cm

14 - 22 cm



I am the therapy chair for when you are young – I will give you security from the start.

I will place you on a level with your friends. My comfortable seat can be easily adjusted in height so that you can sit at any table and then be lowered back down to a height suitable for play. I will offer you just the right amount of support to help you grow and learn through your youngest years.





madita-fun mini.

Push handle

facilitates transfer.

Armrests

can be adapted in height and width.

Pelvic harness

provides the necessary support.

Headrest

is available in a three-dimensionally adjustable or height adjustable version.

Lateral pads

facilitate a stable and upright sitting position.

Seat

ensures optimum development of the hip joint.

Foot support

padded.

Seat height

can be adapted with support from a gas pressure spring.

Castors

facilitate transfer.

As individual as you are: in addition to a standard seat pad, you can select a pad with anatomic leg guidance or a dysplasia seat.



Ideal for
initial care.



	madita-fun mini.
Seat depth	12 - 19 cm
Seat width	14 - 22 cm
Height, short back	28 - 33 cm
Height, long back	30 - 35 cm
Seat height from floor	25 - 62 cm
Lower leg length	12 - 22 cm
Total width	52 cm
Length of base frame	58 cm
Seat angle	(-) 5° - 37° / 0° - 42°
Back angle	(-) 5° - 25°
Max. load	20 kg
Weight	16 kg



You can find more information on fabrics & embroidery
on pages 104 - 107.

You can find more information on accessories & harnesses
in the price List.



GMFCS
Level III-V



I am your integrative therapy chair.

Whether at the breakfast table at home or the play corner in school you can sit at the same level as your friends. I can be adapted to suit your body and your taste and I am available in lots of great colours! Did you know that my backrest can be removed? This helps to train your back muscles and helps you to sit up straight.





Push handle

for easy manoeuvrability.

Back unit

optionally removable, can be adjusted in angle and height.

Foot support

is available in three new versions:

- without lock
- with lock
- with progressive angle adjustment.

Castors

available in different sizes, facilitate movement.

Headrest

can be adjusted in height, depth and angle.

Lateral pads

allow a stable and upright sitting position.

Armrests

are individually adjustable.

Seat surface

is individually adaptable.

Base frame

compact and sturdily built.

In just one single step, my back unit can simply be removed. This gives you more freedom to move.

**Frame colours**

● Anthracite metallic

● Telegrey 4



My headrest with angle-adjustable side guide is flexible, adaptable, practical and offers great therapeutic benefits.

	Size 0	Size 0b	Size 1	Size 1b	Size 1c	Size 2	Size 2a
Seat depth*	15 – 25 cm	21 – 30 cm		26 – 36 cm	30 – 42 cm		
Seat width	16 – 26 cm		20 – 30 cm			28 – 38 cm	
Back height with extra short back	32 – 35 cm			35 – 43 cm			
Back height with short back	35 – 43 cm			42 – 52 cm			
Back height with long back	42 – 52 cm			53 – 63 cm			
Seat height				25 – 62 cm			
Seat height with electric height adjustment				32 – 62 cm			
Seat height with hydraulic height adjustment				27 – 63 cm			
Lower leg length	12** – 25 cm	12** – 25 cm	16** – 33 cm	19** – 36 cm	25 – 43 cm		
Total width	54 cm			64 cm			
Length of base frame	61 cm			75 cm			
Total length with combi subframe	76 cm			86 cm			
Tilt-in-space	(-)8° – 37°						
Back angle	(-)5° – 25°						
Max. load	30 kg			45 kg		50 kg	
Weight	19 kg		20 kg	22 kg		24 kg	
* reduced by 1cm in combination with adjustable pelvic guide		** the minimum lower leg length is achieved by repositioning the screws					

* reduced by 1cm in combination with adjustable pelvic guide

** the minimum lower leg length is achieved by repositioning the screws



You can find more information on fabrics & embroidery
on pages 104 – 107.

You can find more information on accessories & harnesses
in the price list.



Support beyond childhood.

GMFCS
Level III-V

kg
max. 125



(-)5°-25°

(-)20°-20°

30-55 cm

40-55 cm

28-48 cm




I am your grown up therapy chair.

I'm your daily companion for schoolwork and comfort. I ensure you can sit with perfect posture no matter how tall you get. With my discreet design, I blend in seamlessly with everyday chairs, so you never feel out of place.

I keep you positioned just right, allowing you to focus entirely on your tasks. My sturdy central column and extended wheelbase provide exceptional stability, making me strong enough to handle heavy use.

new





Unique design

discreet look which fits perfectly into a teenagers room or school environment.

Growth

sized to accommodate children up to approximately 195 cm and 125 kg.

Electric motor with hand control

for independent adjustment of height, seat and backrest angles.

Headrest

is upholstered and adjustable in height and angle.

Lateral pads

enable a stable and upright sitting position.

Armrests

are width, height and depth adjustable and can also be folded down.

Seat

anatomically shaped ensuring optimal pelvic positioning.

Base frame

is compact and robustly manufactured.

Castors

can also be fitted with an optional user brake for added independence.

Foot rest

can be adjusted in height and knee angle and is firmly attached to the seat unit.

High level of safety

thanks to solid and stable base frame.

Durable cover

in modern slate grey with excellent quality.



	Size 1	Size 2	Size 3
Seat depth	30 - 42 cm	38 - 48 cm	45 - 55 cm
Seat width	28 - 38 cm	36 - 48 cm	
Short back height	42 - 54 cm		
Long back height	53 - 62 cm		
Seat height	40 - 55 cm		
Lower leg length	27 - 45 cm	35 - 52 cm	
Total width	60.5 cm		
Base frame length	84.5 cm		
Seat angle	(-)20° - 20°		
Back angle	(-)5° - 25°		
Max. load	125 kg		
Weight	45 kg		

maxi. offers more independence thanks to electric adjustment of the height and depth as well as the angle of the seat and backrest.



You can find more information on fabrics
on pages 104 – 107.

You can find more information on accessories & harnesses
in the price List.



GMFCS
Level IV-V

kg
max. 100



130°–180°

0°–40°

115–190 cm

new



I am your modular cradle system.

I'll make sure you're always as comfortable as possible. For the most complex of cases my optional vacuum mattress relieves pressure points and pain, while my optional cushioning provides comfort for less demanding needs.

I'm designed for both indoor and outdoor use, giving you the freedom to explore. With my electronic height adjustment and customizable seat, back, and knee angles, you'll always find the perfect position – whether sitting or lying down. And thanks to my smooth-rolling castor wheels, I'm easy to push and maneuver, so you can take me anywhere with ease.







You can recline on me, but also sit comfortably too so that you can always have an overview of the world around you.



Frame colour



White-blue

Side protection
for maximum safety.



Electric pump
for the vacuum mattress.



Cradle subframe	
Height of cradle support from floor	53 - 79 cm
Adjustable back angle	130° - 180°
Seat angle	0° - 40°
Adjustable knee angle	110° - 180°
Base frame width	64 cm
Max. load	100 kg
Weight	70 kg




Cradle system	Size 1	Size 2	Size 3
Width	59 cm (with fall out protection + 12 cm)		
Length of head section	60 cm	75 cm	90 cm
Length seat section	35 - 50 cm		
Length leg section	45 cm	55 cm	65 cm
Total length	140 - 155 cm	165 - 180 cm	190 - 205 cm
Height	6 cm (with fall out protection + 9 cm)		



You can find more information on accessories & harnesses
in the price list.



Transport + Travelling.



Annelie loves going on trips with her family.
Today they are going to the playground together!
She is well protected and supported in her
marie integral.

Our other transport aids also keep you on the move
and get you safely to your destination.

A ride that's fun.



Transport + Travelling.



Finally, she can dig in the sand with her parents. That's what Annelie is looking forward to most today. She enjoys a relaxed ride to the playground in her marie integral car seat. Together, the family builds a great castle!
On the way back, Annelie takes a nap. And mom and dad are relaxed too, because they know their daughter is well protected and supported.

Mobility is essential for children with disabilities and their families. It not only enhances their quality of life but also helps them build social connections. Transportation aids like car seats, rehab buggies, and sub-frames ensure children can travel safely—whether it's to school, nursery, a friend's house, or even a trip to the countryside.

Our specially designed rehab car seats make it easier for parents and caregivers to transport children safely and securely. They meet all relevant regulations and have been independently tested for safety. The best way to travel is with the right harness and the innovative Side Impact Buddy. This smart mechanism allows the car seat to rotate away from the impact point in the event of a side collision.

But for us, it's not just about safety—it's also about the child's comfort and making everyday handling as simple as possible.

Sitting Correctly, Feeling Comfortable

Children should ride as comfortably as possible. That's why the backrest height, seat tilt – both forward and backward – and other adjustments can be tailored to meet their individual needs. These changes are quick and easy to make, usually in just a few steps.

Proper positioning is essential. A child sits well when their pelvis, upper thighs, and feet are securely supported. For good balance, the torso, shoulders, and head must be correctly aligned. The more freedom a child has to move their head, arms, and hands, the more actively they can engage with the world. These essential adjustments can be made using a range of accessories and built-in customisation options. And when standard solutions aren't enough, additional aids are available to accommodate special requirements.

Transport + Travelling.



carrot 3.
38–41



marie integral.
42–45



charlie integral.
46–49



carrot 3.

Always on the move.

GMFCS
Level
III-V

kg
max. 75



cm
25-44

cm
26-38



I am your supportive car seat.

I like car rides. Do you? You can take me on the drive with you. With me you can get anywhere safely and comfortably. I have been crash tested to the highest levels for use in both the front and the back seats of the vehicle.

To make it easier to get into and out of I can turn on a turntable to either 90 or 100 degrees.



With the seat extensions,
the carrot 3. can grow
with the child.





carrot 3.

Seat cover
removable.

Chest pads
fixed with Velcro for stable support.

Back
reclines.

Seat extension
fitted insert for correct leg length support.

Padded headrest
adjustable in height.

Straps
for support and correct positioning.

Anatomical side supports
to help align the body.

Footrest
optional and height adjustable.

Leg rest
optional.

Swivel base
with ISOFIX attachment.

Colour of car
seats

Red

Black



	Size S	Size L	Size XL
Seat depth	25 - 28.5 cm	25 - 28.5 cm	44 cm
Seat width	26 cm	26 cm	38 cm
Back height	30 - 40 cm	45 - 55 cm	50 - 60 cm
Headrest height	52 - 62 cm	67 - 77 cm	73.5 - 83.5 cm
Total length	42 cm	46 cm	46 cm
Footrest height		25 - 46 cm	
Max. load*	10 - 50 kg	10 - 50 kg	40 - 75 kg
Body height	up to 145 cm	up to 160 cm	up to 176 cm
Weight**	from 8 kg	from 8 kg	from 8 kg

* The **carrot 3**. has been tested to ECE R44 / 04 from 15 to 36 kg. In addition, dummies of 10, 49 and 75 kg were tested with success in accordance with the standard.

** In cases with support from GB or consultant.



You can find more information on accessories in the Price List.



Scan this QR code and watch videos about the product directly on our homepage.



GMFCS
Level III-V

kg
max. 16,5



61-105 cm

43 cm

61 cm

360°

34,5°-74,5°



ECE R129/03
CERTIFIED

i-size
universal
ISOFIX

new

I am your safe companion on the go.

And I am your best friend for relaxed journeys.

With my ergonomic seat, you can sit comfortably even on long journeys. In addition, my unique EPP shell structure provides improved impact protection. In combination with the side impact buddy, which absorbs up to 40 % of the energy from an impact, you are optimally protected. So that you can experience car journeys together for as long as possible, you can install me facing forwards or backwards. Thanks to the Isofix function, installing me in the back seat is child's play! Your parents are guaranteed to love my 360-degree rotation! They can easily turn me when they lift you in or out.





The EPP shell structure
of the seat offers improved
impact protection.

Adjustable tilt angle
for 4 different sitting and
reclining positions.

Isofix mount
for use in
in both directions.

360° Rotation
for comfortable
installation and removal.

Head rest
can be adjusted in height.

With the 5-point-harness
the child stays safely in place.

Lateral pads
are adjustable in width.

Abduction wedge
enables optimum positioning.

Foot rest
to relieve the legs,
the back and neck.

Height adjustment of the footrest
the car seat grows with the child.


**Rear facing
travel position.**

0 to
16.5 kg

**Forward facing
travel position.**

9 to
16.5 kg

**Comfortable and
easy to get into
and out of.**



	marie integral.
Recommended body size	61-105 cm
Length	76 cm
Width	43 cm
Height	61 cm
Seat height	12.5 cm (lowest point to the horizontal plane in the seat)
Seat depth	28 cm
Seat width	28 cm
Arm rest height	34 cm
Back rest height	Max. 41 cm
Tilt-in-space	Backwards: 34.5°-49.5°; Forward: 59.5°- 74.5° measured to horizontal plane
Head rest height	19 cm
Footrest height	- 2 cm to +16 cm to the horizontal surface
Weight	15.6 kg
Max. load	16.5 kg
Fastening the safety device	5-point-harness



You can find more information on accessories
in the price list.



charlie integral.

Relax and enjoy the ride!

GMFCS
Level III-Vkg
max. 45

100-150 cm

43-51 cm



I'm your cool car seat. Here to stay with you all the way into your teenage years. I'm your trusty companion on exciting trips, helping you gain more independence along the way.

Operating me is simple and intuitive, so you can quickly learn to fasten your seatbelt on your own. This boosts your confidence and saves your parents from struggling to fasten it in a way that isn't quite right!

My unique EPP shell structure offers enhanced impact protection, while the Side Impact Buddy absorbs up to 40% of impact energy, ensuring you're always optimally protected during your travels.

new



Version with lap-belt and pelvic support pad.



Version with 5-point-harness.



Version with 5-point-harness with armrest & abduction wedge.



charlie integral.





charlie integral.	
Recommended body size	100 - 150 cm
Length	39 cm
Width	43 - 51 cm
Height	59 - 79 cm
Seat height	10.8 - 16 cm (lowest point to the horizontal plane in the seat)
Seat depth	26 - 32 cm (short / long seat)
Seat width	30 - 37 cm (back / front)
Back height	29 - 47 cm
Headrest height	23 cm
Lower leg length when using the footrest	21 - 42 cm
Weight	7.7 kg
Max load.	45 kg

Optional foot rest
is ergonomically shaped and
height-adjustable.



You can find more information on accessories
in the price list.



Walking + Cycling.



A bike trip that inspires.



Marie loves the wind blowing around her nose.
With her **momo tricycle**, she is often on the road.
Our other walking and cycling aids are also
really good assistants, and allow you to
discover the world independently.



Walking + Cycling.



Marie is simply unstoppable with her tricycle! She loves to go on excursions on her own and to explore the world. Meanwhile she can always rely on her momo tricycle. And in addition: cycling trains muscles and balance. It also strengthens psyche and well-being and promotes walking.

Walking is a natural process although it has to be learned. Although their steps are very wobbly at first, healthy children are able to walk properly from the age of approximately three years. Did you know that when walking your whole body is involved – from the skeleton and joints via the muscles and nerves, right up to the lung and heart.

Anterior or posterior?

The process of walking is very complex. Shifting the weight from one leg to the other, rotating the hips, tilting the pelvis, rocking the body from one side to the other – these are processes with which children with disabilities require assistance. Walking aids should already be used during the first attempts at walking.

Depending on the treatment and therapy objective, we distinguish between anterior and posterior walking aids. Anterior walking aids are aids which are pushed in front of the body and which place higher demands on the child. Using these, for example, the torso musculature can be trained. If the torso musculature is particularly weak, posterior walking aids are recommended which are pulled behind the body. Here the weight of the body is supported proportionally by the walking aid.

Becoming independently mobile – not just the right walking aid can assist here. Tricycles or therapy bikes can also support children with disabilities in becoming mobile. An important aspect of this is that muscles and balance are trained, incorrect positioning is improved and musculature tension is controlled. In addition, the feeling of well-being is increased because the cardiovascular system and also the digestive system are stimulated. It therefore does both body and soul good to explore the environment with your own bicycle. And a therapy bicycle looks almost like a conventional bicycle!

Walking + Cycling.



malte.
54–57



malte outdoor.
58–59



marcy.
60–63



momo tricycle.
64–67

GMFCS
Level III-IV

kg
max. 100



46–124 cm

29–113 cm



I am your versatile posterior walking aid.

I am always there for you – even from your first year of life. Walking is my biggest hobby. It does not matter where you want to go; I will follow you. If you support yourself on my handles, you will quickly notice how I take on your body weight. In addition, my special abducted frame shape supports you so that you do not swing outwards when walking. My pelvic support pad ensures that you walk upright and straight. Cool, aren't I?!





Combined pelvic and buttock pad
can become a backrest in one
single easy step.

Aluminium frame
can easily be folded up.

Backwards stop
prevents inadvertent
rolling back.

EVA guide wheels
make the walking aid manoeuvrable.

Forearm rests
with handles.

Padded seat
can easily be folded back.

Direction locks
facilitate
walking straight.



Abducted frame shape
supports you when walking
and prevents you from
swinging outwards.



Aluminium frame
can easily be folded up via
folding joints.

Frame colour
size 0 mini

● Lazurite blue
metallic

Frame colours
size 1

● Telemagenta

● Sky blue

Frame colours
size 2

● Lazurite blue
metallic

● Traffic red

Frame colours
size 3 and size 4

● Anthracite
metallic

○ Pure white



	Size 0 mini*	Size 1	Size 2	Size 3	Size 4
Handle height	29 - 53 cm	34 - 58 cm	51 - 75 cm	68 - 94 cm	87 - 113 cm
Height of forearm rests	-	46 - 70 cm	63 - 87 cm	79 - 103 cm	98 - 124 cm
Handle spacing	22 - 36 cm		30 - 44 cm	38 - 52 cm	43 - 57 cm
Seat height	-	33 cm	45 cm	60 cm	70 cm
Total width	52 cm	51 cm	61 cm	68 cm	75 cm
Total length	58 cm		70 cm	77 cm	87 cm
Frame height (folded)	37 cm	25 cm	32 cm	34 cm	37 cm
Transport dimensions (L x W x H)	58 x 52 x 37 cm	58 x 51 x 25 cm	70 x 61 x 32 cm	77 x 68 x 34 cm	87 x 75 x 37 cm
Max. load	20 kg	30 kg	50 kg	75 kg	100 kg
Weight	3.9 kg	7.1 kg	8.3 kg	8.5 kg	10.2 kg

* Frame cannot be folded up



You can find more information on accessories
in the price list.





I am your outdoor posterior walking aid.

Just like my brother, I am always with you when you really want to get moving, and will offer you the best-possible support to do so. Thanks to my outdoor outfit, I am happiest outdoors, and am the perfect companion for a walk: my large wheels will roll along any surface and my drum and parking brakes provide additional security.



Optional anti-tip mechanism
protects you from tipping back.

Frame colours size 1



Frame colours size 2



Frame colours size 3





	Size 1	Size 2	Size 3
Handle height	34 - 58 cm	51 - 75 cm	68 - 94 cm
Height of forearm rests	46 - 64 cm	63 - 87 cm	79 - 103 cm
Handle spacing	22 - 36 cm	30 - 44 cm	38 - 52 cm
Seat height	33 cm	45 cm	60 cm
Total width	56 cm	66 cm	73 cm
Total length	66 cm	75 cm	82 cm
Frame height (folded)	25 cm	32 cm	34 cm
Transport dimensions (L x W x H)	66 x 56 x 25 cm	75 x 66 x 32 cm	82 x 73 x 34 cm
Max. load	30 kg	50 kg	70 kg
Weight	7.1 kg	8.3 kg	9.8 kg

Did you know that I am one of the few walking trainers which has an outdoor approval? With larger and wider wheels and a drum and parking brake, I am well-equipped for big adventures outdoors. Brilliant!



You can find more information on accessories in the price list.





GMFCS
Level II-IV



I am your anterior walking aid, and will assist you in your very first attempts to walk.

I can hardly wait to start with you! I will support your balance and simultaneously train your torso musculature. You can adapt my wheelbase to suit: if you adjust it so that it is somewhat shorter, I become highly manoeuvrable and you can move me better. I give you more security and support with a long wheelbase. If you find it hard to walk in a straight line, I am happy to assist you with my direction locks.





Handles

can be adapted in height and width.

Aluminium frame

highly compact; sizes 2+3 can be folded.

Friction brakes on the rear wheels

throttle the speed so that you can concentrate completely on walking.

Storage basket

provides space for toys or drinks.

Direction locks

support you in walking straight ahead.

PU tyres

are highly resistant.

**Three dimensionally-adjustable forearm supports**

stabilise your gait pattern.

**Optional parking handbrake**

provide you with more safety and control - also in combination with the forearm supports.

**Frame colours
size 0+1****Frame colours
size 2****Frame colours
size 3**



	Size 0	Size 1	Size 2	Size 3
Handle height	30 - 47 cm	38 - 59 cm	53 - 75 cm	63 - 88 cm
Height of forearm rests	41 - 58 cm	49 - 70 cm	63 - 85 cm	73 - 98 cm
Handle spacing	18 - 32 cm	22 - 36 cm	30 - 44 cm	36 - 50 cm
Seat height	25 cm	30 cm	45 cm	
Total width	44 cm	48 cm	62 cm	71 cm
Total length	36,5 cm	42,5 cm	57,5 cm	61 cm
Total length (telescopic wheelbase)	39 - 47 cm	45 - 55 cm	60 - 74 cm	65 - 79 cm
Max. load	20 kg	30 kg	50 kg	65 kg
Weight	4.2 kg	5 kg	7.4 kg	9.1 kg



You can find more information on accessories in the price list.



momo tricycle.

Let's go for a little ride.

GMFCS
Level I-V

kg
max. 120



30-98 cm

16-23 cm

I am your cool tricycle. I can be there for you whenever you want to take a ride around your streets or a lake. It's fun, but we also train and mobilise your joints, muscles and your circulation. You can also select me in different versions, saddle sizes, handlebar types and colours. Your age does not matter, because I am available in lots of different sizes: for children, adolescents and adults.





momo tricycle.



**You can find a diverse
range of accessories**
for optimum support & individual
equipment in the price list.



**You have the choice! Between lots of
handlebars, saddles, foot pans,
head supports and different drives!**



Frame colours
12", 16" & 16"/12"

-  Window grey / Pure white
-  Telemagenta (Pink) / Pure white
-  Turquoise blue / Pure white

Frame colours
20" & 24/20"

-  Window grey / Pure white
-  Pastel blue / Pure white
-  Kiwisorbent green / Pure white

Frame colours
24", 26" & 26" (XL)

-  Window grey / Pure white
-  Ivory / Pure white
-  Black matt / Pure white



	12"	16/12"	16"	20"	24/20"	24"	26"	26" (XL)
Pedal to top edge of saddle*	30 – 43cm	35 – 48 cm	42 – 62 cm	49 – 71cm	53 – 75 cm	55 – 77 cm	65 – 87cm	72 – 94 cm
Pedal to top edge of saddle**	33 – 40 cm	38 – 45 cm	45 – 60 cm	52 – 67cm	56 – 71 cm	58 – 81cm	68 – 91cm	75 – 98 cm
Saddle support tube to front tube	33cm	40 cm	39 cm	45cm	54 cm	52 cm	58 cm	
Wheelbase	68 cm	75 cm	80 cm	91cm	99 cm	104 cm	112 cm	
Access point height	16 cm		18 cm	20 cm		22 cm	23 cm	
Full length	102cm	110 cm	120 cm	140 cm	153 cm	165 cm	178 cm	
Full width	58cm	65 cm		69 cm	75 cm			
max. load	40 kg	50 kg	60 kg	80 kg	90 kg	100 kg	120 kg	
Tare weight	12 kg	13.5 kg	14.5 kg	17 kg	18.5 kg	19.5 kg	21.5 kg	21.7 kg

* Dimensions with standard saddle post / ** Dimensions with T-saddle post



**You can find further information
on accessories & harnesses
in the price list.**



Standing + Transfer.



Listening to jokes is great fun.



Siphon is so happy with his **till**.
He is really attentive and awake as
physiotherapist Jan makes the fine
adjustments to the flexible knee pads
FKP. and jokes around with him.

Standing + Transfer.



Finally, Siphon can now stand at school and is no longer confined to his nursing bed. This new freedom is thanks to his supine standing system, till.. At last, he can be at eye level with his friends and experience school in a whole new way. What Siphon doesn't know yet is that till. offers much more than just a new perspective - it has so much more to give!

Children and adolescents who are dependent on seating aids and wheelchairs should take a break from long phases of sitting with regular standing exercises. Sitting for too long is not healthy and can cause muscle shortening. Standing also has health benefits for your child; it increases bone density, promotes digestion, reduces joint contractures and facilitates breathing. The shifting of weight onto the legs helps to develop strong and sturdy joints around the hips, knees and ankles.

However, standing isn't just important for health, children become much more attentive when up-right and can perceive their environment better. When standing children can be more involved in their environment and can interact with their friends better. If standing frames provide a very small standing area, it is easier for little or big children to get close to their friends or to reach the school blackboard.

Standing correctly is easy!

During periods of standing, it is important to observe the position of the pelvis. Using the correct adjustments and body support of our standing aids, the pelvis is always kept symmetrical so that each leg takes equal weight. The head and trunk should then be centralised. Depending on the child's condition or treatment objective, the child can be tilted forward or backwards. Standing with the support against the stomach stimulates the child to stretch more actively. Standing with the support against the back is good for a child who is unable to support their own head for long periods.

Standing + Transfer.



till.
72-77



todd.
78-81



tim.
82-85

Positions you well.


GMFCS
Level IV-V

kg
max. 100



75–185 cm

0°–90°



I am your height adjustable supine standing frame. With me, your transfer from the lying to the standing position is safe and comfortable – even if you are not actively involved. With my various adjustment options, I adapt very precisely to your specific body shape – from head to toe. Regular standing has a positive effect on your body structures and functions and promotes your activities. In doing so, I support you and at the same time offer you security. With me you have a new perspective – you can be at eye-level with your friends and can be moved around. This makes participation in different situations and environments during standing possible and enjoyable.


Height adjustment

is easily done via motor or with gas pressure spring support.




Growth adjustments

can easily be made by adjusting the length of the back and legs.







Push bar
for easy transfer.

Headrest
three-dimensionally
adjustable.

Upper body support
offers comfortable support.

Trunk and pelvic support
supports you and facilitates
a stable and upright
standing position.

Therapy tray
with padding offers plenty
of space for activities

Knee angle
adjustable up to 65°.

Flexible knee positioners (FKP)
are individually three-dimensionally
adjustable. Length-adjustable
contact straps adapt optimally to
the shape of the knee.

Footrest
individually height adjustable with adjustment
for dorsi/plantar flexion and eversion/inversion,
with heel edges incl. strap with magnetic closure.

5 Castors
secure, stable and easy to move



You can also easily adjust my height and angle. So nothing stands in the way of playing with my friends.



Precise tilting of the lying surface
is possible from 0-90° via
gas pressure spring.





Headrest with adjustable side guide

is flexible, adaptable, practical and offers great therapeutic benefits.



Flexible knee positioners (FKP)

are individually three-dimensionally adjustable. Length-adjustable contact straps adapt optimally to the shape of the knee.



An optimal foot position

is achieved by individually adjustable 3D footrests. These can be lowered to the floor.



Perfect leg guidance

is made possible by legs adjustable in length and angle.



Precise abduction of the legs
is possible from 0-30°
via adjusting screw.



		Size 1	Size 2	Size 3	Size 4
Body size		75 - 115 cm	100 - 135 cm	120 - 160 cm	145 - 185 cm
Back	Hip/thorax width	14 - 27 cm	17 - 31 cm	20 - 37 cm	27 - 46 cm
	Length of lying surface	30,5 - 43 cm	38 - 48 cm	49 - 61 cm	58 - 69 cm
Legs	Sole to knee	20,5 - 30,5 cm	27,5 - 37 cm	32,5 - 45 cm	42 - 54 cm
	Knee to crotch	12 - 18,5 cm	16,5 - 24 cm	20 - 29,5 cm	24,5 - 34 cm
	Knee angle	0° - 65°			
	Hip angle	0° - 45°			
	Abduction	0° - 30°			
	Level compensation	0 - 5 cm			
Feet	Equinus	0° - 20°			
	Flat foot	0° - 15°			
	Rotation	0° - 15° per direction			
Lying surface height (for transfer)		53,5 - 73,5 cm	60 - 85 cm	71,5 - 93,5 cm	82,5 - 104,5 cm
Tilt		0° - 90°			
Transport dimensions (L x W x H)		95 x 51 x 75 cm	105 x 61 x 85 cm	125 x 76 x 97 cm	152 x 77 x 110 cm
Max. load		30 kg	60 kg	80 kg	100 kg
Weight		30 kg	35 kg	80 kg	90 kg



You can find more information on accessories
in the price list.



todd.

Looks good with you.

GMFCS
Level III-IVkg
max. 80

75–175 cm

50–99 cm

30–84 cm

(-)0° – 15°




I am your moveable standing trainer.


Using me, you can practice standing on your own. Regular standing trains your muscles and improves bone density. My body supports, which can be adapted in height, width and depth, can adapt to your body measurements. My two large wheels round me off perfectly. Using these, you can go anywhere – completely independently.



Flexible knee positioners (FKP)

are individually three-dimensionally adjustable. Length-adjustable contact straps adapt optimally to the shape of the knee.





Pelvic frame
can simply be adjusted in depth using a crank.

Frame colours



Pure white



Melon yellow



Anthracite metallic



Playing cards together at eye level is much more fun!

		Size 1		Size 2		Size 3		Size 4			
Body size		75 - 110 cm		100 - 130 cm			120 - 155 cm		150 - 175 cm		
Wheel size		24"	28"	30"	32"	36"	32"	36"	36"	40"	42"
Grip height: footplate - hand rims at top		40 - 48 cm	50 - 58 cm	57 - 65 cm	64 - 69 cm	72 - 80 cm	60 - 68 cm	71 - 79 cm	82 - 90 cm	87 - 95 cm	
Recommended from elbow height		50 cm	60 cm	70 cm	80 cm	92 cm	82 cm	92 cm		97 cm	99 cm
Width of hand rims at top	0°	-	-	-	-	-	65 cm		67 cm		
	3°	-	-	54 cm	53 cm	52 cm	61 cm	60 cm	63 cm	62 cm	61 cm
	6°	-	-	51 cm	49 cm	47 cm	57 cm	55 cm	59 cm	57 cm	56 cm
	9°	45 cm	44 cm	48 cm	45 cm	42 cm	54 cm	51 cm	55 cm	53 cm	51 cm
	12°	42 cm	41 cm	45 cm	41 cm	38 cm	-	-	-	-	-
	15°	40 cm	-	-	-	-	-	-	-	-	-
Width of hand rims at bottom (min. - max.)		68 - 74 cm	70 - 73 cm	71 - 83 cm		71 - 84 cm	73 - 83 cm	73 - 85 cm	76 - 86 cm		76 - 90 cm
Footplate height - centre of knee support pads		19 - 30 cm		28 - 36 cm			34 - 44 cm		42 - 51 cm		
Footplate height - buttocks support pad support		26 - 42 cm		43 - 58 cm			53 - 69 cm		65 - 80 cm		
Recommended for inside leg length		30 - 47 cm		48 - 63 cm			58 - 74 cm		70 - 84 cm		
Abduction split footplate		0° - 6°					0° - 30°				
Abduction split footplate (with abducted basic model)		0° - 40°					-		-		
Pelvis depth*		11 - 24 cm					15 - 28,5 cm				
Pelvic depth** (adjustment range of the pelvic position by hand crank max. 6 cm)		11 - 20 cm (+6 cm)					15 - 27 cm (+6 cm)				
Sternum height (Footplate height - centre of chest support pad)		54 - 80 cm		68 - 101 cm			92 - 123 cm		111 - 140 cm		
Tray height from footplate		43 - 66 cm		61 - 79 cm			72 - 96 cm		92 - 110 cm		
Max. load		30 kg		45 kg			65 kg		80 kg		
Weight		28 kg		31 kg			38 kg		42 kg		

* with bracket for rear support pad (standard) ** with bracket for rear support pad (flip away)



You can find more information on accessories
in the price list.



GMFCS
Level III-V

kg
max. 80



75-160 cm

(-)15° - 30°



I am your personal standing frame.

Using me, you can practice standing on your own two feet. Together, we will get your muscles fit, strengthen your bones and get your circulation going. Don't be afraid, I will be holding you all the time. Even if you lean forwards or backwards - my tilt angle can be adjusted. All my components can be adjusted in height, depth or angle precisely to suit your requirements.







You can order my wooden tray with or without adjustable arm supports, a bowl or a tray cover for therapy or play purposes.

Frame colour
Size 1

 Melon yellow

Frame colour
Size 2

 Dormant red

Frame colour
Size 3

 Lazurite blue
metallic

The pelvic frame helps you into an optimum upright position. Using a crank, your parents or carers can adjust the frame so that you can slowly stretch out your legs.



Flexible knee positioners (FKP) are individually three-dimensionally adjustable. Length-adjustable contact straps adapt optimally to the shape of the knee.


	Size 1	Size 2	Size 3
Body size	75 – 110 cm	100 – 140 cm	120 – 160 cm
Tray height from footplate	48,5 – 80,5 cm	71 – 95 cm	80 – 112 cm
Pelvis depth	12 – 22 cm	16 – 28 cm	20 – 35 cm
Depth adjustment of the pelvis frame (via hand crank)	6 cm		
Total length	84 cm	90 cm	100 cm
Total width	57 cm	65 cm	71 cm
Abduction split footrests	0°–27°		
Pelvis height	45 – 60 cm	57 – 80 cm	65 – 90 cm
Sternum height	55 – 85 cm	75 – 110 cm	95 – 125 cm
Max. load	40 kg	60 kg	80 kg
Weight	27 kg	29 kg	32 kg



You can find more information on accessories
in the price list.



Bathroom + Toilet.




Theo loves the water. He enjoys feeling it splash onto his skin. And the **boris.** helps him to really relax. Our other bathing aids are also ready and waiting and will accompany you in everything you want to do in the bathroom.

Bring on the water!



Bathroom + Toilet.



Wow, that was fun! All afternoon, Theo has been playing with his sisters on the playground, building sandcastles with them - but now it's time for a shower. Because Theo can't stand alone, Boris assists him in the bathroom with hygiene matters where it can. So, bring on the water!

Body care such as showers, baths, cleaning our teeth and trips to the toilet are some of the things we do every day in the bathroom. Children and teenagers with physical disabilities often cannot deal with their bodily hygiene on their own. They require assistance. To help children to accept their daily care in the bathroom and also to achieve independence, for example when going to the toilet, there are hygiene aids available for different requirements. These aids all have one thing in common; they offer children sufficient support and security, and make everyday life easier for parents or companions.

Sometimes, you simply don't have enough hands to wash a child in the bathtub or to support them safely in the shower. Here aids are really important to give the child a feeling of security. Many adjustments can be carried out on the aids using just one hand. Accessories such as straps or headrests also help to provide a secure hold. Using the right aids, you will always have one hand free to hold the child or adolescent and provide them with security.

Take me with you, please!

Our hygiene aids have been cleverly designed to make things as safe and pleasant for little and big people as possible. For example, the covers are particularly gentle on the skin and do not leave pressure marks on the body. Being showered on the padded seat with its soft materials feels pleasant, and not cold at all.

The best thing about it is that our hygiene aids can be taken with you - on a weekend trip to your grandparents house for example. Some aids are extra-small, compact and foldable.

Bathroom + Toilet.



boris.
90-93



basti.
94-97



bibi.
98-99

GMFCS
Level III-V

kg
max. 150



26–46 cm

31–48 cm

(-)5°–30°



I am your comfortable shower and toilet chair.

Whether cleaning your teeth or showering, my seat tilt can be adapted to any situation. My anatomically shaped seat surface gives you the required support to stabilise your upper body and can be adjusted in height. You can sit comfortably on my soft PU pads which can be easily removed and cleaned. Because I have been produced from highly robust material, you and your parents will have a lot of fun with me for a long time to come. Did you know that I am available in five sizes and can carry up to 150 kg?



You can sit really comfortably on the **removable seat pad**.



Swing away arm supports also help to stabilise you.



Headrest

can be adjusted in height and angle.

Tray

with padded armrest.

Seat pan

Anatomically moulded, provides stability in the pelvic area and stabilises the upper body.

Insertion pad

prevents pressure marks on the body.

Robust sub-frame

made from stainless steel.

Footplate

can simply be folded up for easier access.


Foot pedal

adjusts the seat pan continuously from (-)5° to 30°.

**Toilet equipment**

You can use me as a toilet chair independently of a WC.

I am also available with a larger hygiene opening in sizes 2-4.





		Size 0	Size 1	Size 2	Size 3	Size 4
Seat width (rear)		26cm	31cm	36cm	41cm	46cm
Seat depth		31cm	36cm	41cm	46cm	48cm
Seat shell width (front)		31cm	36cm	41cm	46cm	51cm
Back height		49cm	51cm	56cm	61cm	65cm
Seat angle		(-) 5° - 30°				
Can be positioned over a toilet		47 - 57 cm	45 - 55 cm	42 - 52 cm		
Total width		56cm	61cm	63 cm	71,5cm	74 cm
Total length		65cm			70cm	
Max. lower leg length		38 cm		48 cm		
Front edge of hygiene opening to rear edge of sub-frame		47 cm	49cm	51cm	54 cm	57cm
Front opening width	A	10cm				
Rear opening width	B	12 cm			13,5cm	
Rear opening width*	C	-	-	18,5cm		
Max. load		60 kg	80 kg	100kg		150 kg
Weight		28kg	30 kg	33kg	35kg	38kg

* Only applies to sizes 2 b, 3 b and 4 b



You can find more information on accessories & harnesses in the price list.



GMFCS
Level III-V



I am your versatile bathing and shower cradle.

I will adapt completely to your requirements, whether in the bathroom, in the sea or on the beach. You will always lie safely and comfortably on me! My cover is washable to 30 °C, is particularly skin-friendly and leaves no pressure marks on your body. Because you can sit very low in the bathtub with me, I can save a lot of water. Aren't I practical? What's also practical is that my back angle can be adjusted with one hand during use - leaving the other hand free to support you!



Thanks to the adjustable headrest, no water will run into your eyes.



Head section

can be easily folded.

Head guide pad

can be adjusted in height and width.

Cover

can be removed for cleaning purposes.

Backrest

can be adjusted from 90° to 180°.

Chest strap

provides reliable retention.

Knee angle adjustment

provides a comfortable sitting position.

Optional pelvic mould depression

stabilises and positioned the pelvis, with adjustable seat support for ischiatic bone.

Lower leg support

can be adjusted between 0° and 90°.

Sub-frame

for the shower.

**Pelvic mould depression**

optionally, a cover is available with an integrated pelvic mould depression.

Back angle

can be adjusted even during use with one hand.



Size 4 has two stainless steel joints.

Frame colour



Jet black



	Size 1	Size 2	Size 3	Size 4
Seat width	30 cm	35 cm		40 cm
Seat depth	30 cm	34 cm	38 cm	44 cm
Back height	45 cm	61 cm	71 cm	82 cm
Standing space length	69 cm	74 cm	88 cm	93 cm
Lower leg length	19 cm	20 cm	28 cm	
Back angle	90° - 180°			
Knee angle	0° - 90°			
Headrest angle	0° - 90°			
Headrest angle	0° / 5 cm			
Max. load	35 kg	45 kg	55 kg	75 kg
Weight	4 kg	4,1kg	4,3kg	4,8 kg

* In combination with the cover with pelvic mould depression



You can find more information on accessories & harnesses
in the price List.



GMFCS
Level II-IV



Backrest
can be adjusted
from 90° to 180°.

Pelvic harness
provides reliable retention.

Seat surface
moulds to fit body
proportions.

I am your mobile bathing seat. I would really love to go bathing with you. But not only that I would also be delighted to build castles with you at the beach. Hair full of sand, mucky hands – we have earned ourselves a shower! My backrest can be adjusted so that you can either sit bolt upright or lean back and relax. By the way I am so lightweight and compact in size that you can just fold me up. That's really practical when I accompany you on journeys or to friends' houses.




Back angle adjustment
brings your back into the
required position.

Frame colour



Jet black



	bibi.	bibi. with long back
Seat depth	30 cm	
Seat width	36 cm	
High back	45 cm	55 cm
Total length	60 cm	
Back angle	90° - 180°	
Max. load	30 kg	
Weight	2.3 kg	


Aluminium frame

It can be easily folded up to save space.




You can find more information on accessories & harnesses in the price list.

General guidelines.



Look, the stain has disappeared!



"Now we can wipe down the tray, and **tim.** will be clean soon"! Not just **tim.**, all of our aids can be easily cleaned and dirt free. We explain on the following pages how they should be cleaned. Even more exciting, choose your favourite pattern from our large selection of fabrics and colours.

Introduction to using our products.

Please note that **people with special needs, reading difficulties and impaired vision** must ask a parent, supervisor or technical advisor to read or explain this technical documentation so that safe usage of our products can be guaranteed.

We would like to point out that only **spare parts and accessories** from our company may be used and installed so as not to endanger the safety of the users. In individual cases, the combination of our products with those of other manufacturers may require clarification with us as the manufacturer.

Settings and adjustments to the product or accessories may only be made by people who have been given the necessary instructions by a product advisor. Please ensure that none of the user's extremities are in the respective area when making adjustments of any kind to minimise the risk of injury

The product should be reduced to its smallest dimensions **for transportation**. For example, reduce the seat, the armrests, the back etc. to the smallest dimensions. Remove the parts of the product that can be easily removed, e.g. the wheels of the **todd.** standing frame. Fold up the foot support on products with this feature. Please see the User Manual belonging to the respective product for more detailed information.

Instructions for use on the products are available on www.schuchmann.co.uk/media-library.

If applicable, the **product-specific cleaning guidelines** in the respective instructions for use are to be taken into consideration.

Cleaning guidelines for the different materials.

The care and cleaning of our products should be carried out at regular intervals to ensure that they remain functional over their useful life. The choice of suitable cleaning agents is crucial.

Wood: Wooden parts should be cleaned with a damp cloth or similar. When using detergents, please consult the manufacturer's instructions regarding detergent concentration. Caustic household detergents or corrosive substances are not suitable. Never leave natural wood in a damp or wet condition over long periods, as it will take on water and can change in shape (swelling). Therefore, only wipe and clean damp (i.e. with wrung-out cloths). Test a small, non-visible area for compatibility.

Powder-coated parts: Aluminium parts should be cleaned with a sponge or damp cloth using normal antibacterial household detergents.

Stainless steel: Stainless steel parts should be cleaned using a sponge or damp cloth with a mild soap. Bleach must not be used.

Covers / Fabric: Covers and fabric should be washed at 40° C in the washing machine using normal detergent. Do not tumble dry.

Plastics: Plastics should be cleaned using a sponge or damp cloth with mild soap. Bleach must not be used.

PU (polyurethane): PU parts should only be cleaned using a sponge or damp cloth with a mild soap. Bleach must not be used.

Should you have any
queries, you are welcome
to contact us directly!



Disinfection measures by material type.

Various products can be used for **surface disinfection** of **metal, wood** and **plastic parts**. The disinfectants used should have a wide range of efficacy (at a minimum, bactericidal and levurocidal; if necessary, also virucidal) and be approved for medical devices.

Suitable disinfectants can be found in the List of Disinfectants Certified by the VAH (Association for Applied Hygiene). It must be ensured that only disinfectants classified as medical devices themselves are used. Other disinfectants are not approved for the disinfection of medical aids.

Liquid disinfectants are available as a ready-to-use solution which is sprayed on and evenly applied with a soft cloth. Alternatively cloths soaked with disinfectant can be used and wiped over the whole product. In both cases complete coverage is necessary.

The exposure times may vary and can be found in the manufacturer's instructions for the agents used. Products based on alcohol should not be used.

For occupational safety purposes, disinfection systems may be used. There the contaminated devices are completely disinfected by cold fogging. Methods of disinfection by icing are not permitted.

If necessary, **product-specific disinfection instructions** from the respective instructions for use must be taken into account.



Hygiene measures for all our products.

1. Who should carry out the hygiene measures?

Users, careers, service staff from specialist dealers and other authorised institutions.

2. How regularly?

- Cleaning: As required (maintenance/repair/care)
- Cleaning + disinfection: Regularly and as required (in case of possible contamination and before maintenance and servicing)
- Prior to cleaning and disinfection: Damaged sub-frames or components must first be removed or repaired.
- Before reuse: reprocessing by trained personnel (thorough cleaning and disinfection). It is recommended to replace the cushions and covers!

3. With what should the work be carried out?

We recommend wearing protective gloves (e.g. made of nitrile or butyl rubber)! Avoid the penetration of liquids into Bowden cables, frame tubes or similar. Remove visible soiling with a damp sponge cloth and, if necessary, a cleaning agent. Do not use scouring agents or aids!

4. How should the work be conducted exactly?

For cleaning, all accessible surfaces are cleaned with a damp sponge cloth and, if necessary, a suitable cleaning medium. Then dry the surfaces with a cotton or disposable cloth. Hard-to-access places can also be dried using a hairdryer on its lowest setting. The use of high-pressure cleaners is not permitted, as their use can lead to damage to the products. For disinfection the disinfectant should be used according to the manufacturer's instructions. Ensure **complete coverage**. If necessary, remove excess preparation quantities with unstained disposable tissue. It is particularly important to observe the exposure times of the manufacturers!



Now things are getting colourful.





As round as a spinning top and as colourful as a marble run - let's discover the world through play!

The same applies to our fabrics. With all the usual quality demands that can only be made on the supply of aids.

Good to know - all our fabrics are:

- simple and easy to process
- Highly abrasion-resistant
- flame-retardant in accordance with DIN EN 1021-1+2
- skin-compatible in accordance with DIN EN ISO 10993-5
- fulfil the requirements of the MDR for seat shells and special constructions



Should you request,
we are happy to send you a
fabric sample folder!




Fabric Flexi.



9001043 (No. 43)



9001044 (No. 44)



9001050 (No. 50)

Woven fabric.



9001123 (No. 123)



9001124 (No. 124)



9001127 (No. 127)

Knitted fabric.



9001091 (No. 91)



9001092 (No. 92)



9001105 (No. 105)



9001106 (No. 106)




9001111 (No. 111)*



9001112 (No. 112)*



9001114 (No. 114)*



9001115 (No. 115)*



9001113 (No. 113)**




Even cooler: have your name embroidered on the back pad of your seating aid. All materials - except for flexi - can be embroidered with a maximum of twelve letters in silver-white. Your friends will be amazed!

Leather imitation.



9001 190 (No. 191)




9001 192 (No. 192)**



9001 193 (No. 193)**

Mesh fabric.



9001 046 (No. 46)*



Note:

We reserve the right to carry out design, colour and material modifications according to the latest technology and within the scope of standard market production development. The images of the products and colours merely serve as orientation aids. Due to printing processes, deviations from the original cannot be excluded. All data without guarantee.

* This product does not conform to DIN EN 1021-1+2. ** Discontinued model





schuchmann.co.uk

Schuchmann LLP · 3 Lakeside Court · Kingston Park · Peterborough
01733561885 · info@schuchmann.co.uk