

## CLINICAL STIMULATOR LEADWIRE CHART

MANUFACTURER	MODEL #	UNIT MADE PRIOR TO 2000		COMPLIANT UNITS
		ADAPTOR #	WIRE #	WIRE #
Amrex	HV752	AD4500	W4872 (Ban.-Ban.), W4772 (Ban.-Pin)	W4872 (Ban.-Ban.), W4772 (Ban.-Pin)
	LVG325A	AD4500	W4872 (Ban.-Ban.), W4772 (Ban.-Pin)	W4872 (Ban.-Ban.), W4772 (Ban.-Pin)
	MS Series (all)	AD4500	W4872 (Ban.-Ban.), W4772 (Ban.-Pin)	W4872 (Ban.-Ban.), W4772 (Ban.-Pin)
	Z-Stim IF Series	AD4500	W4872 (Ban.-Ban.), W4772 (Ban.-Pin)	W4872 (Ban.-Ban.), W4772 (Ban.-Pin)
Chattanooga	CPS 200 Combo, CPS 200 Stim	No CMP solution for U.S.A.	W18112 (Cdn. Only)	W7172
	CPS 400 Combo, CPS 400 Stim	No CMP solution for U.S.A.	W18112 (Cdn. Only)	W7172
	Forte 200, Forte 400 Combo	No CMP solution for U.S.A.	W18112 (Cdn. Only)	W7172
	Forte CB 450, Forte ES 450	AD8500	W6572	
	Intelect 120, 170, 775MP, 750MP, 600, 900	AD8500	W6572	
	Intelect 500, 500S, 700, 700C, 775		W2372	
	Intelect 550, 551	AD4500	W4472 or W4372	
	Intelect Combo, Legend Stim	No CMP solution for U.S.A.	W18112 (Cdn. Only)	W7172
	Intelect IFC2 (squared)	AD8500	W6572	
	Intelect Legend XT, Vectra Genisys, Transport			W17672
	Intelect Low Volt 110			W7972
	Vectra Pro, Vectra			W7172
	VMS		W2372	
VMS II	AD8500	W6572		
Dynatronics	Dynatron 50, 125, 525	AD6500	W6672	W6672
	Dynatron Plus 550, 650, 850, 950	AD6500	W6672, W6272 (combo)	W6672, W6272 (combo)
Enraf	Endomed 182, 482			W7772 (DIN-pins) W7872 (DIN-ban.)
	Endomed 582, 982			W8372 (DIN-pins) W8272 (DIN-ban.)
	Endomed CV-405			W8596B
	Sonopulse 464			W2872IVBRB
Enraf (Henley)	Endomed 433	AD8400	W8472 (DIN-pins) W8572 (DIN-ban.)	
Excel	Ultra / Multi / Isotron Series	AD6500	W6572, W6372 (combo)	W6572, W6372 (combo)
Fisher	Model 700	AD4500	W4772 (Ban.-Pin)	
H Wave		AD4500	W4472 or W4372	
Hakomed				W8672
Healthcare Mnfg.	Masterstim	AD6500	W6572	W6572
Kaladar (KMI)	The Stim	AD6500	W6572	
LSI	System II and IV, MS300, MS400	AD7500	W7572	W7272
Medco	Medcollator	AD8500	W6572	
Metron	All models	AD4500	W4772 (Ban.-Pin)	W4772 (Ban.-Pin)
			W4872 (Ban.-Ban.)	W4872 (Ban.-Ban.)

**CLINICAL STIMULATOR LEADWIRE CHART cont'd.**

MANUFACTURER	MODEL #	UNIT MADE PRIOR TO 2000		COMPLIANT UNITS
		ADAPTOR #	WIRE #	WIRE #
Mettler	Sonicator Plus 992, 994, 930	AD9500	W7472, W7372 (combo)	W7472, W7372 (combo)
	Sys*Stim 294, 226, 206	AD9500	W7472	W7472
Monad	Myomatic	AD4500	W4472	
Multiplex		AD8400	W8072	
Nemectron	Edit			W8672
Neurocare	Neurocare 1000			W8972IVBRB
Omnisound	2500C	AD4500	W4472BB	
Orthion Corp.	Sentry Mark V Model 9000		W2372	
P.E.T. Eurotech	Magnum 100	AD6500	W6572	
Phyacton				W4772 (Ban.-Pin) W4872 (Ban.-Ban.)
Physio Tech.	Omnistim	AD6500	W6572	
Physiodynamic	Therastim 1000	AD6500	W6412	W7412
Physiomed		AD8400	W8072	
Precision Microcurrent		AD4500	W4472	W4472
Prostim	ET 3000			W18672
Rebox Unit		AD4500	W4472	
Richmar	HV 1000, 2000	AD8500 with or AD9500 with	W6572 W7572	
	HV Combo	N/A AD4500	W2372 (Active) W4472BB Dispersive	
	HV 20		W2672 (Red Active) W2772 (Blk Active) W2372 (Dispersive)	
	IF 200	AD8400	W8072	
	IF 2000	AD8500	W6572	
	III, IIIG	AD8500	W6472	
	RM6	AD8500 with or AD9500 with	W6572 W7572	
	Theramini 2, 3C, 3P	AD8500 with or AD9500 with	W6572 W7572	W4772 (Ban.-Pin) , W4872 (Ban.-Ban.)
	Theratouch 4.7, 7.7	AD8500 with or AD9500 with	W6572 W7572	W4772 (Ban.-Pin) , W4872 (Ban.-Ban.)
VI		W2972RB		
Scrip	Superstim MF Model 4000	AD8500	W6572	W7272IVRR
Sterne	Acutrain V	AD8500	W6572 W5672 (Cdn. only)	W7272IVRR
	Myotrain V	AD8500	W6572 W5672 (Cdn. only)	W7272IVRR
	NeuroTrain III			W7272
US Medical	Model 450		W2372	
Uniphy				W4872
Vectrodyne		AD8400	W8072	
Verimed	Biofeedback			W2172
Williams	Intertron 4000 etc.	AD6500 with or AD7500 with	W6572 W7572	W8172
	ZC2, ZS2, ZS4	AD6500 with or AD7500 with	W6572 W7572	W8172

**PORTABLE STIMULATOR LEADWIRE CHART**

<b>MANUFACTURER/STIMULATOR NAME</b>	<b>LEADWIRE PART #</b>
<b>Biomedical Life Sciences (all models)</b>	W1848
<b>Bioscope TENS</b>	W1736
<b>Cefar</b>	W1442R, W1462R
<b>Cefar Primo, Basic, Basix 3, Femina</b>	W1342
<b>Classic TENS, Classic NMS</b>	W1342, W17742V
<b>Comfort TENS</b>	W1736
<b>Elpha 2000, Elpha 3000</b>	W1042, W1062
<b>Empi 300 PV</b>	W18440PV
<b>Empi Eclipse, Eclipse Plus</b>	W1342
<b>Empi Focus</b>	W18440F
<b>EMS 5000, 6500</b>	W1636
<b>Epix VT</b>	W18440VT
<b>Epix XL</b>	W18440E
<b>ETS, Verity Medical</b>	W1342, W17742V
<b>IF2, IF4000</b>	W1462C, W18860
<b>Intelect TENS, Chattanooga</b>	W1736
<b>Mantra TENS and Mantra Sports</b>	W1342, W17742V
<b>Medstar 100, 150</b>	W1462C, W18860
<b>Myostim EX, ETD</b>	W1848
<b>Neuromod Selectra 7720</b>	W9748
<b>Neurotrac 3, 4, Verity Medical</b>	W1342, W17742V
<b>Polaris Digital, Spectramed</b>	W1342, W17742V
<b>Procura</b>	W1636
<b>Promed 300 TENS</b>	W1636
<b>Protens</b>	W1636
<b>Rehabilitare GV II, NuWave &amp; NM III</b>	W1442R, W1462R
<b>Rehabilitare SMP-Plus &amp; SX Plus</b>	W1442R, W1462R
<b>Respond II</b>	W9748
<b>Rice AR1-32</b>	W1636
<b>RS Medical TENS</b>	W1848
<b>Simplex, Verity Medical</b>	W1342, W17742V
<b>Spartan</b>	W1442R, W1462R
<b>Spectrum Max-SD &amp; II</b>	W1442R, W1462R
<b>Sport IFC, Care IFC</b>	W1462C, W18860
<b>Stadyn EMS+2, Promax</b>	W1442R, W1462R
<b>Stadyn Maxima II &amp; III</b>	W1442R, W1462R
<b>Stadyn Units - Older Style</b>	W18536
<b>TENS 2500, 3000, 6000</b>	W1636
<b>Vision Quest IFC-Plus</b>	W1462C, W18860