




PROCARBON - Pre-gelled Carbon Conductor

For use with TENS & EMS Applications



Catalogue No.	Size	Shape
C5052PF	1.5" x 1.75"	
C5054PF	1.75" x 4.0"	

NOTE: Each package contains 4 electrode pads.




- Reusable, self-adhesive
- Economical
- Pliable, flexible
- Convenient, easy to use

ELECTRODES NON-GELLED

PROCARBON for TENS - Carbon Rubber Conductor

For use with TENS & EMS Applications



Catalogue No.	Size	Shape
C5005PF	1.0" Round	
C5010PF	1.5" x 1.75"	
C5020PF	1.75" x 4"	

NOTE: Each package contains 4 electrode pads.

- Economical, reusable
- Long lasting
- For use with conductive gel and tape or adhesive tac gel
- Available in three popular sizes




CMP is certified to ISO 13485:2003 by a CMDCAS recognized registrar.




PROCARBON for IFC - Carbon Rubber Conductor

For use with IFC and EMS Applications

NON - INSULATED

Catalogue No.	Size	Shape
C5015PF	1.5" x 2.25"	
C5023PF	2.25" x 3.0"	
C5034PF	3.0" x 4.5"	



NOTE: Each package contains 4 electrode pads.

INSULATED

Catalogue No.	Size	Shape
C5120-B/R	2.0"	
C5130-B/R	3.0"	
C5140-B/R	4.0"	
C5124PF	2" x 4.5/2"	Grey



2.0" @ \$16.95
 3.0" @ \$21.50
 4.0" @ \$23.95
 2"x4.5" @ \$26.95
 (Special Order)

NOTE: Each package contains 2 electrode pads; Black or Red.

- High quality, long lasting reusable electrodes
- Very low resistance
- Soft, flexible carbon rubber



VACUUM SUCTION CUPS & SPONGES

VACUUM SUCTION CUP ELECTRODES

Suction cup electrodes with an internal valve and a horizontal hose connector, a major advancement over the previous classic design with the vertical hose connector. The horizontal connection eliminates hose kinking and allows for towel or hot pack placement over the electrodes. The internal valve feature keeps the electrode sealed until it is applied to the patient. This simplifies and speeds up the application of these electrodes. All of our suction cup electrodes are made of high grade grey silicone rubber with stainless steel contact plates and are autoclavable.



SUCTION CUPS

	Small	(30mm / 1.2")	VCUP10S	\$ 24.00
1/Pkg.	Medium	(60mm / 2.4")	VCUP10M	\$ 34.00
	Large	(90mm / 3.5")	VCUP10L	\$ 39.00

COMPRESSED SPONGES

(16mm / 0.6" thick when expanded)

12/Pkg.	Small	(30mm / 1.2")	VSPG10SC	\$10.00
16/Pkg.	Medium	(60mm / 2.4")	VSPG10MC	\$ 30.00
12/Pkg.	Large	(90mm / 3.5")	VSPG10LC	\$ 30.00

EXPANDED SPONGES

(20mm / 0.8" thick)

4/Pkg.	Medium	(63mm / 2.5")	VSPG10ME	\$10.00
4/Pkg.	Large	(85mm / 3.3")	VSPG10LE	\$ 12.00

SPONGES (Other)



SPONGE POCKETS (Rectangular)

	1.7" x 2.5"	F62-0041	----
4/Pkg.	2.4" x 3.1"	F62-0042 (fit C5120)	\$ 13.00
	3.1" x 4.7"	F62-0043	----
	3.8" x 4.0"	F62-4591 (fit C5130)	\$ 18.00
	4.0" x 4.4"	F62-1471	----
	4.25" x 4.5"	F62-4604 (fit C5140)	\$ 16.00

SPONGE PADS (Round)

4/Pkg.	2" Round	F62-0448	----
	3" Round	F62-4491	\$ 4.00
	3.5" Round	F62-0450	----

