











Standz by Jenx

A versatile standing system offering prone or supine standing that accommodates abducted or neutral leg positioning for children between approximately 1 - 9 years.

Key features

- Tool free adjustments
- Individual leg abduction and adjustment
- Offers Prone or Supine standing
- Antimicrobial covers to aid infection control in multi-user environments

Always Supportive...

Address:









Standz by Jenx

A versatile standing system offering prone or supine standing that accommodates abducted or neutral leg positioning.

Inspired by the clinical benefits of abducted standing and created alongside therapists and carers all over the world, Standz ensures the perfect combination of clinical positioning with simplicity of use. Standz offers Prone or Supine standing and easily accommodates either abducted or neutral leg positioning for children between the ages of 1-9 years. Leg length discrepancy is easily accommodated by the individual control for knee and foot positioning, and all adjustments are tool free to enable quick and effective clinical postural adjustments. With an extensive range of accessories, Standz can accommodate children with a wide range of postural requirements.

height – Ideal for accommodating leg length discrepancy. This tool free adjustment can be made in seconds and securely clamps into place using our clever clamping

marquee design.





Standz Frame

Comes complete with:

- 4 lockable swivel casters
- Tilt mechanism with foot pedal release and safety lock
- 2 leg posts marked with numbers for easy positioning of leg supports
- 0 30 degree independent abduction movement of leg posts
- Angle gauge to read frame angle frame angle range of 5 90 degrees

In addition to the base you will need to add the following accessories to ensure a minimum safe configuration.

Prone – sandals, knee support, pelvic support & trunk support.

Supine – sandals, knee support, knee support straps, pelvic support, trunk suppor & head support.

Always Configurable Accessories



Flat Pad Head Support

A simple yet effective head support pad for use in supine for children with reasonable head control. It can also be used in prone to extend the height of the chest pad.

MH09-XX 7"(w) x 2"(d) x 4"(h)



Multigrip Head Support (supine standing only)

The original and innovative Jenx Multigrip™ head support is available in 3 sizes for all Jenx seating/standing systems.

1 Link MH06 11½"(w) x 2"(d) x 3½"(h) 2 Link MH0101 18½"(w) x 2"(d) x 4½"(h) 3 Link MH05 24"(w) x 2"(d) x 4½"(h)



Molded PU Head Support

Simple headrest for children who require minimal head support. Comes mounted on a flexible bracket.

MH07 XX 11"(w) x 2½"(d) x 8½"(h)



Flat Shoulder Support

A simple yet effective wide pad for supporting the shoulders when using Standz in a supine position. Can also be used in conjunction with the **Pelvic Support Board - back pad** in prone.

STZ 014-XX 2"(d) x 12"(w) x 3½"(h)



Lateral Supports with Strap

Firm lateral support pads offer stability to the trunk, pelvis or both together, and the incorporated support strap can be used to facilitate extension or rotation.

STZ 012-XX 5"(d) x 1"(w) x 4"(h)



STZ 013-XX

Rotation Control Strap

A wide support belt designed to be used to position the pelvis, trunk or both together. Threaded through the support fixings, the strap can be fastened securely around the child and then pulled to facilitate extension or rotation.



STZ 015-XX 2"(d) x 7½"(w) x 6"(h)

Pelvic Support Board

For use in prone positioning only, the pelvic support board offers a full contact surface to facilitate the full extension of the pelvis and prevent pelvic rotation, pelvic obliquity and unwanted hip flexion.



STZ 016-XX 6"(w) x 2"(d) x 6"(h)

Pelvic Support Board -Upper Back Pad

Used in conjunction with the pelvic support board, the back pad offers support for the child's upper back. A head support can also be mounted to the back pad to help keep a neutral head position for those children who hyper extend their neck and head when standing in prone.



Multi-Axis Knee Supports

Use for supine standing only. For use where there is the need to accommodate fixed knee flexion.

STZ 017-XX 5½"(w) x 4½"(h)



STZ 007-XX 3-4"(w) x 5"(d) x 4½"(h)

Standard Knee Support

Individual knee supports offering the ability to accommodate knee flexion when used in supine (requires knee straps when used in supine). Width adjustable to correctly gutter the knee in position.



Standard Knee Strap

Intended for use in conjunction with the standard knee supports, the standard knee strap simply holds the knee in the desired position and prevents knee flexion from occurring.

MU54-XX



MU45-XX

Rotation Control Knee Straps

Intended for use in conjunction with the standard knee supports, the rotation control straps allow for a cuffing of the child's leg and allow for control of internal or external knee rotation.



Sandals (Pair)

Available for Standz in either size 1 or 2, these sandals ensure solid foot placement to give the foundation for effective leg positioning and safe standing.

S16 4"(w) x 6½"(d) **S2** 42½"(w) x 8½"(d)



Sandals with Ankle Snugs

Available in 2 sizes, the ankle snug offers a unique magnetic buckle fastening mechanism, which can be operated single handed.

S19 4"(w) x 6½"(d) **S21** 4.3"(w) x 8½"(d)



Multi-position Tray

Cleverly designed to be used in both prone and supine, the tray is angle adjustable to ensure the correct position in

STZ 008-XX 23"(w) x 15"(d)

Accessories continued



Prone Only Tray
Only suitable for when the product is used in the prone configuration to prevent potential clash issues when fully abducting the leg position.





Tray Extension Wings
Strongly recommended for use in supine and just as effective when used in prone, these wings attach to the rear of the tray to offer a support surface for the child's elbows and prevents arms from slipping off the tray.

STZ 011-XX
6"(w) x 6"(d)

6"(w) x 6"(d)



Tray Holder

A clever accessory that offers a solution of where to store the tray when it is not being used.

STZ 010

Standz by Jenx Technical Specification

	Units	Prone	Supine
Age Range (approx)	yrs	1 - 9	1 - 9
Max User Weight	lbs	99	99
Max User Height	inch	541⁄4	541⁄4
Angle Range (support surface angle range)	degrees	5°- 90°	5°- 90°
Chest Width (min - max)	inch	7 - 11¾	7 - 11¾
Hip Width (min - max)	inch	7 - 11¾	7 - 11¾
Kneeblock Width (centre to centre: min)	inch	4¾	4¾
Leg Abduction (min - max)	degrees	0°- 30°	0°- 30°
Footplate Angle (dorsiflexion - planar flexion)	degrees	10°- 25°	25°- 10°
Height of Stander (floor to top of upper pad: in upright)	inch	42½	421/2
Height of Stander (floor to support board: at horizontal)	inch	N/A	31¾
Height from footplate (to top of chestplate: min - max)	inch	22¾ - 39¼	22¾ - 39¼
Base Footprint (min - max)	inch	25 x 351/4	25 x 351/4



Standz Color Range











Address:

975 Fraser Drive, Unit 16 **Burlington, Ontario L7L 4X8** Phone:

Phone: (905) 847-1983 Toll Free: (866) 664-1983 **Email & Website:** Email: info@RandMHealthCare.ca Web: RandMHealthCare.ca

