



THEIR PROTECTION
YOUR ASSURANCE

ProEdgeDental.com | 888.843.3343 | 7042 S Revere Pkwy #400, Centennial, CO 80112

DENTAL WATER SAFETY

Waterline Test Safety Report

Practice Information:

Larry R. Goldstein DDS
2606 Abbey Ct, Alpharetta, GA 30004
770-475-3210

Practice-Reported Treatment Protocol:

Shock:	Shock Date:	Treatment:	Source Water:
Sterilix ULTRA	12/7/2018	BluTab	Filter/Distiller

Test Date: 12/28/2018 Report Date: 1/4/2019

ProEdge Dental Water Labs approves the following results as accurate and complete utilizing our advanced and proven neutralization methods.


Richard Vigil, MS
Laboratory Director


Kellie Thimmes
Water Safety Specialist

Water Lab Results

ProEdge Sample #	Vial #	Location	Room/Chair/Operator #	Device or Source Water	Microbial Count (CFU/mL)	Waterline Result	*Safety Level
18-57898	1	Operator	1	A/W Syringe	0	PASS	✓
18-57899	2	Operator	2	A/W Syringe	10	PASS	✓
18-57900	3	Operator	3	Handpiece	0	PASS	✓
18-57901	4	Operator	3	Source Water	60	PASS	✓

Questions?

Our team of dental water experts are available to help explain your results and take the next steps to achieve or maintain safe water.

Call or Email
888.843.3343
info@proedgedental.com

TNTC = Too Numerous to Count. QNS = Quantity Not Sufficient (of sample water). Pass / Fail results are based on the ≤ 500 CFU/mL standard set by the CDC and ADA.
*Safety Level: Red = Immediate Shock & Retest | Yellow = Immediate Shock | Green = Continue Treatment Protocol. See attached document for more information about your report.
ProEdge Dental Water Labs neutralization methods have been validated by an independent research lab following ASTM E1054-08. For full reports, visit ProEdgeDental.com.

105-108



THEIR PROTECTION
YOUR ASSURANCE

ProEdgeDental.com | 888.843.3343 | 70425 Revere Parkway #400, Camasville, CO 80112

DENTAL WATER SAFETY Waterline Test Submission Form

FILL OUT THIS FORM COMPLETELY FOR EVERY TEST TO ENSURE YOUR SAMPLES CAN BE PROCESSED. KEEP A COPY FOR YOUR RECORDS.

Practice Information:

New Customer Returning Customer

Practice Name: Larry R. Goldstein, DDS, PC Contact Person: Larry Goldstein

Practice Address: 2600 Abbey Court Contact Title/Position: owner

City, State & Zip: Alpharetta, GA 30004-4023 Contact Email: blackbirdsix@gmail.com

Phone: (770) 475-3210 Fax: (999) 344-4037 Compliance Officer Email: _____
Test results will be sent to this email

Waterline Treatment Protocol Information (Select All That Apply):

A. TREATMENT PRODUCT(S) USED:

BluTab ICX
 Citrisil/Citrisil Blue DentaPure Straw/Cartridge
 Sterisil Straw Other: _____
 Install Date: _____
 Shock Prior to Install: Yes No

B. SHOCK PRODUCT(S) USED:

Bleach Citrisil Shock
 Sterilix Ultra
 Other: _____
 Date of Last Shock: 7 DEC 19
 Never Shocked:

C. OTHER TREATMENT SYSTEMS:

Sterisil System O-So Pure
 Vista System
 Other: _____

D. USING WATER BOTTLES:

Yes
 No

E. SOURCE WATER (We fill our water bottles from...):

City Water/Tap Other (e.g. Delivered) _____
 Bottle Distilled In-Office Filter, Distiller or R/O Unit (Brand: Tuflwater)

Test Sampling Information:

Sampling Date: _____

Do not sample on Fridays or before a holiday. Ship same day as sampled.

FOLLOW CLOSELY TO ENSURE ACCURATE TEST RESULTS

- Fill out this form completely every time you test.
- Do not sample or ship on Fridays or the day before a holiday.
- Ship samples on the same day they are taken.
- Do not include more than 16 vials in each shipment.

FOR LAB USE ONLY ProEdge Sample #	Vial Number	Location (Circle One)	Room/Chair/Operator #	Device or Source Water (Circle One)	Notes or Other Information
105	1	Room Chair <u>Operator</u>	1	<u>AW Syringe</u> Handpiece Scaler Source Water	
106	2	Room Chair <u>Operator</u>	2	<u>AW Syringe</u> Handpiece Scaler Source Water	
107	3	Room Chair <u>Operator</u>	3	AW Syringe <u>Handpiece</u> Scaler Source Water	
108	4	Room Chair <u>Operator</u>	3	AW Syringe Handpiece Scaler <u>Source Water</u>	<u>FAUCET</u>
	5	Room Chair Operator		AW Syringe Handpiece Scaler Source Water	
	6	Room Chair Operator		AW Syringe Handpiece Scaler Source Water	
	7	Room Chair Operator		AW Syringe Handpiece Scaler Source Water	
	8	Room Chair Operator		AW Syringe Handpiece Scaler Source Water	
	9	Room Chair Operator		AW Syringe Handpiece Scaler Source Water	
	10	Room Chair Operator		AW Syringe Handpiece Scaler Source Water	
	11	Room Chair Operator		AW Syringe Handpiece Scaler Source Water	
	12	Room Chair Operator		AW Syringe Handpiece Scaler Source Water	
	13	Room Chair Operator		AW Syringe Handpiece Scaler Source Water	
	14	Room Chair Operator		AW Syringe Handpiece Scaler Source Water	
	15	Room Chair Operator		AW Syringe Handpiece Scaler Source Water	
	16	Room Chair Operator		AW Syringe Handpiece Scaler Source Water	

The guidelines outlined by the Center for Disease Control (CDC) for bacteria in water used as a coolant/irrigant for non-surgical dental procedures should be as low as reasonably achievable, and at a minimum ≤ 500 CFU/mL. Your practice's water samples will be tested against this CDC standard or dental water safety.

3/16/17/SE/05

Questions?

Call us today at 888.843.3343

41726