


Axonics Bladder Therapy Information

What is the Device

An implanted wire that is attached to the sacral nerve. The wire sends constant electrical signals through your sacral nerve to your brain to help you feel sensation in your bladder before it becomes urgent.



How Do I Use the Device?

Leave the device on at all times - 24/7!
The device will send constant signals to your brain so you can feel your bladder being full before it becomes urgent.

Do I Adjust the Device?

It is normal to become unaware of the tingling from the device. **Only adjust the strength when you notice an increase in leakage or in urgency.**