Intelligent Textiles[®]

for Medicine

Therapy while you sleep!



Drug Free Pain Relief

FDA Indications

- Pain Relief
- Increase Microcirculation
- Increase Range of Motion
- Retard Disuse Atrophy

Medicare Approved for:

- Acute Post-OP Pain Relief
- Chronic Pain Relief
- Retard Disuse Atrophy (Reimbursed by most private insurance)

Evidence Based Medicine

- 50 years of clinical research with 109 studies on HVPG
- 4 published research on The Silver-Thera Stocking
- Medicare approved procedure codes for wound healing
- U.S. Dept. of Health & Human Services Clinical Practice Guidelines recommend for wound healing



Silver-Thera[™] Stocking



Silver-TheraTM Sock

Nocturnal below Sensory Stimulation: High Volt Pulsed Galvanic




Micro-Z Mini[™] Kit

- Micro-Z Mini[™]
- Velcro[®] Wrap
- Battery and Manual

Studies and Reviews: View these studies at www.prizm-medical.com

- "The Benefit of Electrical Stimulation to Enhance Perfusion in Patients with Diabetes Mellitus"
- "Is Electrical Stimulation Effective in Reducing Pain in Patients with Diabetes?" (1)
- "Electric Stimulation as an Adjunct to Heal Diabetic Foot Ulcers: A Randomized Trial"
- "Electrotherapy in Tissue Repair"

Painful Diabetic Neuropathy Study (1)



Notice the Dramatic Pain Relief



800.447.4422 • 770.622.0933 www.prizm-medical.com contact@prizm-medical.com

CE registered and approved Micro-Z Mini™ Copyrighted Software

Precautions:

- Clinician to determine protocols
- When there is a tendency to hemorrhage
- Over the menstruating uterus
- Specifications subject to change without notice

Contradictions:

- Patients with cardiac demand pacemakers
- Should not be used on cancer patients

Warnings:

- Keep out of reach of children
- Safety has not been established for use during pregnancy
- Persons with suspected or diagnosed epilepsy
- Not to stimulate over the carotid sinus nerves

Caution:

• Federal Law restricts this device for sale by or on the order of a licensed physician or any other practitioner licensed by the law of the state in which s/he practices to use or order the use of this device.

Micro-Z Mini[™] Specifications:

- One AAA battery
- Treatment time 15 minutes up to 24 hours
- Pulse rate 1-120 pps in 5v increments
- Programmable
 Delewiter manifi
- Polarity positive or negative
 Intensity amplitude 0, 100 welt
- Intensity amplitude 0-100 volts
 Wave form: pulsed galvanic monophasic exponential spike
- wave form: pursed galvanic monophasic exponential spike with a fixed interpulse interval of 100 microseconds
 With voltage law exponence pulsed galvanic structure
- High voltage, low-amperage pulsed galvanic stimulation
 Dimensions 2" x 2" x .5"
- Dimensions 2" x 2" x .5"
 Unit weight approximately 2 oz.
- Warranty: 1 year of defective parts and labor

