


# Tested and Approved: Level 2 Face Mask MEDICAL GRADE LEVEL TWO 3-PLY MASK



Brand / Model: TOL GLOBAL HEALTH L2 MEDICAL 3-PLY Adult Size Mask / Non-Sterile / 50 Per Box

Case Pack: 40  
Gross Weight: 19 lbs  
Pallet Count: 11 x 6

EN 14683:2019 + AC : 2019(E) ANNEX C



LEVEL 2 PROTECTION



NO LATEX - NON IRRITATING



BACTERIAL FILTRATION EFFICIENCY  $\geq 98\%$



MEDICAL GRADE

## PRODUCT SPECIFICATIONS

## PRODUCT UPC

K202211
 

 ASTM F2101-19  
 ASTM F1862  
 ASTM F1862M-17  
 ASTM F2100-19e1  
 ASTM F2299

16 CFR PART 1610-2012

**WHITE** 8 46146 05053 6  
**BLUE** 8 46146 05054 3  
**BLACK** 8 46146 05055 0  
**PINK** 8 46146 05056 7

<b>NAME</b>	TOL GLOBAL HEALTH Medical Grade Level 2 Face Mask
<b>DESIGN</b>	Face Mask With Ear Loops
<b>SURFACE</b>	Smooth
<b>COLOR</b>	White, Blue, Black, Pink. Custom Colors & Designs Available
<b>MATERIAL</b>	Inner & Outer layer: Spun-bound Polypropylene Middle Filter Layer: Melt-Blown Polypropylene
<b>STORAGE</b>	Product Shall Be Stored Under Room Condition, Avoid Direct Sunlight
<b>SHELF LIFE</b>	2 Years From Date Of Manufacture With The Above Storage Condition

## TOL GLOBAL HEALTH PACKAGING DIMENSIONS

