


Tested and Approved: Level 3 Face Mask MEDICAL GRADE LEVEL THREE 3-PLY MASK



Brand / Model: TOL GLOBAL HEALTH L3 MEDICAL 3-PLY Adult Size Mask / Non-Sterile / 50 Per Box

Case Pack: 40
Gross Weight: 20.28 lbs
Pallet Count: 11 x 6

EN 14683:2019 + AC : 2019(E) ANNEX C



LEVEL 3 PROTECTION



NO LATEX - NON IRRITATING



BACTERIAL FILTRATION EFFICIENCY ≥ 98%



MEDICAL GRADE

PRODUCT SPECIFICATIONS

PRODUCT UPC



K202211



ISO 9001
ISO 13485



ASTM F2101-19
ASTM F1862
ASTM F1862M-17
ASTM F2100-19e1
ASTM F2299



16 CFR PART 1610-2012

WHITE



BLUE



BLACK



PINK



NAME	TOL GLOBAL HEALTH Medical Grade Level 3 Face Mask
DESIGN	Face Mask With Ear Loops
SURFACE	Smooth
COLOR	White, Blue, Black, Pink. Custom Colors & Designs Available
MATERIAL	Inner & Outer layer: Spun-bound Polypropylene Middle Filter Layer: Melt-Blown Polypropylene
STORAGE	Product Shall Be Stored Under Room Condition, Avoid Direct Sunlight
SHELF LIFE	2 Years From Date Of Manufacture With The Above Storage Condition

TOL GLOBAL HEALTH PACKAGING DIMENSIONS

