


# Tested and Approved: Level 1 Face Mask

## MEDICAL GRADE LEVEL ONE 3-PLY MASK



Brand / Model: TOL GLOBAL HEALTH L1 MEDICAL 3-PLY Adult Size Mask / Non-Sterile / 50 Per Box

Case Pack: 40  
Gross Weight: 19 lbs  
Pallet Count: 11 x 6

EN 14683:2019 + AC : 2019(E) ANNEX C

**LEVEL 1 PROTECTION**

**NO LATEX - NON IRRITATING**

**BACTERIAL FILTRATION EFFICIENCY ≥ 95%**

**MEDICAL GRADE**

### PRODUCT SPECIFICATIONS

### PRODUCT UPC

ASTM F2101-19  
ASTM F1862  
ASTM F1862M-17  
ASTM F2100-19e1  
ASTM F2299

✓ 16 CFR PART 1610-2012

WHITE

8 46146 05045 1

BLUE

8 46146 05046 8

BLACK

8 46146 05047 5

PINK

8 46146 05048 2

<b>NAME</b>	TOL GLOBAL HEALTH Medical Grade Level 1 Face Mask
<b>DESIGN</b>	Face Mask With Ear Loops
<b>SURFACE</b>	Smooth
<b>COLOR</b>	White, Blue, Black, Pink. Custom Colors & Designs Available
<b>MATERIAL</b>	Inner & Outer layer: Spun-bound Polypropylene Middle Filter Layer: Melt-Blown Polypropylene
<b>STORAGE</b>	Product Shall Be Stored Under Room Condition, Avoid Direct Sunlight
<b>SHELF LIFE</b>	2 Years From Date Of Manufacture With The Above Storage Condition

### TOL GLOBAL HEALTH PACKAGING DIMENSIONS

