


HY106

- *Lightweight construction aids in helping caregivers avoid back injuries. Hydraulic pump handle can rotate from side to side for caregiver convenience.
- *Easy disassembly for transport and set-up
- *780mm to 1740mm range allows patient to be lifted from the floor
- *Locking casters available
- *Economical with a 180 kg weight capacity



90 cm



TECHNICAL DATA

1200 mm



Length

620 mm



Width

860 mm



Max fork Distance

900 mm



Lifting Distance

39 Kg



Net Weight

180 Kg



Max Loading Weight