








# HY118

The Easy HY118 is the perfect choice when a strong, sturdy floor lift is required. It features a lighter Aluminum which makes maneuvering it extremely easily, even in tighter spaces. The legs open and close electronically and the compact size make it ideal for use in both a home or a facility.



## TECHNICAL DATA

1424 mm  Length	656 mm  Width	931 mm  Max fork Distance
1740 mm  Lifting Distance	65 mm  Ground Clearance	24V 7000n  Motor
60-80 Times  Battery Capacity	38 Kg  Net Weight	180 Kg  Max Loading Weight