

HS480

Stable, Trusty, Durable: With a focus on safety and reliability, the HS480 scooter is equipped with multiple features to ensure a safe and easy ride. Imagine being able to undertake long journeys in complete safety and comfort. Whether you enjoy a daily trip to the shops or love to travel further afield, the HS480 scooter gives you freedom of independence. Combining ergonomics with style, HS480 scooter is the obvious choice for a relaxing, reliable drive.

Easy Disassembly

HS480 Scooters provide fast and simple disassembly. From our compact travel scooters to our larger luxury models, innovative design allow for convenient storage and easy transport.



TECHNICAL DATA

360 mm Backrest Height	430 mm Seat Width	410 mm Seat Depth	580 mm Driving Unit Width	1200mm Overall Length	500-560 mm Ground Clearance of Seat
65 mm Ground Clearance	60 mm Spanning Ability	12° Max. Safe Slope	25 Km Travel Distance	1500 mm Turning Radius	280 W Motor
22 Ahx2 Lead acid battery Battery Capacity (Option)	24V 2.0A Charger	19 Kg The Heaviest Part	59 Kg Net Weight	135 Kg Weight Capacity	7 Km/h Max. Speed
9" Solid PU wheels Wheels					