

# HS280

## Easy Disassembly

HS280 Scooters provide fast and simple disassembly. From our compact travel scooters to our larger luxury models, innovative design allow for convenient storage and easy transport.

## Easy Portability

Our scooters feature easy to use, excellent controls to make operation smooth and worry free.

## Elegant Style

Our scooters are designed to be a visually striking as they are high performing. Our sporty designs allow your scooter to be an extension of your unique personal style.



## TECHNICAL DATA

320 mm	440 mm	390 mm	580 mm	1200mm	500-560 mm
Backrest Height	Seat Width	Seat Depth	Driving Unit Width	Overall Length	Ground Clearance of Seat
50 mm	35 mm	10°	15 Km	1350 mm	180 W
Ground Clearance	Spanning Ability	Max. Safe Slope	Travel Distance	Turning Radius	Motor
12 Ahx2 Lead acid battery	24V 2.0A	17.4 Kg	45 Kg	125 Kg	7 Km/h
Battery Capacity (Option)	Charger	The Heaviest Part	Net Weight	Weight Capacity	Max. Speed
8" Solid PU front wheels 9" Solid PU rear wheels					
Wheels					