

HS300

In Performance

HS300 Scooters are engineered and factory tested high performance operation that you can count on. Our high quality construction and proven durability means you can depend on it for a long time.



TECHNICAL DATA

320 mm	440 mm	390 mm	570 mm	1090 mm	500-560 mm
Backrest Height	Seat Width	Seat Depth	Driving Unit Width	Overall Length	Ground Clearance of Seat
75 mm	35 mm	10°	25 Km	1400 mm	180 W
Ground Clearance	Spanning Ability	Max. Safe Slope	Travel Distance	Turning Radius	Motor
17 Ahx2 Lead acid battery	24V 2.0A	32.3 Kg	52 Kg	125 Kg	7 Km/h
Battery Capacity (Option)	Charger	The Heaviest Part	Net Weight	Weight Capacity	Max. Speed
9" Solid PU wheels					
Wheels					