



BestShape[®] Care
AI Non-contact Care Series

BestShape Care -The Main Functions & Features

User

- No-image / Non-wearable
- Ease & comfortable

Caregiver

- 24/7 continuous monitoring
- Auto notification (Line/WhatsApp)



BestShape Care _ Competitive Analysis

	Technology and application	Alarm	Function	Strengths and weaknesses
BestShape Care Non-contact care system	mmWave Radar Motion	<ul style="list-style-type: none"> Notification and Reminder Personal settings SNS notification 	<ul style="list-style-type: none"> Location detection Prolonged motionless warning Not-in-bed warning 	<ul style="list-style-type: none"> 24/7 considerate care No image and no wearable, easy to clean Whole space care (connect multi-devices) Water resistance and can be used in the bathroom
Pressure sensing Mattress / Vibration sensor	Fiber-optic sensor Motion	Abnormal warning	<ul style="list-style-type: none"> Not-in-bed warning 	<ul style="list-style-type: none"> High cost Should be cleaned frequently
Wearable device	Motion sensor Motion	Abnormal warning	<ul style="list-style-type: none"> Not-in-bed warning 	<ul style="list-style-type: none"> Should be cleaned frequently forget to wear
Thermal camera/ Video monitor	Camera Motion	Abnormal warning	<ul style="list-style-type: none"> Real-time and recorded video 	<ul style="list-style-type: none"> Privacy issues Ambient light and temperature trigger false alarms Blind or Block Cameras

Application - Sites

Accessible toilets



Single wards



Senior apartments


























Two types of Real-time status - Senior apartments

No.	Sub No.	Name	Bed status	HR	RR	Room status				
						Bedroom	Livingroom	Kitchen	Toilet	Not in room
204	01	零〇零		--	--					
217	01	零〇六		--	--					
310	01	零〇二		--	--	Pause reminder				
403	01	零〇六		--	--	Pause reminder				
503	01	零〇七		85	18					
507	01	零〇八		--	--	Pause reminder				
606	01	零〇八		73	12	Pause reminder				
607	01	零〇六		--	--	Pause reminder				
909	01	零〇三		66	22					
911	01	零〇九		62	15	Pause reminder				
1003	01	零〇一		63	11	Pause reminder				
1011	01	零〇一		70	18					
1103	01	零〇六		--	--	Pause reminder				
1106	01	零〇五		68	18	Pause reminder				
1205	01	零〇零		--	--	Pause reminder				

This page shows 1-15 records, a total of 15 records.

Sleep + motion

Dashboard – Status and Icons

	Bedroom	Living Room	Kitchen	Toilet	Not in room
Occupied					
Alarm					
Fall					
Disconnect					
Prolong motionless					
No device					
Vacant	