

Drawback of the traditional Nasogastric Tube fixation method (1)

Frequently using adhesive tape to fix the Nasogastric Tube in the same position on the skin can cause skin tearing injuries due to inadequate rest time for the skin in the nasal area.



Drawback of the traditional Nasogastric Tube fixation method (2)

Using adhesive tape to fix the Nasogastric Tube creates a pressure point in the nasal cavity, which can lead to nasal mucosal damage due to long-term compression of small blood vessels.



Nasogastric Tube Fixation Patch



Product feature(1)

- ◆ Changing the adhesive position → **Reduces the risk of skin tearing.**



- ◆ Suspending the tube in the nasal cavity




→ **Reduces the risk of nasal pressure injury.**

We offer two kind of Nasogastric Tube:

1. Two Pieces Nasogastric Tube
2. Connector

Nasogastric Tube Combination Connector

【 A two-piece Nasogastric Tube 】 can be customized to meet the individual needs of patients.




When wearing a mask, the Nasogastric Tube is not visible, enhancing the patient's appearance and preserving their dignity.



| | |
|---|---|
| 1 | Customizing the one-piece Nasogastric Tube into a two-piece tube based on individual patient height provides cost-effective access to the benefits of a two-piece tube, with the patented anti-reflux valve design effectively resolving the limitations of existing one-piece Nasogastric Tubes. |
| 2 | The exposed Nasogastric Tube is shorter and lighter, reducing the risk of dislodgement from patient movement or accidental pulling. |
| 3 | The Nasogastric Tube does not protrude outside the mask, enhancing the patient's appearance and preserving their dignity. |
| 4 | Long-term use of a two-piece Nasogastric Tube provides greater patient comfort due to its lighter weight. |

REDUCE-BACKFLOW DETACHABLE NASOGASTRIC TUBE



- **Reduce the chance of pulling out the Nasogastric Tube**
- **Maintain patient dignity**
- **Enable easy feeding**
- **Reduce food reflux from Nasogastric Tube**