



# L I F E C A R E H E A R I N G C E N T R E

2190 Warden Avenue (at Sheppard)  
Scarborough, Ontario M1T 1V6

Tel.: (416) 496-3200  
Fax: (416) 496-3833  
lifecarehearingcentre@gmail.com

## TIPS – for Rechargeable Hearing aids

### Please clean hearing aid weekly.

RIGHT HEARING AID - **RED** or **R**

LEFT HEARING AID - **BLUE** or **L**



- Please charge hearing aids every night
- Hold bottom button for 5 seconds to turn on/off
- Charging hearing aids will automatically turn it off
- Green light while in charger = fully charged
- Red light while in charger = charging
- Press upper button to increase volume
- Press lower button to decrease volume

---

## Hearing Aid Not Working?

- Charge the hearing aid when you do not hear anything for at least 1 hours
- Look for wax buildup on the hearing aid and clean with the pick/brush

For help Videos: [www.lifecarehearingcentre.ca/hearing-aid-help](http://www.lifecarehearingcentre.ca/hearing-aid-help)

Call **(416)496-3200** for help, or to book an Appointment