



American Medical Equipment, Inc.

29450 Gratiot Ave, Roseville, MI 48066

PH: (586) 298-6161 FX: (866) 593-7841




LIFT CHAIRS HERITAGE COLLECTION MADE BY PRIDE MEDICAL



PETITE LC-358P 3-Position, Full Recline, Chaise Lounge	PETITE WIDE LC-358PW 3-Position, Full Recline, Chaise Lounge	SMALL LC-358S 3-Position, Full Recline, Chaise Lounge	MEDIUM LC-358M 3-Position, Full Recline, Chaise Lounge	LARGE LC-358L 3-Position, Full Recline, Chaise Lounge	X-LARGE LC-358XL 3-Position, Full Recline, Chaise Lounge	XX-LARGE LC-358XXL 2-Position, Partial Recline, Chaise Lounge
Specifications: <ul style="list-style-type: none"> Weight Capacity: 375 lbs. User Height Range: 5'2" & Below Seat to Floor: 18" Seat Depth: 18" Seat Width: 18" Top of Back to Seat: 24" Back Style: Button Back Fully Padded Chaise: Yes Recline Positions: 3-Position Heat & Massage Option: Yes Distance From Wall: 17.5" Overall Chair Width: 34.5" Chair Weight: 112 lbs. 	Specifications: <ul style="list-style-type: none"> Weight Capacity: 375 lbs. User Height Range: 5'3" & Below Seat to Floor: 19" Seat Depth: 19" Seat Width: 24" Top of Back to Seat: 23.5" Back Style: Button Back Fully Padded Chaise: Yes Recline Positions: 3-Position Heat & Massage Option: Yes Distance From Wall: 17" Overall Chair Width: 39" Chair Weight: 118 lbs. 	Specifications: <ul style="list-style-type: none"> Weight Capacity: 375 lbs. User Height Range: 5'3" & Below Seat to Floor: 19" Seat Depth: 19" Seat Width: 18" Top of Back to Seat: 24" Back Style: Button Back Fully Padded Chaise: Yes Recline Positions: 3-Position Heat & Massage Option: Yes Distance From Wall: 17" Overall Chair Width: 35" Chair Weight: 118 lbs. 	Specifications: <ul style="list-style-type: none"> Weight Capacity: 375 lbs. User Height Range: 5'4" - 6'0" Seat to Floor: 20" Seat Depth: 20" Seat Width: 20" Top of Back to Seat: 24.5" Back Style: Button Back Fully Padded Chaise: Yes Recline Positions: 3-Position Heat & Massage Option: Yes Distance From Wall: 17" Overall Chair Width: 35" Chair Weight: 117 lbs. 	Specifications: <ul style="list-style-type: none"> Weight Capacity: 375 lbs. User Height Range: 5'9" - 6'2" Seat to Floor: 21" Seat Depth: 22" Seat Width: 22" Top of Back to Seat: 25.5" Back Style: Button Back Fully Padded Chaise: Yes Recline Positions: 3-Position Heat & Massage Option: Yes Distance From Wall: 17.5" Overall Chair Width: 37" Chair Weight: 120 lbs. 	Specifications: <ul style="list-style-type: none"> Weight Capacity: 500 lbs. User Height Range: 5'6" - 6'1" Seat to Floor: 20" Seat Depth: 22" Seat Width: 26" Top of Back to Seat: 29" Back Style: Button Back Fully Padded Chaise: Yes Recline Positions: 3-Position Heat & Massage Option: Yes Distance From Wall: 18" Overall Chair Width: 39.5" Chair Weight: 159 lbs. 	Specifications: <ul style="list-style-type: none"> Weight Capacity: 600 lbs. User Height Range: 5'6" - 6'1" Seat to Floor: 20" Seat Depth: 22" Seat Width: 30" Top of Back to Seat: 29" Back Style: Button Back Fully Padded Chaise: Yes Recline Positions: 2-Position Heat & Massage Option: Yes Distance From Wall: 18" Overall Chair Width: 43" Chair Weight: 168 lbs.
Retail Price: \$1165	Retail Price: \$1285	Retail Price: \$1165	Retail Price: \$1230	Retail Price: \$1230	Retail Price: \$1,570	Retail Price: \$1,800

**Pricing above could change due to manufacturer pricing increases

Color Options:



** IMPORTANT **

Please be sure to choose the appropriate chair for your needs because chairs are non-returnable & non-refundable

*Colors on printed paper may vary slightly on chair