











Article	Description
ZBH-4D	Tool holder + fixation screw
	standard delivered, 3 pieces
31199701AS	Fixation screw (20 Pcs)
	
1000009895	Detection pin
	standard delivered
1000010985	Detection jig
	standard delivered
ZS-4D	Spindle replacement kit
	


ZSB-140	Spindle belt
	


1000015824	Spanner	standard delivered
		

ZAT-4D	4 position ATC
	


ZV-4D	Multiple slot clamp
	

6701609030	Maintenance grease (10ml)	standard delivered
		


21675102	Maintenance grease for collet (10ml)
	

ZC-4D	Collet replacement kit (Collet grease 21675102 included)
	


Extraction

D0844A8201	Dust Pro Universal - DWXrange
A1030170	Replacement filter bag - DWX
A1030161	Replacement HEPA filter - DWX
	

Software


Millbox-4	Roland Dental CAM software for DWX-4
	

ExoCAD	Exocad DentalCam standard version
ExoCAD modules	
Exoplan viewer	Exocad dicom viewer
Exoimplant	Exocad implant module
Exomodel creator	Exocad creator module
Exoprovisional	Exocad provisional crown creator
Articulator	Exocad virtuell articulator
TruSmile Techn.	Exocad TruSmile technology
Tooth library	Exocad tooth library



Milling tools


Ball end cutter




ZCB-40D	R0,4 x 2(l) x 10(EL) x 50(L) x 4(D) x 2NT
ZCB-50D	R0,5 x 2(l) x 12(EL) x 50(L) x 4(D) x 2NT
ZCB-100D	R1,0 x 10(l) x 16(EL) x 50(L) x 4(D) x 2NT
ZDB-30D	R0,3 x 2(l) x 8(EL) x 50(L) x 4(D) x 2NT
ZDB-50D	R0,5 x 1,5(l) x 10(EL) x 50(L) x 4(D) x 2NT
ZDB-100D	R1,0 x 3(l) x 16(EL) x 50(L) x 4(D) x 2NT




SF06-10	0,6(R) x 10(l) x 50(L) x 4(d) x 2NT
SF1-20	1,0(R) x 20(l) x 50(L) x 4(d) x 2NT
SF2-20	2,0(R) x 20(l) x 50(L) x 4(d) x 2NT
SF3-20	3,0(R) x 20(l) x 50(L) x 4(d) x 2NT
SF1-12	1,0(R) x 12(l) x 50(L) x 4(d) x 2NT
SF2-16	2,0(R) x 16(l) x 50(L) x 4(d) x 2NT

	
04.003.11DIA	Dental milling bur, ball end D=0.3mm, l = 2 mm
04.006.11DIA	Diamond coated bur, ball end D=0.6mm, l = 11mm
04.010.18DIA	Diamond coated bur, ball end D=1.0mm, l = 18mm
04.020.20DIA	Diamond coated bur, ball end D=2.0mm, l = 20mm

Square end cutter

	
PT06-10	0,6(D) x 10(l) x 50(L) x 4(d) x 2NT
PT1-20	1,0(D) x 20(l) x 50(L) x 4(d) x 2NT
PT2-20	2,0(D) x 20(l) x 50(L) x 4(d) x 2NT

Lollipop

	
Lollipop	R1,0 x 14(l) x 50(L) x 4(D)

Scanners

A855	Scanner activity 855
A885	Scanner activity 885



9002	Fixator system Artex	For Activity A885 scanner
9003	Fixator system Adesso Split	For Activity A885 scanner
9004	Fixator system SAM	For Activity A885 scanner



Ovens

HT-S	High Temperature Furnace HT-S
HT-S-SPEED	High Temperature Furnace HT-S SPEED
TABEO	High Temperature Furnace TABEO



Milling material

Zirconium



252455	Zirkon BioStar HT Smile, 76x40x16mm
252456	Zirkon BioStar HT Smile, 76x40x22mm
R-76-16-NT-WT-02	Nexx Zirkon Translucent 76x40x16mm, 2PC
R-76-22-NT-WT-02	Nexx Zirkon Translucent 76x40x22mm, 2PC



252459	Zirkon BioStar HT S, 40x19x15,5mm -10PC	* ZV-4D needed
R-40-19-NT-WT-05	Nexx Zirkon Translucent 40x19x15mm, 5PC	* ZV-4D needed
R-20-15-NT-WT-10	Nexx Zirkon Translucent 20x15x14mm, 10PC	* ZV-4D needed
R-20-19-NT-WT-10	Nexx Zirkon Translucent 20x19x15mm, 10PC	* ZV-4D needed

Wax



R-76-16-WX-BL-02	Sagemax Milling wax 76x40x16mm, 2PC
R-76-16-WX-WT-02	Sagemax Milling wax 76x40x16mm, 2PC
R-76-22-WX-BL-02	Sagemax Milling wax 76x40x22mm, 2PC
R-76-22-WX-WT-02	Sagemax Milling wax 76x40x22mm, 2PC