

Prescribing Information:

Federal and State Worker's Compensation:

Prescription stating: Neurolumen Daily
Diagnosis Codes
Patient Demographics
Patient's Case File Number

Third Party Insurance:

Prescription stating: Neurolumen Daily
Diagnosis Codes
Patient Demographics
Copy of Insurance Card (front & Back)

Veterans Administration: Over 50 Facilities
are prescribing

Approved through the Central Office
FSS/GSA : V797D-40311
Vendor: Neurolumen LLC Prime
Prime Vendor Contract:
AMD: 999269626

Tax ID: 74-3260316
SAMS: 5VJR9
CPT Code: E0761
DUNS: 957741148
FDA: K142025

Prosthetics Consult suggestion:

Neurolumen PN-1000-1
Patient has tried progressive evolving
therapeutic methods and requires
Neurolumen: 12 Lasers, 24 LEDs, 8
electrical stimulation electrodes.

Fax to 800-547-1296

" I was a paratrooper in Vietnam and broke my back in a jump into hostile territory. After many spinal and neck surgeries, I was in a wheel chair for a long time which resulted in two rotator cuff surgeries. I've had fibromyalgia for as long as I can remember and have tried just about every therapy to give me some comfort and mobility. I just hurt everywhere. I started using the Neurolumen system for pain and was amazed! For the first time in years I didn't hurt" -L.P.

"I injured my right shoulder March of 2011 which led to a more recent diagnosis of arthropathy of the AC joint along with PTSD, after years of physical therapy and pain medications I finally found a real solution for my pain. With in a few treatments my pain levels went from an 8 to a 2! My PTSD symptoms have also drastically improved! With my physicians help I went from 11 different medications to now just 3! My quality of life has drastically improved, I am actually able to sleep at night and play with my children." -C.C.

"What a relief to finally find help for my chronic foot, knee, back and neck pain. After many years of suffering, my pain is finally becoming tolerable after only a few Neurolumen treatments. The treatments have eased the pain considerably. My family and coworkers have noticed a positive difference in my posture and demeanor after starting the treatments. I had abandoned exercise because of the extra pain, and am now gradually working back into a more active life. What a blessing after thinking for so long that no help was available except medication." -S.P.

855-855-4648

Your Pain Stops Here



NEUROLUMEN



Pain relief when you need it!

405-463-6525



TREATABLE AREAS:

- BACK
- NECK
- SHOULDERS
- HIPS
- KNEES
- FOOT
- ANKLE
- HAND
- ELBOW
- WRIST

RELIEVE DISCOMFORT IN THE FOLLOWING CONDITIONS:

- ARTHRITIS
- NEUROPATHY
- CHRONIC/ACUTE PAIN
- LYMPHEDEMA
- WOUND CARE
- MULTIPLE SCLEROSIS
- SCIATICA
- PAD
- CVI (VENOUS INSUFFICIENCY)
- DRUG ADDICTION
- SCOLIOSIS
- SHINGLES
- TRIGEMINAL NEURALGIA

**Three Proven Technologies,
One Powerful Therapy**

PAIN RELIEF WHEN AND WHERE YOU NEED IT!

Neurolumen utilizes low-level LASERS, electrical stimulation and LEDs combined into an easy to apply flexible wrap. The restorative power of the light therapy not only decreases inflammation, but also increases circulation to offer an effective alternative for your pain relief.

Neurolumen has taken the most effective wavelengths of light and strategically placed them to provide the best therapeutic benefit.

Treatments with this device will administer light therapy while simultaneously delivering electrical stimulation. This will create muscle movement that will assist the removal of toxins in the area while delivering nutrients more effectively.

The combination of these powerful technologies come together to create a synergy that proves to be more effective than each modality alone.

Technical Specifications

Laser Diodes

Wavelength..... 808 nm
Output Power 60mW

Electrical Stimulation

Frequency..... 2 Hz
Output Current 100 mA max
Pulse Voltage..... 75VDC max
Pulse Width 100us biphasic

LEDs RED

Wavelength..... 660 nm
Output Power 15 mW

LEDs IR

Wavelength..... 904 nm
Output Power 22 mW



- Battery operated to travel with you anywhere.
- Hands-free
- Use directly over hardware
- Up to 12 Lasers, 24 LEDs and 8 electrical stimulation pads per unit.
- You can treat up to six areas at one time.
- Insurance coverage.
- Case included
- Gel pads last up to three months
- Can be modified to be used with pacemaker
- Cleared for home use.
- Training provided

We are here to help! Please call us today to find out more information.



855.855.4648 405.463.6525
Fax: 800.547.1296

9632 N May Avenue
Oklahoma City, OK 73120
www.neurolumen.com