

OneKnot® Polyester is a sterile single-use closure system with a pre-made knot. Contains polyester suture with a silicone coating of medical silicone elastomer MED2174. Nonabsorbable multifilament, braided (green). Stainless steel needle AISI 302.

OneKnot Polyester				
USP	Needle type	Length	Unit/Box	Product REF
0	17mm Reverse cutting 1/2 circle	65 cm	12 devices	1RH0CP
0	17mm Reverse cutting 1/2 circle	90 cm	12 devices	1RH0AP
0	17mm Reverse cutting 1/2 circle	120 cm	12 devices	1RH0BP
0	22mm Taperpoint 1/2 circle	65 cm	12 devices	3HH0CP
0	22mm Taperpoint 1/2 circle	90 cm	12 devices	3HH0AP
0	22mm Taperpoint 1/2 circle	120 cm	12 devices	3HH0BP
0	22mm Reverse cutting 1/2 circle	65 cm	12 devices	3RH0CP
0	22mm Reverse cutting 1/2 circle	90 cm	12 devices	3RH0AP
0	22mm Reverse cutting 1/2 circle	120 cm	12 devices	3RH0BP
0	26mm Taperpoint 1/2 circle	65 cm	12 devices	5HH0CP
0	26mm Taperpoint 1/2 circle	90 cm	12 devices	5HH0AP
0	26mm Taperpoint 1/2 circle	120 cm	12 devices	5HH0BP
0	26mm Reverse cutting 1/2 circle	65 cm	12 devices	5RH0CP
0	26mm Reverse cutting 1/2 circle	90 cm	12 devices	5RH0AP
0	26mm Reverse cutting 1/2 circle	120 cm	12 devices	5RH0BP
0	26mm Tapercut 1/2 circle	65 cm	12 devices	5YH0CP
0	26mm Tapercut 1/2 circle	90 cm	12 devices	5YH0AP
0	26mm Tapercut 1/2 circle	120 cm	12 devices	5YH0BP
0	36mm Taperpoint 1/2 circle	65 cm	12 devices	6HH0CP
0	36mm Taperpoint 1/2 circle	90 cm	12 devices	6HH0AP
0	36mm Taperpoint 1/2 circle	120 cm	12 devices	6HH0BP
0	36mm Reverse cutting 1/2 circle	65 cm	12 devices	6RH0CP
0	36mm Reverse cutting 1/2 circle	90 cm	12 devices	6RH0AP
0	36mm Reverse cutting 1/2 circle	120 cm	12 devices	6RH0BP
0	40mm Reverse cutting 1/2 circle	65 cm	12 devices	7RH0CP
0	40mm Reverse cutting 1/2 circle	90 cm	12 devices	7RH0AP
0	40mm Reverse cutting 1/2 circle	120 cm	12 devices	7RH0BP

OneKnot Polyester				
USP	Needle type	Length	Unit/Box	Product REF
2/0	17mm Reverse cutting 1/2 circle	65 cm	12 devices	1RH2CP
2/0	17mm Reverse cutting 1/2 circle	90 cm	12 devices	1RH2AP
2/0	17mm Reverse cutting 1/2 circle	120 cm	12 devices	1RH2BP
2/0	22mm Taperpoint 1/2 circle	65 cm	12 devices	3HH2CP
2/0	22mm Taperpoint 1/2 circle	90 cm	12 devices	3HH2AP
2/0	22mm Taperpoint 1/2 circle	120 cm	12 devices	3HH2BP
2/0	22mm Reverse cutting 1/2 circle	65 cm	12 devices	3RH2CP
2/0	22mm Reverse cutting 1/2 circle	90 cm	12 devices	3RH2AP
2/0	22mm Reverse cutting 1/2 circle	120 cm	12 devices	3RH2BP
2/0	26mm Tapercut 1/2 circle	65 cm	12 devices	5YH2CP
2/0	26mm Tapercut 1/2 circle	90 cm	12 devices	5YH2AP
2/0	26mm Tapercut 1/2 circle	120 cm	12 devices	5YH2BP
2/0	26mm Taperpoint 1/2c	65 cm	12 devices	5HH2CP
2/0	26mm Taperpoint 1/2c	90 cm	12 devices	5HH2AP
2/0	26mm Taperpoint 1/2c	120 cm	12 devices	5HH2BP
2/0	36mm Tapercut 1/2 circle	65 cm	12 devices	6YH2CP
2/0	36mm Tapercut 1/2 circle	90 cm	12 devices	6YH2AP
2/0	36mm Tapercut 1/2 circle	120 cm	12 devices	6YH2BP
2/0	36 mm Reverse cutting 1/2 circle	65 cm	12 devices	6RH2CP
2/0	36 mm Reverse cutting 1/2 circle	90 cm	12 devices	6RH2AP
2/0	36 mm Reverse cutting 1/2 circle	120 cm	12 devices	6RH2BP
2/0	36 mm Reverse cutting 5/8 circle	65 cm	12 devices	6RF2CP
2/0	36 mm Reverse cutting 5/8 circle	90 cm	12 devices	6RF2AP
2/0	36 mm Reverse cutting 5/8 circle	120 cm	12 devices	6RF2BP
2/0	40mm Reverse cutting 1/2 circle	65 cm	12 devices	7RH2CP
2/0	40mm Reverse cutting 1/2 circle	90 cm	12 devices	7RH2AP
2/0	40mm Reverse cutting 1/2 circle	120 cm	12 devices	7RH2BP

OneKnot Polyester				
USP	Needle type	Length	Unit/Box	Product REF
3/0	22mm Tapercut 1/2 circle	65 cm	12 devices	3YH3CP
3/0	22mm Tapercut 1/2 circle	90 cm	12 devices	3YH3AP
3/0	22mm Tapercut 1/2 circle	120 cm	12 devices	3YH3BP
3/0	26mm Tapercut 1/2 circle	65 cm	12 devices	5YH3CP
3/0	26mm Tapercut 1/2 circle	90 cm	12 devices	5YH3AP
3/0	26mm Tapercut 1/2 circle	120 cm	12 devices	5YH3BP
3/0	26mm Taperpoint 1/2 circle	65 cm	12 devices	5HH3CP
3/0	26mm Taperpoint 1/2 circle	90 cm	12 devices	5HH3AP
3/0	26mm Taperpoint 1/2 circle	120 cm	12 devices	5HH3BP
3/0	26mm Reverse cutting 3/8 circle	65 cm	12 devices	5RE3CP
3/0	26mm Reverse cutting 3/8 circle	90 cm	12 devices	5RE3AP
3/0	26mm Reverse cutting 3/8 circle	120 cm	12 devices	5RE3BP
3/0	36mm Tapercut 1/2 circle	65 cm	12 devices	6YH3CP
3/0	36mm Tapercut 1/2 circle	90 cm	12 devices	6YH3AP
3/0	36mm Tapercut 1/2 circle	120 cm	12 devices	6YH3BP
4/0	17mm Tapercut 3/8 circle	65 cm	12 devices	1YE4CP
4/0	17mm Tapercut 3/8 circle	90 cm	12 devices	1YE4AP
4/0	17mm Tapercut 3/8 circle	120 cm	12 devices	1YE4BP
4/0	24mm Reverse cutting 1/2 circle	65 cm	12 devices	4RH4CP
4/0	24mm Reverse cutting 1/2 circle	90 cm	12 devices	4RH4AP
4/0	24mm Reverse cutting 1/2 circle	120 cm	12 devices	4RH4BP
4/0	26mm Tapercut 1/2 circle	65 cm	12 devices	5YH4CP
4/0	26mm Tapercut 1/2 circle	90 cm	12 devices	5YH4AP
4/0	26mm Tapercut 1/2 circle	120 cm	12 devices	5YH4BP
4/0	26mm Reverse cutting 1/2 circle	65 cm	12 devices	5RH4CP
4/0	26mm Reverse cutting 1/2 circle	90 cm	12 devices	5RH4AP
4/0	26mm Reverse cutting 1/2 circle	120 cm	12 devices	5RH4BP

OneKnot Polyester				
USP	Needle type	Length	Unit/Box	Product REF
5/0	17mm Taperpoint 1/2 circle	65 cm	12 devices	1HH5CP
5/0	17mm Taperpoint 1/2 circle	90 cm	12 devices	1HH5AP
5/0	17mm Taperpoint 1/2 circle	120 cm	12 devices	1HH5BP
5/0	17mm Tapercut 3/8 circle	65 cm	12 devices	1YE5CP
5/0	17mm Tapercut 3/8 circle	90 cm	12 devices	1YE5AP
5/0	17mm Tapercut 3/8 circle	120 cm	12 devices	1YE5BP
5/0	19mm Reverse cutting 3/8 circle	65 cm	12 devices	2RE5CP
5/0	19mm Reverse cutting 3/8 circle	90 cm	12 devices	2RE5AP
5/0	19mm Reverse cutting 3/8 circle	120 cm	12 devices	2RE5BP

EU MDR CE Class IIb medical product. Intended for use by healthcare professionals only. For indications, contraindications, and technical specifications, consult Instructions for Use (IFU).

Manufactured in Germany.  
© 2026 Clouz GmbH. All rights reserved.