

SHARK™ 

®VIAL HOLDER
FOR INJECTABLE PRODUCTS



FEATURES

A SLEEK DESIGN COMBINING FUNCTIONALITY AND EFFICIENCY

NARROW CLAMP HOLDER FOR
BETTER VISIBILITY

FLEXIBLE, REMOVABLE VIAL CLAMPS

LOCKABLE ARTICULATED HEAD

LOCKABLE TELESCOPIC TUBE

RIGID ANODIZED ALUMINUM TUBE

VERSATILE INSTALLATION BASE



FEATURES



- FLEXIBLE, REMOVABLE VIAL CLAMPS



SUPPORT HEAD FOR 6 CLAMPS

CLAMP AVAILABLE IN DIFFERENT
MODELS (A, B, C, D, E)


COMPATIBLE WITH MOST STANDARD
VIALS FROM 1 TO 100 ML



QUICK RELEASE BUTTON

FEATURES


- LOCKABLE ARTICULATED HEAD



FEATURES

- LOCKABLE TELESCOPIC TUBE

LOCKING SLEEVE




TELESCOPIC TUBE (580 - 980 MM)




FEATURES

- VERSATILE INSTALLATION BASE

STANDARD BASE FOR FLAT SURFACE OR
SQUARE TUBE (MAX. THICKNESS 80 MM)




ROUND TUBE ADAPTER (Ø 22 - 32 MM)



FEATURES

- US/EU RAIL COMPATIBLE MOUNT



NX MED

NEXT INNOVATION. NEXT SOLUTION.



[HTTPS://WWW.LINKEDIN.COM/COMPANY/NXMED/](https://www.linkedin.com/company/nxmed/)



[HTTPS://WWW.INSTAGRAM.COM/NXMED.INC/](https://www.instagram.com/nxmed.inc/)

