

Clinical Lead Wire Compatibility Chart

Manufacturer	Model	Unit made prior to 2000		Complaint Unit Wires
		Adaptor	Wire	
AMREX	HV752	AD4500	W47 (Banana-Pin)	W47 (Banana-Pin)
	LVG325A	AD4500	W47 (Banana-Pin)	W47 (Banana-Pin)
	MS SERIES (ALL)	AD4500	W47 (Banana-Pin)	W47 (Banana-Pin)
	Z-STIM IF SERIES	AD4500	W47 (Banana-Pin)	W47 (Banana-Pin)
BETTER LIFE	AP-439 PLUS	NA	NA	W80 (5-Pin DIN-Pin)
CHATTANOOGA	CPS 200 Combo, CPS 200 Stim	No CMP Solution for USA	W18112 (CDN Only)	W71 (4-Pin DIN)
	CPS 400 Combo, CPS 400 Stim	No CMP Solution for USA	W18112 (CDN Only)	W71 (4-Pin DIN)
	Forte 200, Forte 400 Combo	No CMP Solution for USA	W18112 (CDN Only)	W71 (4-Pin DIN)
	Forte CB 450, Forte ES450	AD8500	W65 (Stereo-Pin)	
	Intelect 120, 170, 775MP, 600, 900	AD8500	W65 (Stereo-Pin)	
	Intelect 500, 500S, 700, 700C, 775		W2372 (3 Pin Amphenol)	
	Intelect 550, 551	AD4500	W44 (Banana-Pin)) W43 (Banana-Pin)	
	Intelect Combo, Legend Stim	No CMP Solution for USA	W18112 (CDN Only)	W71 (4-Pin DIN)
	Intelect High Volt Stim			W8860BLBRB (Mini 3-Pin DIN-Pin)
	Intelec IFC2 (squared)	AD8500	W65 (Stereo-Pin)	
	Intelect Legend XT (#2760), Vectra Genisys, Transport, US 2738			W176 (2 Cond Socket-Pin)
	Intelect Low Volt 110			W79 (4-Pin DIN-Pin)
	Low Volt 120			W71 (4-Pin DIN-Pin)
	Vectra Pro, Vectra			W71 (4-Pin DIN-Pin)
	VMS		W2372 (3-Pin Amphenol)	
	VMS II	AD8500	W65 (Stereo-Pin)	
DYNATRONICS	Dynatron 50, 125, 525	AD6500	W66 (Gold Stereo-Pin)	W66 (Gold Stereo-Pin)
	Dynatron Plus 550, 650, 850, 950	AD6500	W66 (Gold Stereo-Pin)	W66 (Gold Stereo-Pin) W62 (Gold Stereo combo-Gold Banana)
	Dynatron Solaris 708	1X W4772IVRR + 2X AD4500 and 1X W4772IVBB + 2X AD4500		
DIMEQ				W72 (4-Pin DIN-Pin)
ENRAF	Endomed 182, 482			W77 (4-Pin DIN-Pin)
	Endomed 582, 982			W83 (R/A 4-Pin DIN-Pin W82 (R/A 4-Pin-Banana)
ENRAF(HENLEY)	Endomed 433	AD8400	W84 (5-Pin DIN-Pins)	
EXCEL	Ultra/Multi/Isotron Series	AD6500	W65 (Stereo-Pin)	W65 (Stereo-Pin)
FISHER	Model 700	AD4500	W47 (Banana-Pin)	
H WAVE		A4500	W42 (Banana-Banana) W43 (Banana-Pin)	
HAKOMED				W86 (5-Pin DIN-Pin)
HEALTHCARE MNFC	Masterstim	AD6500	W65 (Stereo-Pin)	W65 (Stereo-Pin)
HILL LABS	Most Machines			W75 (3-Pin DIN-Pin)
JOHARI	AP-439 PLUS (2 channel)			W80 (5-Pin DIN-4 Pins)
KALADAR (KMI)	The Stim	AD6500	W65 (Stereo-Pin)	W65 (Stereo-Pin)
LSI	System II and IV, MS300, MS400	AD7600	W72 (4-Pin DIN-Pin)	W72 (4-Pin DIN-Pin)
MEDCO	Medcollator	AD6500	W65 (Stereo-Pin)	W65 (Stereo-Pin)
METRON	All models	AD4500	W47 (Banana-Pin)	W47 (Banana-Pin)
METTLER	Sonicator Plus 992, 994, 930	AD9500	W74 (3-Pin DIN-Pin) W73 (R/A 3-Pin DIN combo)	W74 (3-Pin DIN-Pin) W73 (R/A 3-Pin DIN combo)
	Sonicator 940	No current CMP solution for 8-Pin DIN		
	Sys*Stim 294, 226, 206, 220	AD9500	W74 (3-Pin DIN-Pin)	W74 (3-Pin DIN-Pin)
MONAD	Myomatic	AD4500	W44 (Banana-Pin)	
MULTIPLEX		AD8400	W80 (5-Pin DIN-4 Pins)	
NEMECTRON	Edit			W86 (5-Pin DIN-Pin)
NEUROCARE	Neurocare 1000			W89 (Mini DIN-Pin)
OMNISOUND	2500C	AD4500	W47 (Banana-Pin)	
ORTHION CORP	Sentry Mark V Model 9000		W2372 (3-Pin Amphenol)	
P.E.T. EUROTECH	Magnum 100	AD6500	W65 (Stereo-Pin)	
PHYACTON				W47 (Banana-Pin)
PHYSIO TECH.	Omnistim	AD6500	W65 (Stereo-Pin)	
PHYSIODYNAMIC	Therastim 1000	AD6500	W64 (Stereo-Pin)	
PHYSIOMED		AD8400	W80 (5-Pin DIN-4 Pins)	
PRECISION MICROCURRENT		AD4500	W44 (Banana-Pin)	W44 (Banana-Pin)
ProActive	Thera3		W26842WHWRB	
PROSTIM	ET 3000			
PureCare	EMS Primo Interferential 960			W71 (eliminates small leads)
Quattro	Tri-Socket to Dual Pins			W194 (Triple Socket-Pin)

Clinical Lead Wire Compatibility Chart

Manufacturer	Model	Unit made prior to 2000		Complaint Unit Wires
		Adaptor	Wire	
Quattro II	4 Channel with U/S			W80 (5-Pin DIN-4 Pins)
REBOX UNIT				
RICHMAR	HV 1000, 2000	AD8500 AD9500	W65 (Stereo-Pin) W75 (3-Pin DIN-Pin)	
	HV Combo	NA AD4500	W23 (3-Pin Amphenol-Pin) W44 (Banana-Pin)	
	HV 20		W26 (Red Active) W27 (Blk Active) W23 (Dispersive)	
	EVO Winner Series			W47 (Banana-Pin)
	IF 200	No Solution		
	IF 2000	AD8500	W65 (Stereo-Pin)	
	III, IIG	AD8500	W64 (Stereo-Pin)	
	RM6	AD8500 AD9500	W65 (Stereo-Pin) W75 (3-Pin DIN-Pin)	
	Theramini 2, 3C, 3P	AD8500 AD9500	W65 (Stereo-Pin) W75 (3-Pin DIN-Pin)	W47 (Banana-Pin)
	Theratouch 4.7, 7.7	AD8500 AD9500	W65 (Stereo-Pin) W75 (3-Pin DIN-Pin)	W47 (Banana-Pin)
	VI		W29 (Socket-Pin)	
ROSCOE MEDICAL	QUATTRO Combo Care			W80 (5-Pin DIN-4 Pins)
SCRIP	Super Stimm MF Model 4000	AD8500	W65 (Stereo-Pin)	W72 (4-Pin DIN-Pin)
	Super Stimm Plus			
STERNE	AcuTrain V	AD8500	W65 (Stereo-Pin) W56 (CDN Only)	W72 (4-Pin DIN-Pin)
	MyoTrain V	AD8500	W65 (Stereo-Pin) W56 (CDN Only)	W72 (4-Pin DIN-Pin)
	NeutoTrain III			W72 (4-Pin DIN-Pin)
SUN MEDISYS		AD6500	W66 (Gold Stereo-Pin)	
SUPERSTIM	NEUROTRAIN III			W72 (4-Pin DIN-Pin)
THERASTIM		AD6500	W65 (Stereo-Pin)	
US MEDICAL	Model 450		W23 (Amphenol 3 Pin-2 Pin)	
UNIPHY				W48 (Banana-Banana)
VECTRODYNE		No Solution		
WILLIAMS	Intertron 4000 etc.	AD6500	W65 (Stereo-Pin)	W81 (5-Pin DIN-Pin)
			W75 (3-Pin DIN-Pin)	

Notes:

Adaptors convert machine plugs to compliant wires

Convertors change patient end (ex. pins to banana plugs) for electrodes