

Adaptor Compatibility Chart

Original Machine Jack	Compliant Wire Machine End	Adaptor Needed	Compatible Wire	Common Machine Compatibility
Non-compliant Stereo	FDA Compliant Stereo	AD6500	W65	Dynatronics, Excel
Mono	FDA Compliant Stereo	AD8500	W65	Chattanooga, Richmar, Sterne
Mono	3 or 5 Pin DIN	AD9500	3-PIN W74 (R/A), W75 5-PIN W81	Mettler, Richmar
Non-compliant Stereo	4 Pin DIN	AD7600	W72	LSI
Non-compliant Stereo	3 or 5 Pin DIN	AD7500	3-PIN W74 (R/A), W75 5-PIN W81	Williams
Unsheathed Banana	Sheathed Banana	AD4500	W47	Amrex, Metron (sell adaptor with every W47 for a good fit)

*** TENSExperts@gmail.com ***

*** 613.800.2525 ***