

#	Order code	Compatible Devices	plug	Product Description	NTC	Hose material	user	Diameter	Appearance
1	PT400D-ES07	General YSI400 module		Disposable Temperature Probes 32in.2.3mm.Rectal/Esophageal, 25°C, 2.252KΩ,24pcs/pk	NTC/YSI400/2.252K@25 °C	PVC	Rectal/Esophageal	7fr	
2	PT400D-ES09	General YSI400 module		Disposable Temperature Probes 32in.3.0mm.Rectal/Esophageal, 25°C, 2.252KΩ,24pcs/pk	NTC/YSI400/2.252K@25 °C	PVC	Rectal/Esophageal	9fr	
3	PT400D-ES10	General YSI400 module		Disposable Temperature Probes 32in.3.2mm.Rectal/Esophageal, 25°C, 2.252KΩ,24pcs/pk	NTC/YSI400/2.252K@25 °C	PVC	Rectal/Esophageal	10fr	
4	PT400D-ES12	General YSI400 module		Disposable Temperature Probes 32in.4.0mm.Rectal/Esophageal, 25°C, 2.252KΩ,24pcs/pk	NTC/YSI400/2.252K@25 °C	PVC	Rectal/Esophageal	12fr	
5	PT10KD-ES07	Comen/10k module		Disposable Temperature Probes 32in.3.2mm,Rectal/Esophageal, 25°C, 10KΩ,24pcs/pk	NTC/10K@25°C	PVC	Rectal/Esophageal	7fr	
6	PT10KD-ES09	Comen/10k module		Disposable Temperature Probes 32in.3.2mm,Rectal/Esophageal, 25°C, 10KΩ,24pcs/pk	NTC/10K@25°C	PVC	Rectal/Esophageal	9fr	
7	PT10KD-ES10	Comen/10k module		Disposable Temperature Probes 32in.3.2mm,Rectal/Esophageal, 25°C, 10KΩ,24pcs/pk	NTC/10K@25°C	PVC	Rectal/Esophageal	10fr	
8	PT10KD-ES12	Comen/10k module		Disposable Temperature Probes 32in.4.0mm,Rectal/Esophageal, 25°C, 10KΩ,24pcs/pk	NTC/10K@25°C	PVC	Rectal/Esophageal	12fr	

9	PTPHD-ES07	Philips module		Disposable Temperature Probes 32in.2.3mm.Rectal/Esophageal, 25°C, 2.252KΩ,24pcs/pk	NTC/YSI400/2.252K@25 °C	PVC	Rectal/Esophageal	7fr	
10	PTPHD-ES09	Philips module		Disposable Temperature Probes 32in.3.0mm.Rectal/Esophageal, 25°C, 2.252KΩ,24pcs/pk	NTC/YSI400/2.252K@25 °C	PVC	Rectal/Esophageal	9fr	
11	PTPHD-ES10	Philips module		Disposable Temperature Probes 32in.3.2mm.Rectal/Esophageal, 25°C, 2.252KΩ,24pcs/pk	NTC/YSI400/2.252K@25 °C	PVC	Rectal/Esophageal	10fr	
12	PTPHD-ES12	Philips module		Disposable Temperature Probes 32in.4.0mm.Rectal/Esophageal, 25°C, 2.252KΩ,24pcs/pk	NTC/YSI400/2.252K@25 °C	PVC	Rectal/Esophageal	12fr	
13	PTOHD-ES07	GE Ohmeda Giraffe module		Disposable Temperature Probes 32in.3.2mm.Rectal/Esophageal, 25°C, 10KΩ,24pcs/pk	NTC/10K@25°C	PVC	Rectal/Esophageal	7fr	
14	PTOHD-ES09	GE Ohmeda Giraffe module		Disposable Temperature Probes 32in.3.2mm.Rectal/Esophageal, 25°C, 10KΩ,24pcs/pk	NTC/10K@25°C	PVC	Rectal/Esophageal	9fr	
15	PTOHD-ES10	GE Ohmeda Giraffe module		Disposable Temperature Probes 32in.3.2mm.Rectal/Esophageal, 25°C, 10KΩ,24pcs/pk	NTC/10K@25°C	PVC	Rectal/Esophageal	10fr	
16	PTOHD-ES12	GE Ohmeda Giraffe module		Disposable Temperature Probes 32in.4.0mm.Rectal/Esophageal, 25°C, 10KΩ,24pcs/pk	NTC/10K@25°C	PVC	Rectal/Esophageal	12fr	
17	PT400D-SK1	General YSI400 module		Disposable Temperature Probes 32in.Skin-surface, 25°C, 2.252KΩ 24pcs/pk	NTC/YSI400/2.252K@25 °C	PVC	Skin-surface	30mm	

18	PT400D-SK2	General YSI400 module		Disposable Temperature Probes 32in.Skin-surface, 25°C, 2.252KΩ 24pcs/pk	NTC/YSI400/2.252K@25 °C	PVC	Skin-surface	40mm	
19	PT400D-SK3	General YSI400 module		Disposable Temperature Probes 32in.Skin-surface, 25°C, 2.252KΩ 24pcs/pk	NTC/YSI400/2.252K@25 °C	PVC	Skin-surface	28.8*30mm	
20	PT10KD-SK1	Comen/10k module		Disposable Temperature Probes 32in.Skin-surface, 25°C,10KΩ 24pcs/pk	NTC/10K@25°C	PVC	Skin-surface	30mm	
21	PT10KD-SK2	Comen/10k module		Disposable Temperature Probes 32in.Skin-surface, 25°C,10KΩ 24pcs/pk	NTC/10K@25°C	PVC	Skin-surface	40mm	
22	PT10KD-SK3	Comen/10k module		Disposable Temperature Probes 32in.Skin-surface, 25°C,10KΩ 24pcs/pk	NTC/10K@25°C	PVC	Skin-surface	28.8*30mm	
23	PTPHD-SK1	Philips module		Disposable Temperature Probes 32in.Skin-surface, 25°C, 2.252KΩ 24pcs/pk	NTC/YSI400/2.252K@25 °C	PVC	Skin-surface	30mm	
24	PTPHD-SK2	Philips module		Disposable Temperature Probes 32in.Skin-surface, 25°C, 2.252KΩ 24pcs/pk	NTC/YSI400/2.252K@25 °C	PVC	Skin-surface	40mm	
25	PTPHD-SK3	Philips module		Disposable Temperature Probes 32in.Skin-surface, 25°C, 2.252KΩ 24pcs/pk	NTC/YSI400/2.252K@25 °C	PVC	Skin-surface	28.8*30mm	
26	PTDGAD-SK1	Drager Siemens module		Disposable Temperature Probes 32in.Skin-surface, 25°C, 2.252KΩ 24pcs/pk	NTC/YSI400/2.252K@25 °C	PVC	Skin-surface	30mm	

27	PTDGAD-SK2	Drager Siemens module		Disposable Temperature Probes 32in.Skin-surface, 25°C, 2.252KΩ 24pcs/pk	NTC/YSI400/2.252K@25°C	PVC	Skin-surface	40mm	
28	PTDGAD-SK3	Drager Siemens module		Disposable Temperature Probes 32in.Skin-surface, 25°C, 2.252KΩ 24pcs/pk	NTC/YSI400/2.252K@25°C	PVC	Skin-surface	28.8*30mm	
29	PTOHD-SK1	GE Ohmeda Giraffe module		Disposable Temperature Probes 32in.Skin-surface, 25°C,10KΩ 24pcs/pk	NTC/10K@25°C	PVC	Skin-surface	30mm	
30	PTOHD-SK2	GE Ohmeda Giraffe module		Disposable Temperature Probes 32in.Skin-surface, 25°C,10KΩ 24pcs/pk	NTC/10K@25°C	PVC	Skin-surface	40mm	
31	PTOHD-SK3	GE Ohmeda Giraffe module		Disposable Temperature Probes 32in.Skin-surface, 25°C,10KΩ 24pcs/pk	NTC/10K@25°C	PVC	Skin-surface	28.8*30mm	