

## EtCO2 Airway Adapter



### Mainstream Parts- EtCO2 Airway Adapter

| Ref#         | CO2 Airway Adapter        | Applicable patients | ET Tube Size | Deadspace |
|--------------|---------------------------|---------------------|--------------|-----------|
| <b>WK-01</b> | Single-Patient Disposable | Adult/Pediatric     | >4.0mm       | <5cc      |
| <b>WK-03</b> | Single-Patient Disposable | Infant/Neonatal     | ≤4.0mm       | <1cc      |

### Compatible Brands

Philips/ Respironics/Mindray/Drager/ZOLL

## Parts-Med Tech

C,2/F, Lilan Technology Building, Jiazitang Community,518107,

Fenghuang Street, Guangming District, Shenzhen, China

Email: [Marketing@parts-med.com](mailto:Marketing@parts-med.com)

Web: [www.parts-med.com](http://www.parts-med.com)

E-Mail: [Marketing@parts-med.com](mailto:Marketing@parts-med.com)

Web: [www.parts-med.com](http://www.parts-med.com)