

Holter ECG



Supports OEM, ODM customization

We are the source manufacturer

PSM+

Parts-Med Tech



Shenzhen Parts-Med Tech Co., Ltd

Address:


C,2/F, Lilan Technology Building, Jiazitang Community,518107,
Fenghuang Street, Guangming District, Shenzhen, China

Email: Marketing@parts-med.com

Web: www.parts-med.com

WhatsApp: +86 152 2023 7501

2025 catalogue



Compliance with

 AAMIEC53:2013

 IEC 62366-1:2015

 ISO 10993-5:2009

 ISO 10993-23:2021

 ISO 10993-10:2021

Wide range of compatible device models



PRODUCT SPECIFICATIONS

Materials

Conductor: Silver plated alloy

Shield: Silver plated alloy

Jacket: TPU

Materials Composition

Not made with BPA, natural rubber latex, PVC, DEHP or other phthalates

Operating environment

Temperature: 5° to 45°C (41° to 113°F)

Humidity: 5% to 95%, non-condensing

Storage

Temperature: -25° to 70°C (-13° to 158°F)

Humidity: 0% to 90%, non-condensing

Pressure: -700 kPa - 1060 kPa

Disposal considerations

Dispose as clinical waste, in accordance with hospital policy, local guidelines and regulations.

Packaging materials

Primary: PE (radiation protection)

Cleaning solutions

Ethyl alcohol: 70-80%

Isopropyl alcohol: 70-80%

Gutaraldehyde: 2% (pH 7.5-8.5)

Hydrogen peroxide: 1.5%

Sodium hypochlorite:

(5% household bleach): Minimum 1:500 dilution (minimum 100 ppm free chlorine) and maximum dilution 1:10

Do not sterilize or autoclave









Compliance











RoHS (Regulation of Hazardous Substances) compliant American Heart Association (AHA) standards







CE Mark

Product compatibility


see label

Models	OEM P/N	Ref No.	Connector	Leads	Patient end	Length/ Quantity
Burdick Cardiac Science X12		ZEH129		7	Snap/AHA	41.3 in/1piece
		15-048-0003		ZEL145	5	Snap/AHA
Biolight 850		15-048-0002		3	Snap/IEC	41.3 in/1piece
				ZEH010	5	Snap/AHA
Biomedical BI9800 / BI9000 EDAN SE-1201 Nihon Kohden Rac 3012		ZEH011		5	Snap/IEC	41.3 in/1piece
		ZEH031		7	Snap/AHA	41.3 in/1piece
		ZEH152		7	Snap/IEC	41.3 in/1piece
		ZEH047		10	Snap/AHA	41.3 in/1piece
		ZEH169		10	Snap/AHA	41.3 in/1piece
		ZEH170		10	Snap/IEC	41.3 in/1piece
		ZEH032		5	Snap/AHA	41.3 in/1piece
BI(Biomedical Instruments)		ZEH035		7	Snap/AHA	41.3 in/1piece
		ZEH036		10	Snap/AHA	41.3 in/1piece
		ZEH132		7	Snap/AHA	41.3 in/1piece
BI(Biomedical Instruments)		ZEH130		10	Snap/AHA	41.3 in/1piece
	14Pin BI(Biomedical Instruments)			ZEH033	7	Snap/AHA
		ZEH049		10	Snap/AHA	41.3 in/1piece
7PIN Beneware CT-08		ZEH043		7	Snap/IEC	41.3 in/1piece
10PIN Beneware		ZEH162		5	Snap/AHA	
		ZEH077		7	Snap/AHA	
		ZEH076		10	Snap/IEC	
13PIN		ZEH182		5	Snap/AHA	

Models	OEM P/N	Ref No.	Connector	Leads	Patient end	Length/ Quantity
Beneware		ZEH180		7	Snap/AHA	
		ZEH181		7	Snap/IEC	
		ZEH179		10	Snap/IEC	
HDMI19P Biocare IH- 12PLUS		ZEH025		5	Snap/AHA	41.3 in/1piece
		ZEH062		7	Snap/IEC	41.3 in/1piece
		ZEH054		10	Snap/IEC	41.3 in/1piece
17P Biocare IE10		ZEH055		10	Snap/IEC	41.3 in/1piece
		1				
Bainuodai, Creative		ZEH005		3	Snap/IEC	41.3 in/1piece
Biolight		ZEH069		5	Snap/AHA	41.3 in/1piece
		ZEH029		5	Snap/IEC	41.3 in/1piece
		ZEH104		10	Snap/IEC	41.3 in/1piece
CARDIOLINE		ZEH184		10	Banana/IEC	41.3 in/1piece
		ZEH186		10	Snap/IEC	41.3 in/1piece
		ZEL042		10	Snap/IEC	41.3 in/1piece
		ZEL043		7	Snap/IEC	41.3 in/1piece
		ZEL044		5	Snap/IEC	41.3 in/1piece
		ZEL050		5	Snap/IEC	41.3 in/1piece
	ZEL051		5	Cable - 1.5mm Din	41.3 in/1piece	
Choice		ZEE10001A D		3	Snap/IEC	41.3 in/1piece
DC3.5		ZEC066		3	Snap/IEC	41.3 in/1piece

Models	OEM P/N	Ref No.	Connector	Leads	Patient end	Length/ Quantity
Normal DIN 1.5		ZEL001		3	Clip/AHA	41.3 in/1piece
		ZEL002		3	Snap/AHA	41.3 in/1piece
		ZEL003		3	Clip/IEC	41.3 in/1piece
		ZEL004		3	Snap/IEC	41.3 in/1piece
		ZEL085		5	Clip/AHA	41.3 in/1piece
		ZEL086		5	Snap/AHA	41.3 in/1piece
		ZEL087		5	Clip/IEC	41.3 in/1piece
		ZEL088		5	Snap/IEC	41.3 in/1piece
		ZEL181		7	Clip/AHA	41.3 in/1piece
		ZEL182		7	Snap/AHA	41.3 in/1piece
		ZEL183		7	Clip/IEC	41.3 in/1piece
		ZEL184		7	Snap/IEC	41.3 in/1piece
		ZEL185		10	Snap/AHA	41.3 in/1piece
		ZEL188		10	Clip/IEC	41.3 in/1piece
		ZEL189		10	Snap/IEC	41.3 in/1piece
HDMI DMS300-4A		ZEH013		3	Snap/IE	41.3 in/1piece
		ZEH012		5	Snap/AHA	41.3 in/1piece
		ZEH037		7	Snap/AHA	41.3 in/1piece
		ZEH038		7	Snap/IEC	41.3 in/1piece
		ZEH052		10	Snap/IEC	41.3 in/1piece
15PIN DMS		ZEH198		4	Snap/IEC	41.3 in/1piece
		ZEH017		5	Snap/IEC	41.3 in/1piece
		ZEH200		7	Snap/AHA	41.3 in/1piece
		ZEH199		10	Snap/IEC	41.3 in/1piece
Philips		ZEH020		5	Snap/AHA	41.3 in/1piece
		ZEH021		5	Snap/IEC	41.3 in/1piece
Fukuda		ZEH03001A		3	Snap/IEC	41.3 in/1piece
		ZEC03001A S		4	Snap/IEC	41.3 in/1piece
		ZEC03001E S		5	Snap/IEC	41.3 in/1piece
HDMI Cardiocare		ZEH116		10	Snap/IEC	41.3 in/1piece
HDMI Cardiocare		ZEH117		10	Snap/IEC	41.3 in/1piece
GE 7pin		ZEH183		7	Snap/IEC	41.3 in/1piece









Models	OEM P/N	Ref No.	Connector	Leads	Patient end	Length/ Quantity	
GE 10pin		ZEH095		10	Snap/IEC	41.3 in/1piece	
		ZEH014		5	Snap/AHA	41.3 in/1piece	
GE 16P Kenz-Cardy 301		ZEH015		5	Snap/IEC	41.3 in/1piece	
		ZEH155		5	Clip/AHA	41.3 in/1piece	
		ZEH027		5	Clip/IEC	41.3 in/1piece	
		ZEH039		7	Snap/AHA	41.3 in/1piece	
		ZEH040		7	Snap/IEC	41.3 in/1piece	
		ZEH042		7	Clip/AHA	41.3 in/1piece	
		ZEH009		7	Clip/IEC	41.3 in/1piece	
	GE 8P			ZEH155		5	Clip/AHA
		ZEL114		5		Snap/AHA	41.3 in/1piece
		ZEL115	6	Clip/AHA		41.3 in/1piece	
		ZEL116	6	Snap/AHA		41.3 in/1piece	
GE 20PIN				5	Snap/IEC	41.3 in/1piece	
		ZEH106		7	Snap/IEC	41.3 in/1piece	
Vales & Hills, Borsam, Zoncare		ZEH176		5	Snap/AHA	41.3 in/1piece	
		ZEH177		5	Snap/IEC	41.3 in/1piece	
		ZEH165		7	Snap/AHA	41.3 in/1piece	
		ZEH166		7	Snap/IEC	41.3 in/1piece	
		ZEH067		10	Snap/AHA	41.3 in/1piece	
		ZEH175		10	Snap/IEC	41.3 in/1piece	
Vales & Hills		ZEH071		10	Snap/AHA	41.3 in/1piece	
Vales & Hills, Borsam		ZEH161		5	Snap/AHA	41.3 in/1piece	
		ZEH159		7	Snap/AHA	41.3 in/1piece	
		ZEH073		10	Snap/AHA	41.3 in/1piece	
Evergrande TYPE-C		ZEH134		7	Snap/IEC	41.3 in/1piece	
		ZEH135		10	Snap/IEC	41.3 in/1piece	
6PIN HL medical		ZEH001		3	Snap/AHA	41.3 in/1piece	
		ZEH016		5	Snap/AHA	41.3 in/1piece	



Models	OEM P/N	Ref No.	Connector	Leads	Patient end	Length/ Quantity
7PIN HL medical		ZEH002		3	Snap/AHA	41.3 in/1piece
		ZEL045		3	Clip/AHA	41.3 in/1piece
HRD		ZEL046		3	Snap/AHA	41.3 in/1piece
		ZEL047		3	Clip/IEC	41.3 in/1piece
		ZEL048		3	Snap/IEC	41.3 in/1piece
		ZEL120		4	Snap/AHA	41.3 in/1piece
		ZEL129		5	Clip/AHA	41.3 in/1piece
		ZEL130		5	Snap/AHA	41.3 in/1piece
		ZEL131		5	Clip/IEC	41.3 in/1piece
		ZEL132		5	Snap/IEC	41.3 in/1piece
HSM		ZEH019		10	Snap/AHA	41.3 in/1piece
HSM HDMI		ZEH089		10	Snap/AHA	41.3 in/1piece
HanWei、Datawe		ZEH090		10	Snap/IEC	41.3 in/1piece
HX HDMI		ZEH124		10	Snap/AHA	41.3 in/1piece
CHANCEFINE		ZEH109		10	Snap/IEC	41.3 in/1piece
CHANCEFINE		ZEH114		10	Snap/IEC	41.3 in/1piece
Contec		ZEH193		10	Snap/AHA	41.3 in/1piece
		ZEH112		10	Snap/IEC	41.3 in/1piece
Comen		ZEH024		5	Snap/AHA	41.3 in/1piece
		ZEH034		5	Snap/IEC	41.3 in/1piece
MINI HDMI KeMeiSi		ZEH098		10	Snap/IEC	41.3 in/1piece

Models	OEM P/N	Ref No.	Connector	Leads	Patient end	Length/ Quantity
Carewell、Lepu		ZEH150		10	Snap/AHA	41.3 in/1piece
16PIN Edan		ZEH026		5	Snap/AHA	41.3 in/1piece
		ZEH045		7	Snap/AHA	41.3 in/1piece
		ZEH044		7	Snap/IEC	41.3 in/1piece
		ZEH093		10	Snap/AHA	41.3 in/1piece
		ZEH094		10	Snap/IEC	41.3 in/1piece
20PIN Edan		ZEH056		10	Snap/IEC	41.3 in/1piece
		ZEH096		10	Snap/AHA	41.3 in/1piece
IT20 Edan		ZEH147		3	Snap/IEC	41.3 in/1piece
		ZEH146		5	Snap/IEC	41.3 in/1piece
Edan	SA-20	ZEE10001E S		10	Snap/IEC	41.3 in/1piece
Labtech Kft.	EC-3H			7	Snap/IEC	41.3 in/1piece
	EC-12H			10	Snap/IEC	41.3 in/1piece
HDMI19P Landcom		ZEH028		5	Snap/AHA	41.3 in/1piece
		ZEH102		10	Snap/AHA	41.3 in/1piece
10PIN Landcom		ZEH097		10	Snap/AHA	41.3 in/1piece
7PIN Mbelec		ZEH006		3	Snap/AHA	41.3 in/1piece
HDMI MEDEX		ZEH068		10	Snap/IEC	41.3 in/1piece
26PIN MEDEX		ZEH158		10	Snap/IEC	41.3 in/1piece
Meditech		ZEH174		7	Snap/IEC	41.3 in/1piece
		ZEL039		10	Snap/IEC	41.3 in/1piece
H4 MGY		ZEH105		10	Snap/AHA	41.3 in/1piece

Models	OEM P/N	Ref No.	Connector	Leads	Patient end	Length/ Quantity
		ZEH163		10	Snap/IEC	41.3 in/1piece
MGY		ZEH160		10	Snap/AHA	41.3 in/1piece
		ZEH022		5	Snap/AHA	41.3 in/1piece
		ZEH023		5	Snap/IEC	41.3 in/1piece
10PIN Mortara		ZEH066		10	Snap/IEC	41.3 in/1piece
10PIN Mortara		ZEH143		10	Snap/IEC	41.3 in/1piece
10PIN Mortara		ZEL03001A C		10	Snap/IEC	41.3 in/1piece
MILLION HDMI		ZEH070		10	Snap/IEC	41.3 in/1piece
美浩 MH		ZEH072		10	Snap/IEC	41.3 in/1piece
7PIN Mindray		ZEH007		3	Snap/AHA	41.3 in/1piece
		ZEH008		3	Snap/IEC	41.3 in/1piece
		ZEH030		5	Snap/AHA	41.3 in/1piece
12PIN Mindray		ZEH145		3	Snap/AHA	41.3 in/1piece
		ZEH189		3	Snap/IEC	41.3 in/1piece
		ZEH144		5	Snap/AHA	41.3 in/1piece
		ZEH123		6	Clip/IEC	41.3 in/1piece
		ZEH127		6	Snap/IEC	41.3 in/1piece
North East		ZEH122		5	Snap/AHA	41.3 in/1piece
		ZEH156		7	Snap/AHA	41.3 in/1piece
Norav		ZEL03001E C		3	Snap/IEC	41.3 in/1piece

Models	OEM P/N	Ref No.	Connector	Leads	Patient end	Length/ Quantity
		ZEL03001E S		4	Snap/IEC	41.3 in/1piece
		ZEH048		5	Snap/IEC	41.3 in/1piece
		ZEH053		10	Snap/IEC	41.3 in/1piece
		ZEH108		10	Snap/AHA	41.3 in/1piece
15P PengYang		ZEH191		9	Snap/IEC	41.3 in/1piece
		ZEH103		10	Snap/AHA	41.3 in/1piece
11P PengYang		ZEH125		9	Snap/IEC	41.3 in/1piece
		ZEH164		10	Snap/IEC	41.3 in/1piece
4+6 Shanghai Quntian		ZEH088		10	Snap/IEC	41.3 in/1piece
HDMI Shanghai Quntian		ZEH082		10	Snap/IEC	41.3 in/1piece
BH-812D Nihon Kohden		ZEH196		10	Snap/AHA	41.3 in/1piece
		ZEH107		10	Snap/IEC	41.3 in/1piece
BH-853D Nihon Kohden		ZEH197		5	Snap/IEC	41.3 in/1piece
RuiDi、FuSiTe		ZEE10001E B		3	Snap/IEC	41.3 in/1piece
Rencare		ZEH188		5	Snap/AHA	41.3 in/1piece
		ZEH172		7	Snap/IEC	41.3 in/1piece
		ZEH099		10	Snap/IEC	41.3 in/1piece
3RAY				10	Snap/IEC	41.3 in/1piece
11PIN JincoMed、 Contec medwinghealthc are		ZEH058		10	Snap/IEC	41.3 in/1piece
11PIN JincoMed ECGMACMedical		ZEH060		10	Snap/IEC	41.3 in/1piece

Models	OEM P/N	Ref No.	Connector	Leads	Patient end	Length/ Quantity
14PIN JincoMed		ZEH100		10	Snap/IEC	original
Times、3RAY、 JinHao		ZEH085		10	Snap/AHA	41.3 in/1piece
		ZEH086		10	Snap/AHA	41.3 in/1piece
Safezx		ZEH126		10	Snap/IEC	41.3 in/1piece
TJJY		ZEH187		3	Snap/AHA	41.3 in/1piece
		ZEH192		4	Snap/AHA	41.3 in/1piece
TR		ZEH092		10	Snap/AHA	41.3 in/1piece
		ZEH138		7	Snap/AHA	41.3 in/1piece
TianZe		ZEH120		5	扣式, Snap	41.3 in/1piece
		ZEH119		7	扣式, Snap	41.3 in/1piece
		ZEH118		10	扣式, Snap	41.3 in/1piece
		ZEE097		10	Banana/AH A	41.3 in/1piece
HDMI Zoncare		ZEH167		10	Snap/IEC	41.3 in/1piece
10P Welch Allyn Protocol		ZEH087		5	Snap/AHA	41.3 in/1piece
DB15 Welch Allyn Protocol		ZEC03001A C		10	Snap/AHA	41.3 in/1piece
We for Health		ZEH171		10	Snap/AHA	41.3 in/1piece
		ZEH190		10	OT	41.3 in/1piece

Models	OEM P/N	Ref No.	Connector	Leads	Patient end	Length/ Quantity
Avantgrade		ZEH064		10	Snap/IEC	41.3 in/1piece
10PIN Schiller		ZEH046		10	Snap/AHA	41.3 in/1piece
		ZEL03001A S		10	Banana/IEC	41.3 in/1piece
7PIN Schiller		ZEL119		5	Clip/IEC	41.3 in/1piece
PCB-8P Schiller		ZEC273		4	Snap	41.3 in/1piece
		ZEC274		6	Snap	41.3 in/1piece
8P Schiller		ZEH153		5	Snap	41.3 in/1piece
		ZEC03001E C		7	Snap	41.3 in/1piece
XM		ZEH133		10	Snap/AHA	41.3 in/1piece
Yocaly, LEPU		ZEH168		6	Snap/AHA	41.3 in/1piece
HDMI YiKang		ZEE10001A B		10	Snap/AHA	41.3 in/1piece
		ZEE10001A S		5	Snap/AHA	41.3 in/1piece

Models	OEM P/N	Ref No.	Connector	Leads	Patient end	Length/ Quantity
YTYZ				10	Snap/IEC	41.3 in/1piece
		ZEH173		10	Snap/AHA	41.3 in/1piece
26P BeiJingYiSi		ZEH195		6	Snap/IEC	41.3 in/1piece
HDMI Biox		ZEH079		10	Snap/AHA	41.3 in/1piece
		ZEH080		7	Snap/AHA	41.3 in/1piece
		ZEH194		5	Snap/AHA	41.3 in/1piece
		ZEH083		10	Snap/AHA	41.3 in/1piece
中健科仪 HDMI Biox		ZEH110		10	Snap/AHA	41.3 in/1piece
3PIN Careshine		ZEH003		2	Snap/AHA	41.3 in/1piece
26PIN Careshine		ZEH136		10	Snap/AHA	41.3 in/1piece
ZKXHL		ZEH018		5	Snap/AHA	41.3 in/1piece
		ZEH081		10	Snap/AHA	41.3 in/1piece
中网心电 Zwxd		ZEH148		10	Snap/AHA	41.3 in/1piece
Кабель		ZEL041		10	Snap/IEC	41.3 in/1piece
		ZEL040		7	Snap/IEC	41.3 in/1piece